



Hu-Friedy's

CLEANING MONITORS

Quality Assured

Proper instrument cleaning is the most important step to instrument reprocessing. If an instrument is not properly cleaned, it cannot be effectively sterilized. Take the steps to assure that your cleaning process achieves the highest quality results by incorporating Hu-Friedy Cleaning Monitors into your instrument reprocessing protocols. Using multi-parameter testing, Hu-Friedy Cleaning Monitors are designed to give you a clear interpretation of your cleaning process.



POINTS OF PERFORMANCE



Colored test soil on strip is designed to parallel the removal of blood and bioburden from an instrument surface

Cleaning Monitor Holder acts as a hinged instrument measuring impingement in washer-disinfectors

Reusable stainless steel holder secures monitor in place in order to provide accurate result



Ideal Cleaning Monitor Usage

WASHER-DISINFECTOR

Run a monitor test in the morning to release the equipment for use and minimum of one monitor per load for load release.

ULTRASONIC CLEANING

Monitor once daily with an empty load for a machine release and periodically throughout the day with instruments for a load release, if needed.

How the best perform

Hu-Friedy

CLEANING MONITORS

CLEAR & SIMPLE RESULTS INTERPRETATION

Removal of all test soil indicates effective cleaning

	Washer-Disinfecter	Ultrasonic
Fail (Unexposed)		
Fail		
		
		
		
Pass		

See Trouble Shooting Charts (N1736 and N1692) for possible reasons for results and corrective actions.



Washer-Disinfecter Cleaning Monitor, 50 PCS
| [IMS-1200W](#)



Ultrasonic Cleaning Monitor, 50 PCS
| [IMS-1200U](#)



Cleaning Monitor Holder, 1 EA
| [IMS-1200H](#)

PARAMETERS TESTED & COMMON FAILURES

	Washer-Disinfecter	Ultrasonic
Parameters tested by the cleaning monitors	• Time • Temperature • Detergent concentration • Spray arm function • Enzyme soak	• Cavitation • Time • Temperature • Detergent
Common Cleaning Process Failures	• Inadequate water spray/impingement • Clogged spray arms • Overloading • Instrument shadowing • Inadequate detergent dosing • Improper detergent dosing • Poor water quality	• Ineffective cavitation • Overloading • Insufficient time, temperature and/or detergent

Visit us online at [Hu-Friedy.com](#)

Call 1-800-Hu-Friedy or contact your authorized Hu-Friedy representative for more information.
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. IMS-689/1015

