🇺🇵 MÉTODOS RECOMENDADOS PARA O PROCESSAMENTO DE IMS

LIMPEZA

Por ultra-sons

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🇧🇷 RECOMMENDED IMS PROCESSING METHODS

CLEANING

Ultrasonic

- Follow manufacturer’s directions for proper ultrasonic cleaner/disinfection operation.
- Hu-Friedy recommends a minimum run time of 16 minutes. Actual cleaning/disinfection time may vary.
- Clean Hu-Friedy instruments or cassettes in pH neutral enzymatic cleaner (such as Hu-Friedy’s Enzymax™) or multipurpose detergent/disinfectant.
- After ultrasonic cleaning, rinse cassette thoroughly with demineralized water to remove detergent or particles.

Instrument/Automated Washer

- Follow manufacturer’s directions for proper cleaning in an automated washer.
- Clean Hu-Friedy instruments or cassettes in a pH neutral detergent/disinfectant only.

Thermal Disinfectior

- Follow manufacturer’s directions for proper cleaning and disinfection operation.
- Clean Hu-Friedy instruments or cassettes in a pH neutral detergent/disinfectant only.

INSPECTION/PREPARATION

Open cassette and inspect instruments, making sure all debris has been removed from instruments and cassette during ultrasonic cleaning. Check for broken or missing instruments. Instruments should be dry before they are wrapped or sterilized to help prevent wet packing. Add an IMS™ Indicator or IMS™ Integrator Strip if you are using steam sterilization. Wrap cassettes with IMS™ Autoclave Wrap (depending on cassette size). Seal with IMS™ Monitor Tape. IMS™ Monitor Tape allows quick visual identification of wrapped cassette and confirmation that items have passed through steam sterilization. Other packaging options include Hu-Friedy Bagette™ Self-Sealing Sterilization Pouches or Hu-Friedy IMS™ Sterilization Containers. Do not use Hu-Friedy sterilization packaging products (IMS™ Wrap, Bagette pouches, or Sterilization Containers) in dry heat sterilization.

STERILIZATION

For all Hu-Friedy Cassettes - Steam autoclave. Do not exceed 350°F/177°C. Steam sterilize for at least 4 minutes at 273°F/134°C. Recommend 30 minute dry time after sterilization cycle. Follow manufacturer’s directions for proper sterilizer operation. Place wrapped cassettes in sterilizer making sure to keep cassettes slightly separated. Follow sterilizer manufacturer’s directions for “wrapped instruments” to determine length of cycle, temperature, and drying instructions. Remove cassettes. Confirm that monitor tape on the outside of the wrap or indicator/Integrator inside the package has changed color (if using a steam autoclave).

	www.hufriedygroup.com/Reprocessing/1-800-Hu-Friedy; info@hu-friedy.com
----------------------------	--

PRIOR TO USE

Inspect the IMS™ Indicator or IMS™ Integrator Strip to verify cassette has been exposed to proper autoclave parameters of time and temperature. Stripes should equal or exceed color standard on the strip.

AFTER USE

Place lid on cassette and latch shut. Note: Align hinges at the back of the cassette for proper closing, with latch(es) in the vertical/open position. When locking the latch(es), apply light pressure to the lid of the cassette then turn latch(es) until locked. Do not overload the cassette. Resin cassettes will take little or no pressure to close. The rail system in the Infinity Series and Signature Series™ Cassettes will require slight pressure to close the cassette. If you encounter resistance when turning the latch to the closed position, stop. Check to see that all instruments are properly placed. Forcing the cassette to close under stress may hinder future performance.

RECOMMENDATIONS

Systematic biological monitoring of the sterilizer with sporicidal tests is recommended by hospital standards. This is the only monitoring method that ensures items being processed in a sterilizer are sterile upon completion of a cycle. Most sterilizer manufacturers offer IMS™ Cassette Racks for their units. To ensure proper sterilization of cassettes, ask your dealer or sterilizer manufacturer for a cassette rack.

PRECAUTIONS

Be sure to clean and maintain ultrasonic and sterilizer per manufacturer’s recommendations. Do not expose cassettes/instruments to temperatures above 350°F/177°C. After use at chairside, wipe off cement residue, glass ionomer, blood and/or other foreign matter to promote faster and more thorough cleaning. Use only with Hu-Friedy IMS related products. Certain chemicals and detergent/disinfectants can adversely affect the cassette causing discoloration and/or corrosion. If you are using other solutions, such as film processing in your ultrasonic, remember to change back to a multipurpose detergent prior to processing instruments.

CAUTION

All Hu-Friedy Stainless Steel Cassettes should not come in contact with the following chemicals:

Activated Chlorous Acid/Chlorine Dioxide	Aluminum Chloride
Barium Chloride	Bichloride of Mercury
Calcium Chloride	Carbolic Acid
Chlorinated Lime	Citric Acid
Dakin’s Solution	Ferrous Chloride
Lysol	Mercuric Chloride
Mercury Salts	Phenol
Potassium Permanganate	Potassium Thiocyanate
Sodium Hypochlorite	Stannous Chloride
Tartaric Acid	Aqua Regia
Ferric Chloride	Sulphuric Acid
Hydrochloric Acid	Iodine

For more detailed cleaning and care information about Hu-Friedy instruments and IMS™ Cassettes, please call Hu-Friedy Customer Care at 1-800-Hu-Friedy - for Europe: Tel. 00800-48374339 or view the Instrument Reprocessing document in the Instrument Care section under Resources of the Hu-Friedy website, HuFriedyGroup.com or for Europe: Hu-Friedy.eu.

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