STERILE CONTAINER SYSTEMS

INSTRUCTIONS FOR USE FOR NORTH AMERICA ONLY

With the purchase of this sterilization container you have chosen a high quality product. The correct handling and use are described in the following.

In order to keep risks and unnecessary burdening as low as possible for the patient, the user and third parties, we request that prior to first use you carefully read through these instructions for use and keep them for further reference.

Article numbers and descriptions

IMDINCO0M	DIN-Container for max. 1 container cassette Dimensions: 312 x 190 x 50mm
IMDINCO1M	DIN-Container for max. 1 container cassette Dimensions: 312 x 190 x 65mm
IMDINCO2M	DIN-Container for max. 2 container cassettes Dimensions: 312 x 190 x 100mm
IMDINCO3M	DIN-Container for max. 3 container cassettes Dimensions: 312 x 190 x 135mm

1 FIELD OF APPLICATION

The Sterilization Container System is a sterile packaging system for dental instruments that utilizes an established filter technology, proved materials and design properties which make our Container Systems reliable. It is a reusable device and features an assortment of sizes that provide an effective packaging method for sterilization, storage and transportation of instruments by healthcare providers. This container system is suitable for use in pre-vacuum Steam Sterilizers. All Containers are available with perforated bottom and perforated lids. The containers must only be used with appropriate single-use paper filters provided by Hu-Friedy.

2 INTENDED USE

The Sterilization Container System is intended to allow sterilization of the enclosed dental device and also maintain sterility during transport and storage for 30 days. The system consists of different dental container sizes. All models are available with perforated lids and perforated bottoms. All models are intended to be used with Hu-Friedy single use paper filters. The containers are reusable devices designed to be used with the following sterilization cycle parameters:

Pre Vacuum cycle:

4 minutes
132°C (270°F)
Drying Time minimum 20 minutes

3 minutes
135°C (273°F)
Drying Time minimum 16 minutes

Loading: Metal dental instruments (scalars, scissors, clamps, forceps, etc).

2.1 Contraindications

- Sterilization Containers have not been validated for sterilization of instruments with channels or lumens, compressed air-driven instruments or power systems.
- An application of gas sterilization is not suitable because it is carried out with Formaldehyde or Ethylene oxide. The containers are not made for it.
- Do not externally stack the containers!

3. HANDLING

3.1 General

The container is made of aluminum alloy with an anodized oxide surface which resists corrosion. Abrasive cleaners, metal brushes or abrasive cleaning pads can cause permanent damage to the container surface and therefore must not be used. Warranty exclusions will result if these instructions are not followed.

The product and accessories must be used only by qualified or trained and experienced personnel, in order to prevent damage to the containers, closing devices, seals and sterile filters.

The closing device can be locked with a security seal, which has to be broken when opening.

3.2 First use

- The new container needs proper cleaning before first use.
- Use included automated cleaning program prior to first use.
- In the machine to use there should be a neutral detergent (pH 7).
- After processing in the cleaning equipment, we recommend sterilization of the container according to the sterilization parameters defined in this instruction.
- We also recommend the maintenance of all moving parts on the container with instrument oil.
- After cleaning, a new filter must be inserted (see 3.6 filter change)

3.3 Preparation before cleaning

- Separate lid and bottom
 - Remove the devices from the inside of the container (baskets, instruments, etc.)
 - Take off the filter holders inside the lid and, if existent, in the bottom (in case of perforated bottom containers)
 - Discard the disposable filter(s)
 - Remove disposable locks and indicators
- Note: all paper-filters are for single use only and should be disposed after each processing-cycle

3.4 Cleaning

This device is only validated for the automated cleaning procedure described in section 3.4.2.

Note: Improper cleaning and disinfection can lead to corrosion and stress cracks. Therefore follow the specifications the cleaning agent's manufacturer recommends. Detergents must be alkaline-, sodium- and carbonate-free, of neutral pH (7) and mild. Only fully deionized water is recommended for reprocessing the product.

- Before every use, the container has to be cleaned.
- Removal containers have to be cleaned each time after use.

3.4.1 Mechanical cleaning

- Mechanical cleaning of the containers is only recommended if the washer/disinfector has a special program for aluminum containers. When mechanically cleaning the containers, it is absolutely essential to observe the specifications of the washer/disinfector manufacturers and their instructions for the cleaning of containers.
- Only use neutral cleaners for cleaning. Do not use any cleaning solutions containing soda or caustic soda.
- Use neutral cleaners only with aluminum containers, which have been expressly approved by the manufacturer for the cleaning of aluminum containers. The exact dose of the cleaner is to be used as specified by the manufacturing company. With an optimization of the program these products are also suitable for cleaning surgical instruments.
- No additional use of acidic neutralizers.
- Fully demineralized water has to be used for final rinsing, as the salts in the water during the subsequent sterilization can also lead to spotting
- The cleaning system must be constructed for the cleaning of sterilization containers. That is specially indicated for the safe placing in the washing baskets and the positioning of the spray nozzles or -units.
- The container lid and bottom need to be cleaned separated from each other.
- The container bottom needs to be placed upside down in the washer/disinfector to avoid the collection of water.
- The inside of the container lid should face the machine-bottom with the latching mechanism folded into the inner part of the lid.
- All component parts of the container (bottom, lid, filter holders...) should be disassembled and placed in the washer/disinfector baskets that are specially designed for containers and accessories.
- After the washing/disinfection program is finished, the container and accessories need to be dried with a soft dry cloth or by (medical) compressed air.
- Contaminations that cannot be removed in a normal cleaning cycle independent of the process (sticky labels, indicator strips, writing), can be removed with eloxal cleaning agent. After this special treatment, the containers should be cleaned as usual.

3.4.2 Recommended Cleaning procedure

Based on validation studies, we recommend the following cleaning procedure by using qualified equipment for automated cleaning:

- 1 min prewash with cold water
- 3 minutes cleaning with Mediclean 0,5% (Dr. Weigert) at 45° C
- Neutralization with demineralized water

3.5 Inspection

The sterilization containers have to be checked every time before use to ensure their correct functionality. If the user detects any visible signs of wear, such as discoloration, cracking, peeling or flaking, the containers must be discontinued and prevented from being reused.

3.6 Filter change

After placing the filters over the perforated areas on the inside of the lid and on the bottom of the container, the filter holders have to be pressed into place until they snap into position. Bottoms and lids should only be assembled with Hu-Friedy filter holders.

- The disposable paper filters have to be changed before every new sterilization cycle.
- Only Hu-Friedy filters can ensure suitability and exactness of fit.
- Warranty can only be accepted when original Hu-Friedy filters have been used.

WARNING!

Only original Hu-Friedy component parts such as lids, bottoms, filters, seals and filter holders should be combined with each other so that the container sealing and function as a microbial barrier is not compromised.

4 LOADING

The overall weight for loading containers should not exceed the following loads:

MODEL	DIMENSION	Maximum Recommended Load in kg (lb)	
		INSTRUMENTS	TEXTILES
IMDINCO0M	312 x 190 x 50 mm	0.6 (1.3)	0.5 (1.1)
IMDINCO1M	312 x 190 x 65 mm	0.9 (1.9)	0.7 (1.5)
IMDINCO2M	312 x 190 x 100 mm	1.5 (3.3)	1.2 (2.6)
IMDINCO3M	312 x 190 x 135 mm	1.8 (3.9)	1.4 (3.0)

PRECAUTION!

Complex instruments like Endoscopes, instruments with a lumen, compressed air-driven instruments or power systems and instruments with cannulae are to be prepared according to the manufacturer's instructions for sterilization and to be sterilized accordingly.

Small baskets, trays, other types of accessories, especially with cover or lids, should only be used with the sterilization container, if the sterilization container has been specifically designed and tested for that purpose.

WARNING!

When using non-absorbent tray liners (e.g., plastic/silicone-fingered organizing mats) for the containers this can cause condensation to pool. Instead of this use moisture absorbing mats if necessary.

Always attach the lid to the bottom via the latching mechanism before placing the container in the sterilizer. Otherwise the content of the container will become unsterile as soon as the sterilizer door will be open.

5 PLACEMENT IN THE STERILIZER

PRECAUTION!

Never wrap the container in any kind of outer packaging. Never cover the perforated area of the container with any kind of foil packaging or similar during sterilization, because this will block the air and steam flow through the perforation. The result will be a vacuum damage due to insufficient pressure venting and the container content will not be sterilized.

During loading and unloading of the sterilizer and during transport, the sterile container must always be carried by the bottom and never by the lid.

6 PROCESSING

- Follow the time and temperature recommendations of the sterilizer manufacturer for every chosen sterilization cycle.
- To minimize condensation inside the container, leave container on container carts until cool enough to handle.
- Always check if the sterilization has been successful by observing indicators and the container latch (must be closed).
- Wet filters or visible moisture within a load may compromise the sterile barrier and should not be considered sterile when removed from the sterilizer.

7 STORAGE

- Temperature 15 – 26°C / 59-79°F
- Air humidity 30 – 50%
- Air pressure, normal atmospheric pressure

The sterilization containers were tested for a storage period of 30 days by spraying *Bacillus subtilis* – Spore suspension in an unventilated room.

8 MAINTENANCE

- Cover seals should not be treated with spray, oil or solvents. For cleaning and maintenance simply wipe off occasionally with a damp cloth.
- If any damage on the container or its parts is noticed the container must be discontinued and prevented from being reused.
- Replacement parts can be obtained from Hu-Friedy:
 - Disposable Paper Filters (100 per pack) - IMCO-FIL3M
 - Seals for Container (100 per pack) - IMCOM-SEAL
 - Indicator Card for Container - (250 per pack) - IMCO-INDNA

9 MATERIALS

Sterilization containers are made of anodised aluminium alloy.

10 CONTAINER DESIGN

In order to ensure sterile safety, tests were carried out by an independent and accredited test laboratory. The purpose of these trials was to validate a sterilization process for the reusable sterilization container with steam.

HU-FRIEDY DISCLAIMS LIABILITY FOR CUSTOMERS FAILURE TO FOLLOW INSTRUCTIONS FOR USE.

 Charge number, identification assigned to a particular quantity or lot of material from the manufacturer

 Read the instruction for use

Manufactured for:
Hu-Friedy Mfg. Co., LLC
2501 Sanders Rd.
Northbrook, IL 60062 | USA
1-800-Hu-Friedy | HuFriedyGroup.com

 **Hu-Friedy Europe LLC & Co. KG**
Bogenstrasse 6 - 8
78576 Emmingen - Liptingen
Germany

 **MDSS CH GmbH**
Laurenzenvorstadt 61
5000 Aarau
Switzerland

 **MDSS-UK RP Ltd.**
Parkway House, Palatine Rd,
Northenden, Wythenshawe,
Manchester M22 4DB
United Kingdom

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved.

HuFriedyGroup
The Best In Practice

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©2026 Hu-Friedy Mfg. Co., LLC. Tous droits réservés.