# Accutron™ Axess™ Scavenging Circuits

## INSTRUCTIONS FOR USE

Axess™ RFS™ Scavenging Circuits (52006, 53006)


Figure 1

Figure 2

**A** 53034-9, 53034-12, 53034-16, 53035-9, 53035-12, 53035-16, 53037-9, 53037-12, 53037-16, **B** 28624-FRU, **C** 28417-FRU, **CC** 28463-FRU, **DD** 28559-FRU


Figure 3


Figure 4


Figure 5

**Axess™ Inline Circuits (52005, 53005)**



Figure 6

Figure 7

**A** 53034-9, 53034-12, 53034-16, 53035-9, 53035-12, 53035-16, 53037-9, 53037-12, 53037-16, **DD** 28625-FRU, **C** 28417-FRU, **D** 28530-FRU, **E** 25591 **F** 26254


Figure 8


Figure 9

Figure 10


Figure 11

## System Setup

1. Remove the reservoir bag from the flowmeter and install the plug  per Figures 3, 5 and 8. **Note:** Plug will pop off if the gas flow line is kinked or blocked.
2. Install the mask on the circuit per Figure 7. It does not matter which line goes into which side of the mask.
3. Connect right angle adapter  to bagtee outspout per Figures 3, 5 and 8.
4. **RFS™ Scavenging Circuits:** Connect the Axess™ RFS Circuit to the RFS Outlet per Figure 3. Fittings are sized to only fit one way.  
**Accu-Vac™ Vacuum Controller:** Connect the Axess RFS Circuit to Accu-vac Controller per Figure 5. Note that the Accuvac Adapter  (28559-FRU) is required. Do not use with standard fitting.  
**In-line:** Connect the Axess In-line Circuit per Figure 8. **Note:** the vacuum connection must be made on the gauge side. Connect the spiral vacuum tube from the controller to a vacuum source in the operator. Some options are shown in Figures 9 and 10.  
**Caution:** Only one vacuum controller can be used at one time. Do not connect an In-line controller to the RFS Circuit or Accu-vac Controller.
5. Place the circuit and mask on the patient per Figure 2 making sure that the slide is snug under the chin to minimize mask movement.

## Scavenging Vacuum Operation

6. **RFS™ Scavenging Circuits:** Adjust the flow on the RFS Vacuum Controller using knob  so that the needle is in the green zone of the gauge  per Figure 4.  
**Accu-Vac™ Vacuum Controller:** Adjust the flow on the Accu-vac Controller to fully open (horizontal) using knob  per Figure 5. Disregard light colors as they do not function with Axess™ Circuits.  
**In-line:** Adjust the flow on the In-line controller using knob  until the needle is in the green zone per Figure 4.  
 It is recommended that the gauge is checked periodically during the procedure as facility vacuum pressure can change from time to time.

## Gas Flow Operation

**Digital Ultra™ Flowmeters:** Turn on the unit. It should start at TOTAL FLOW of 5 lpm with 0% N<sub>2</sub>O. Gradually increase TOTAL FLOW and %N<sub>2</sub>O until desired effect is reached. Use the Titration method when adding N<sub>2</sub>O to avoid negative patient reactions. See note 8 below. At the end of the procedure, reduce %N<sub>2</sub>O to 0 for five minutes.

**Ultra PC™ Flowmeters:** Turn the unit on and adjust the Total Flow to 5 lpm. Gradually increase TOTAL FLOW and %N<sub>2</sub>O until desired effect is reached. Use the Titration method when adding N<sub>2</sub>O to avoid negative patient reactions. See Note 8 below. At the end of the procedure, reduce %N<sub>2</sub>O to 0, increase Total Flow (oxygen) to 9 lpm.

**Ultra DC™ Flowmeters:** Turn unit on. O<sub>2</sub> flow should start at 3 lpm. Gradually increase N<sub>2</sub>O flow using the Titration method until desired effect is reached. If needed increase O<sub>2</sub> flow. See note 8 below. At the end of the procedure, reduce N<sub>2</sub>O flow to 0 and increase O<sub>2</sub> flow to 9 lpm.

**NOTE:** Axess™ Systems do not operate in the same fashion as standard reservoir bag systems because of the small tubing and lack of bag. A reservoir bag system can provide up to three times the amount of gas flow during inhalation as the bagless system. Because of this, the flows must be higher and room air must often be pulled in through special ports to compensate. Whenever room air is pulled into the gas flow, the N<sub>2</sub>O is diluted. To achieve the same effect on most patients, higher flows and max N<sub>2</sub>O percentages must often be used. As a precaution, the max %N<sub>2</sub>O from the flowmeter is limited to 70% to keep the N<sub>2</sub>O concentration at a safe, reasonable level for patients of all ages. The scavenging vacuum flow is capable of removing these flows so that the staff remains safe.

*Axess™ Circuits can be reprocessed per 26613 or 26613-INT. Do not exceed 134° C for 4 minutes.  
 UDI Product Tag included for product traceability purposes.*

## Indications for use:

To be used in nitrous oxide/oxygen sedation systems for delivering to a patient a mixture of nitrous oxide and oxygen gases with a maximum nitrous oxide concentration of 70%.

Axess™ Nasal Masks are single use only.

## Contraindications:

Contraindications for use of nitrous oxide/oxygen inhalation may include patients with:

- Some chronic obstructive pulmonary diseases
- Severe emotional disturbances or drug-related dependencies
- First trimester of pregnancy
- Treatment with bleomycin sulfate
- Methylenetetrahydrofolate reductase deficiency
- Severe asthma

Whenever possible, appropriate medical specialists should be consulted before administering analgesic/anxiolytic agents to patients with significant underlying medical conditions (e.g., severe obstructive pulmonary disease, congestive heart failure, sickle cell disease, acute otitis media, recent tympanic membrane graft, acute severe head injury)! Operator must review patient history with regard to these health issues.

<sup>1</sup> American Academy of Pediatric Dentistry. Clinical Guidelines: Guidelines on Use of Nitrous Oxide for Pediatric Dental Patients. *Pediatr Dent* 2009;31(6):148-151. Available at: "http://www.aapd.org/media/policies.asp". Accessed February 23, 2010.

Do not use this device for the administration of general anesthesia or as part of, or in conjunction with, a general anesthesia administration system.

## Warnings:

To be used only by a licensed and trained professional in the use of nitrous oxide, using titration method. Patient should always be closely monitored during nitrous oxide use. If patient has an adverse reaction, reduce or stop the flow of nitrous oxide as needed. The O<sub>2</sub> flush button can be used to rapidly purge the lines of N<sub>2</sub>O. If patient does not show signs of quick recovery, remove nasal hood and treat with pure oxygen from either the O<sub>2</sub> resuscitator fitting or an auxiliary oxygen tank using a demand valve, oxygen assisted manual resuscitator, or equivalent. Call for emergency assistance if rapid response is not achieved.

## Cautions:

Federal (U.S.) law restricts this device to sale by or on order of a dentist or physician.


Accutron, Inc.  
1625 W. Pinnacle Peak Rd., Phoenix, AZ 85027 USA  
800.531.2221 | HuFriedyGroup.com/Accutron

[EU REP] MT Promedt Consulting GmbH  
Ernst-Heckel-Straße 7  
66386 St. Ingbert  
Germany

[CH REP] MDSS CH GmbH  
Laurenzenvorstadt 61  
5000 Aarau  
Switzerland

[UK REP] MDSS-UK RP Ltd.  
6 Wilmslow Road, Rusholme  
Manchester M14 5TP  
United Kingdom

CE 0482

No expiration date required.  
Not made with natural rubber latex.

Made in USA  
28904 / REV 5 - 06/2025

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.  
©2025 Hu-Friedy Mfg. Co., LLC. All rights reserved.

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EU REP | MT Promed Consulting GmbH
CH REP | MDSS-CH GmbH
UK REP | MDSS-UK RP Ltd.
©2025 Hu-Friedy Mfg. Co., LLC. Toate drepturile rezervate.