
# VistaCool™

EN

Direct-to-Drain Cooling System for Autoclave Wastewater  
Installation Guide & Owner's Manual

**Model CVV7601**

**QAI**<sup>®</sup>  
C US  
Includes QAI-Certified  
VistaCheck™ Dual Check  
Backflow Preventer

**PATENT PENDING**

**VistaCool model CVV7601  
can be used with:**

SciCan STAT/M® 2000 Autoclave  
SciCan STAT/M 5000 Autoclave


**HuFriedyGroup**  
The Best In Practice

# Introduction


The VistaCool™ Direct-to-Drain Cooling System for Autoclave Wastewater is a patent-pending system designed to eliminate the need for autoclave condenser/ waste bottles and condensation tanks by cooling the exhaust water and automatically sending it directly down the drain.

## IMPORTANT INFORMATION

Please read this entire manual before proceeding with installation and operation, and always follow local plumbing codes. This system has been verified for installation on SciCan STAT/M® 2000 and 5000 Autoclave models only.


**x2** **VS9125**  
Wall-Mount C-Clip (3.5")

**VS7507**  
Drain Line Assembly


**CVVC250-A**  
VistaCheck™ Dual Check  
Backflow Preventer


**VS7508**  
Inline Thermal  
Sensor Assembly


**x4** **VS7539**  
Cabinet Grommet


**x2** **VS6112**  
Swivel Plug-In Elbow


**S1056**  
High-Temperature  
Swivel Plug-In Elbow


**VS6119**  
3/8" x 3/8" x 1/4" Union Tee

**VS6135**  
Angle Stop Adapter


**VS6155**  
Inline Shut-Off Valve


**12'** **VTU902**  
Cold Water Feed Line  
1/4" OD x 1/8" ID  
Blue Transparent Polyurethane

**12'** **VTU944**  
Condensate Drain Line  
1/4" OD x 0.170" ID  
Black Opaque LLDPE


**12'** **VTU942**  
Coolant Overflow Drain Line  
1/4" OD x 0.170" ID  
Blue Opaque LLDPE

**9'** **VTU748**  
Teflon™ High-Temperature  
Tubing Feed Line  
1/4" OD x 0.185" ID

# Typical Installation

For **1**  
Autoclave

**1** Drain Line

**2** VS7507  
Drain Line Adapter

**3** Sink Drain

**4** VS6119 or VS6135  
Cold Water Line Fitting

**5** VS7508  
Inline Thermal Sensor Assembly

**6** VS6155  
Inline Shut-Off Valve

**7** CVVC250-A  
VistaCheck™ Backflow Preventer

**8** VTU942  
Coolant Overflow Drain Line

**9** VTU902  
Cold Water Feed Line




**10** VTU944  
Condensate Drain Line

**11** VS7503  
Thermal Regulator Assembly

**12** VTU748  
High-Temperature Feed Line

## Plumbing

- 1 Turn off water supply. Provide for a 1/4" connection to the cold water supply. This is usually accomplished using a compression fitting or one of the supplied fittings (VS6119 or VS6135) on the cold water line. **Be sure any possible construction debris is flushed completely from the supply line.**
- 2 Cut a 3" to 4" piece of the blue poly tubing and connect the installed fitting (from step 1) to the inlet of the inline water valve provided. Be certain the valve is in the closed position (blue lever 90-degrees to the valve body).
- 3 Cut another 3" to 4" piece of blue poly tubing and connect one end to the outlet end of the inline valve, and the other end to the inlet of the supplied QAI-certified VistaCheck™ Dual Check Valve Backflow Preventer.
- 4 Measure the drain tee adapter assembly for placement into the drain piping. Cut and remove a section (usually 1-1/2") in the appropriate place above the trap as shown. Keep the entire assembly as low as possible while maintaining a minimum of 1" offset between the horizontal drain line and the drain line adapter inputs.



- i** Connect the system to cold water supply ahead of foot switch (if present) for uninterrupted water supply.

## Cabinetry

The VistaCool™ Direct-to-Drain Cooling System for Autoclave Wastewater is typically installed in the base/sink cabinet of sterilization centers or at the sterilization tower if plumbed with separate water and drain connections. The 1/4" O.D. Teflon™ High-Temperature Tubing will need to run from the output port of the STAT/M® Autoclave directly to the system. This exhaust line should run in a progressively descending trajectory, per STAT/M Autoclave requirements.



In many cases, holes will need to be drilled through one or more partitions in order to get the tubing from the autoclave to the VistaCool Direct-to-Drain Cooling System for Autoclave Wastewater. In order to protect the tubing from kinking or damage, grommets have been included. Simply drill a 7/8" pilot hole through the metal or wood partition then push the smooth plastic grommet into the holes. The hole entering the cabinet where the system is installed should be slightly above the top of the system so the tubing has a natural gravity angle to help prevent pools of liquid in the tubing.

-----  
**Install the system  
along dashed path**





**The system may be installed horizontally or vertically. If installed horizontally, the coolant overflow port must be UP in the 12 o'clock position. Do not install the system upside down or at an angle other than shown here as acceptable. Affix both included C-clips securely to the desired mounting surface using appropriate screws to ensure the system does not move during operation.**

Coolant overflow  
always at 12 o'clock


## Connecting the System

- 1** Connect the cold water supply to the system manifold using the blue, semi-transparent (1/4" OD x 1/8" ID) polyurethane tubing. Continue the tubing from the outlet of the VistaCheck™ Backflow Preventer to the top stainless hose barb on top of the brass water valve on the thermal regulator. (Optional: warm the end of the tubing with warm water before pushing it over the barb. Do not expose the tubing to direct, uncontrolled heat, such as a soldering torch, as this could cause material damage that might lead to leaks that may not appear immediately.) Make sure the tubing is pushed down completely over the stem of the barb and is flush with the hex stainless base of the barb.
- 2** Connect the high-temp Teflon™ 1/4" OD tubing from the STAT/M® Autoclave to the milky white KYNAR® Fitting on the end of the system manifold. Make sure the tubing does not kink when moving the autoclave or pulling the tubing through the cabinet partition holes. If necessary, use the high-temperature swivel elbow included with the system (S1056) to prevent the high-temp line connected to the STAT/M Autoclave from kinking.
- 3** Connect the black 1/4" OD LLDPE condensate line tubing to the fitting opposite the milky white KYNAR Fitting, then connect the other end to one of the male adapter fittings on the plumber-provided drain connection. If necessary, use an included swivel elbow (VS6112) to prevent the tubing from kinking.
- 4** Cut the black condensate line (from step 3) as close to the drain adapter assembly as possible and install the in-line thermal sensor. Make certain the flow direction on the in-line sensor is oriented towards the drain.
- 5** Connect one end of the 1/4" OD solid blue LLDPE tubing to the coolant overflow fitting on the system—which must be oriented UP in the 12 o'clock position if the system is installed horizontally—and the other end to the drain connection. If necessary, use an included swivel elbow (VS6112) to prevent the tubing from kinking.

## Starting & Using the System

- 1 Turn on the main water supply and check for leaks.
- 2 Slowly open the inline water supply valve fully and check for leaks.
- 3 Prime the system by inserting a small screwdriver in the viewport hole on the brass stem of the thermal regulator and pushing the piston away from the system body to activate cold water flow. The system is primed when water can be heard flowing through the 1/4" solid blue overflow tubing and running to the drain.
- 4 Ensure the system is installed in an elevated position above the drain connection as described previously to ensure proper flow of the condensate line, as required by STATIM® Autoclave instructions. **The coolant overflow fitting should always be UP in the 12 o'clock position if the system is installed horizontally.**
- 5 Run three consecutive autoclave cycles to ensure everything is working properly.

## Troubleshooting

<p><b>System is not draining</b></p>	<ul style="list-style-type: none"> <li>• Ensure water supply valve is in the ON position.</li> <li>• Check for blockages in the tubing.</li> </ul>
<p><b>Autoclave is displaying a cycle fault</b></p>	<ul style="list-style-type: none"> <li>• Exhaust port line is blocked. Check for kinks in the tubing from the autoclave to the system manifold.</li> <li>• Check that water is running from the tank coolant overflow line to the drain adapter.</li> <li>• Make certain the water supply valve is on. If there is no water to cool the coil in the tank, the inline thermal sensor located in the condensate line will respond by closing the pathway to alert you of the problem and protect the plumbing from steam.</li> </ul>
<p><b>The black condensate line appears to be blocked</b></p>	<ul style="list-style-type: none"> <li>• The inline thermal sensor has been activated, blocking condensate flow to the drain, caused by a lack of coolant water in the tank that overheated the condensate. Overheated condensate will cause the wax motor to expand, closing the outlet. To reset the inline sensor and enable flow to resume, allow the system to rest and cool naturally. The inline sensor's wax-motor piston will retract, opening the outlet port.</li> </ul>
<p><b>Water leaking under counter</b></p>	<ul style="list-style-type: none"> <li>• Check all tubing and fitting connections</li> <li>• Push in on all push-to-connect fittings to ensure the tubing is fully engaged in the fitting "O" ring.</li> </ul>

# System Information

Thank you for choosing the VistaCool™ Direct-to-Drain Cooling System for Autoclave Wastewater for your practice! The following chart is for necessary information for future reference. Please fill it out completely and keep this manual in a convenient place for ready access and reference.

<b>VistaCool™ Direct-to-Drain Cooling System for Autoclave Wastewater</b>		<b>System Model #</b>	<b>CVV7601</b>
<b>Install Date</b>		<b>System Serial #</b>	
<b>Dealer</b>		<b>Installer Name</b>	
<b>Dealer Phone</b>		<b>Installer Phone</b>	
<b>Dealer Address</b>		<b>Notes</b>	

## 1-Year Limited Warranty

Vista Systems are warranted as provided in our Warranty Policy available at <https://www.hufriedygroup.com/documents/62903> or by scanning the QR Code below:

### The following apply notwithstanding the Warranty Policy:

1. Filtration elements subject to varying types of water conditions are warranted for defects in materials and workmanship but not for performance due to fouling by local water conditions.
2. Crosstex disclaims all responsibility and liability for any cost of field labor or other charges incurred by customer in removing and/or re-affixing any Crosstex or third party product, part or component thereof.

3. All warranties are disclaimed for any product that is not installed in accordance with applicable local codes, ordinances and good trade practices, and any product that is moved from its original installation location.
4. Within ten (10) days of installation, please register the product by mailing in the completed warranty card that is provided with the product or registering online at <https://www.hufriedygroup.com/en/services/vista-warranty>.
5. In the absence of suitable proof of installation date, the warranty period will commence on the date of manufacture plus one hundred eighty (180) days.

### Direct all notices, etc. to:

Service Department:  
(800) 483-7433

**CROSSTEX**

Manufactured for: Crosstex International, Inc.  
1625 W. Pinnacle Peak Rd., Phoenix, AZ 85027 USA

Made in USA

VM1065 / REV 5 - 04/2026

CROSSTEX™ is a trademark of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies. All other product and company names referenced are trademarks of their respective owner.

©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved.

**HuFriedyGroup**  
The Best In Practice

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Fabriqué pour : Crosstex International, Inc.
CROSSTEX™ est une marque déposée de Hu-Friedy Mfg. Co., LLC, de ses sociétés affiliées ou apparentées.
©2026 Hu-Friedy Mfg. Co., LLC. Tous droits réservés.