ConFirm™ Rapid 20 Minute

STEAM Biological Indicator Instructions for Use

Intended Use

The ConFirm™ Rapid 20 Minute **STEAM** Biological Indicator is used for routine monitoring, qualification testing, load monitoring and product testing of the following steam sterilization cycles:

- Dynamic Air Removal 270°F (132°C) 4 minutes
- Dynamic Air Removal 275°F (135°C) 3 minutes
- Gravity 250°F (121°C) 30 minutes
- Gravity 270°F (132°C) 15 minutes.

When used in conjunction with the ConFirm™ Rapid 20 Minute **STEAM** Incubator, the Incubator provides a fluorescent result within 20 minutes.

Performance Characteristics

The ConFirm™ Rapid 20 Minute **STEAM** BI has a validated incubation time of 20 minutes. Dynamic Air Removal includes precuum (both standard and IUSS) and SFPF sterilization cycles.

Instructions for Use:

- As applicable, load the chamber with items to be sterilized.
- Before use, examine the BI to ensure that the media-filled cap is not fully sealed on the top of the vial and that the foil seal on the cap remains intact. Without removing the cap, examine the BI for evidence of media in the vial. If media is detected in the vial, the seal has been compromised, do not use the BI.
- Observe the process indicator located on the BI cap (1) and verify that the color has NOT changed from pink to brown, indicating exposure to steam.
- Place the BI in an appropriate representative instrument pack/challenge pack in a manner consistent with items being sterilized. The BIs can also be placed within loads for periodic quality assurance testing per ANSI/AAAMI S779.
- For routine monitoring or for loads with implants:
 - For standard wrapped goods, place the instrument pack/challenge pack in a loaded chamber in the most challenging location for steam to penetrate, typically on the bottom shelf near the door and over the drain.
- Initiate the sterilization cycle.
- Upon cycle completion, review the sterilizer cycle printout. If all parameters were met, remove the instrument pack/challenge pack from the sterilizer and continue to Step 8. If any parameter is not met, remove all items, dispose of the BI and, if applicable, reprocess the entire load with a new BI in an appropriate instrument pack/challenge pack.

CAUTION: Instrument pack/challenge packs and test configurations processed through steam sterilization may still be hot.

- Remove the instrument pack/challenge pack from the Sterilizer.
- Remove the BI from the instrument pack/challenge pack.

NOTE: Per ANSI/AAAMI S779, Items containing implants should be quarantined until the results of the BI are available.

- Examine the BI for evidence of media in the vial.
 - If there is any evidence of media in the vial, the seal has been compromised and the test must be repeated using a fresh BI.
 - If no damage is evident, proceed to Step 11.
 - Observe the process indicator located on the BI cap and verify that the color has changed from pink to brown.
 - If the process indicator has changed from pink to brown, the BI has been exposed to steam.
 - If the process indicator did not change correctly, steam did not come into contact with the BI. Follow departmental procedures for reporting sterilization failures. Do not use items processed in the load. These items must be reprocessed.
- NOTE:** If the BI will not be incubated immediately, do not activate. Store the BI at room temperature for up to 72 hours following completion of the sterilization cycle.
- To activate the BI, twist the cap clockwise while holding the vial to close the BI cap. (2) Transfer the media from the cap to the vial by holding the BI firmly by its cap and flicking the twist of the cap. (3) The BI is properly activated when the cap is sealed tightly on the vial and the media has come in contact with the spores on the bottom of the vial. (4) The media in the bottom of the vial will be purple.
 - Immediately place the activated BI into a well of the ConFirm™ Rapid 20 Minute **STEAM** Incubator. The incubator will display results when incubation is complete.
 - If the BI is positive for the detection of remaining organisms, the conditions for sterilization were not achieved. Follow departmental procedures for reporting sterilization failures.
 - If the BI is negative for the detection of remaining organisms, the conditions for sterilization were achieved. Record the results per facility process.

NOTE: See the ConFirm™ Rapid 20 Minute **STEAM** Incubator operator manual.

- Biological controls should be performed each day a BI is used and whenever a new BI lot is used.
- Repeat Steps 12 – 13 using an unprocessed BI with the same lot number as that used in the sterilization cycle. The unprocessed BI is a positive control used to confirm that the BI contains viable organisms and will be positive for growth.

NOTE: The process indicator will remain the pink start color as it has not been exposed to steam.

- If the unprocessed biological control does not demonstrate a positive fluorescent growth response, all BI tests performed since the previous positive control are invalid. Follow departmental guidelines for reporting sterilization failures.
- Document and retain the final results.

Storage Conditions

Prior to use, the ConFirm™ Rapid 20 Minute **STEAM** BI should be stored at 60-75°F (16-24°C), 30-60% relative humidity (RH). Avoid contact with sterilants or chemicals.

Expiration Date

Do not use beyond the expiration date provided on the BI packaging and labeling. The lot number is the expiration date in the format YYYYMMDD. NOTE: Due to space constraints the lot number/expiration date is expressed in the format YYMMDD on the BI label. Go is used to denote that *Geobacillus stearothermophilus* is present for **STEAM** applications.

Performance Limitations

- Single use only. Do not reuse BIs.
- ConFirm™ Rapid 20 Minute **STEAM** Biological Indicators are part of a quality control system for sterility assurance. They cannot be used as a sole means for validating the sterilization process.
- ConFirm™ Rapid 20 Minute **STEAM** Biological Indicators should not be used by anyone with color blindness specific to the colors used in the throughput indicator. Ensure that users are able to correctly distinguish processed and unprocessed indicator colors.

Do not mark or place any identifying objects (i.e. stickers) on the side of the vial prior to completion of incubation and interpretation of the results.

- All personnel involved must be trained in the proper use of the BI and its corresponding indicator.
- BI is not designed to be exposed to freezing conditions for transport/ storage as the potential for freezing and leaking of media may result.
- If there are any problems interpreting the results from this product, contact HuFriedyGroup.

Safety Related Precautions

The BI contains a dye that is safe under normal conditions of use.

The process indicator on the BI contains an indicator formulation that is safe under normal conditions of use.

Disposal Instructions

The ConFirm™ Rapid 20 Minute **STEAM** Biological Indicator should be treated as standard microbiological waste, non-pathogenic species. Use a steam sterilization cycle at 250°F (121°C) for no less than 30 minutes, incineration, or other suitable means to safely dispose.

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- Biologiske kontroller skal udføres hver dag a BI is used and whenever a new BI lot is used.
- Repeat Steps 12 – 13 using an unprocessed BI with the same lot number as that used in the sterilization cycle. The unprocessed BI is a positive control used to confirm that the BI contains viable organisms and will be positive for growth.
- If the unprocessed biological control does not demonstrate a positive fluorescent growth response, all BI tests performed since the previous positive control are invalid. Follow departmental guidelines for reporting sterilization failures.
- Document and retain the final results.

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