INSTRUCTIONS FOR USE

This diamond-coated PWR Piezo Tip has been designed for use in combination with PWR Pair Combo and PWR Piezo ultrasonic scalers distributed by HuFriedy/Group and in accordance with the intended use.

Using the specific Piezo Tips, the following treatments can be performed:

- Scaling: all procedures for the removal of bacterial plaque deposits and supragingival, subgingival and interdental calculus and the removal of stains;
- Periodontology: periodontal therapy for scaling and root-planing/debridement, including cleaning and irrigation of the periodontal pocket;
- Implant surface cleaning (PWR Pair Combo);
- Endodontics: all treatments for canal preparation, irrigation, filling, gutta-percha condensation, endodontic reprocessing and retrograde cavity preparation;
- Restorative and prosthetics: cavity preparation and removal of carious tissue, removal of prosthetics and excess restoration materials, amalgam condensation, finishing of prosthetic abutment.

The torque wrench supplied with the Piezo Tip, or sold separately, must only be used to secure and remove the Piezo Tips from PWR Pair Combo and PWR Piezo ultrasonic scalers distributed by HuFriedy/Group.

⚠ **WARNING:** The Piezo Tip must be used only by specialised and properly trained personnel, such as Physician, Dentist and/or Dental Hygienist, able-bodied adults of any weight, age, height, gender and nationality. This Piezo Tip must be used in accordance with the correct settings described in the annex "PWR tips setting table".

**NOTE:** Regarding the intended patient group, patient selection criteria, indications for use, contraindications and warnings, refer to the PWR Pair Combo and PWR Piezo devices' "USE & MAINTENANCE MANUAL".

The terms **"WARNING"**, **"CAUTION"** and **"NOTE"** identify conditions that require particular care and must, therefore, be carefully examined. Specifically:

⚠ **WARNING:** Identifies conditions and practices that present a risk of death or serious injury for the patient and/or user.

⚠ **CAUTION:** Identifies conditions and practices that may cause minor injuries and/or damage to the devices.

**NOTE:** A warning not associated with a condition of risk or imminent danger.

⚠ **CAUTION:** The information contained in this document is limited to that strictly necessary for the first use of the product and must be consulted as a supplement to the PWR Pair Combo and PWR Piezo devices' "USE & MAINTENANCE MANUAL". Read the PWR Pair Combo and PWR Piezo's "USE & MAINTENANCE MANUAL" paying particular attention to the "Safety Requirements" chapter, before performing any type of operation on the system and commencing any type of treatment.

⚠ **WARNING:** If a serious accident attributable to the device occurs during correct and intended use, it is recommended to report it to the Competent Authority and to the company indicated on the product label.

⚠ **WARNING: SINGLE USE ONLY**

The diamond-coated PWR Piezo tip is intended to be used on an individual patient for single use only.

- Discard it after each use.
- Do not attempt to clean or sterilize it after use.

⚠ **WARNING: BEFORE FIRST USE**

The Piezo Tip and torque wrench are supplied "non-sterile", therefore, before use, they must be reconditioned in accordance with the procedures summarily described below and reported in detail in the "Cleaning and Sterilisation" section of the PWR Pair Combo and PWR Piezo devices' "USE & MAINTENANCE MANUAL". **After each use**, the torque wrench must be immediately cleaned and sterilised.

1- Cleaning procedures

**IMPORTANT:** Submerge the Piezo Tip or the torque wrench in demineralised water or in an enzymatic cleaning solution immediately after use. Do not leave residue or blood deposits on the Piezo Tips and instruments, eliminate larger impurities with a disposable cloth or paper towel.

The cleaning procedure can be performed using two different methods: manual cleaning and automatic cleaning. These two methods are alternative and only one or the other must be performed.

Repeated reconditioning has a minimal effect on these devices and instruments. The end of service life is generally determined by wear or damage due to use.

⚠ **CAUTION:** Do not use hydrogen peroxide to clean the Piezo Tips.

⚠ **WARNING:** To avoid injury to operators, treat the Piezo Tips separately.

Particular attention must be paid while handling sharp and potentially infected instruments to avoid the risk of wounds and infections. Health personnel in charge of use and cleaning of the devices must use suitable protection and always use care when handling instruments with sharp edges and pointed tips.

1.1- Manual Cleaning

**NOTE:** The tip, and torque wrench must be processed separately.

Do not leave the tip inserted into the torque wrench.

Do not place the tip and the torque wrench in the same container.

Do not use the same toothbrush for the tip and torque wrench.

**NOTE:** Process independently validated with Enzymec enzymatic detergent and 0.8% dilution.

1 Prepare a pH-neutral (6-9) enzymatic detergent solution, following the manufacturer's instructions.

⚠ **CAUTION:** Once used, dispose of the enzymatic cleaner solution properly, do not recycle.

2 Place the tip or torque wrench in a clean container, in a horizontal position. Add an adequate amount of enzymatic solution to fully cover the device to be cleaned.

3 Soak the device in the enzymatic detergent solution for 10 minutes at room temperature. This process reduces the amount of blood, protein, and mucus on the instrument.

4 After 10 minutes of immersion in the enzymatic solution, gently scrub the device internal and external surfaces with a soft-bristled brush until any visible dirt is removed. For each part to be cleaned use a brush that is suitable to its dimension. Use different brush for the tip and for the torque wrench.

**NOTE:**  
**For the tip:** Thoroughly clean hard-to-clean areas such as sharp edges and especially gaps between cutting tips.  
**For the torque wrench:** Brush thoroughly, for about 20 seconds, all the following parts:

- through holes and internal channels;
- external metal ring;
- internal cavities, grooves and fissures.

5 Remove the device being cleaned from the enzymatic detergent solution.

6 Thoroughly rinse and brush all the surfaces of the device (see point 4) under running water for:

- at least 10 minutes for the wrench;

• at least 1 minute for the tip.

7 **For the tip**, add the following steps: use a 20 ml disposable syringe to flush the enzymatic detergent solution into hard-to-reach areas (such as inner channel, through-holes or / canulais). Repeat this step three times to ensure the effective removal of dirt from the internal surfaces. Each time use freshly prepared solution.

8 Place the device to be cleaned in a stainless-steel instrument tray with a mesh bottom. Place the tray in the ultrasonic cleaning tub and make sure that in the bath, at 75.2°F ±3.6°F (24°C ±2°C), the instruments are submerged in fresh enzymatic solution (see step 1).

9 Ultrasonicate the torque wrench for at least 20 minutes and the tip for 10 minutes, or as instructed by the manufacturer of the enzymatic cleaner and ultrasonic tank.

10 Remove the device from the ultrasonic cleaning tank and repeat step 6.  
**For the tip:** flush three times the inner channel of the device to be cleaned with injected demineralized water using a 20 ml syringe in order to remove any dirt residues.  
**Soak the device** to be cleaned in demineralized water for 1 minute.  
**For the torque wrench:** Soak the device in demineralized water for 5 minutes.

1.2- Automatic Cleaning

**NOTE:** The Piezo tip, and the torque wrench must be processed separately. Do not leave the tip inserted into the torque wrench.

Automatic cleaning involves the use of a thermodisinfecteur and the following materials:

- Alkaline detergent. For example, neodisher® FA (0.2 % v/v);
- Neutralizing liquid. For example, neodisher® Z (0.1 % v/v);
- Metal basket;
- Adaptors;

**NOTE:** Make sure that the accessories are appropriately blocked in the basket and cannot move during washing. Any blows could damage them. Position the instruments in a way that the water can flow through all the surfaces, even internal.

⚠ **WARNING:** Avoid overloading the thermal disinfecteur as this could compromise cleaning effectiveness.

- 1 Place the accessories in a metal basket using the special adaptors (supplied as optional).
- 2 Set the sequence and parameters for the cycle in the thermal disinfecteur:
  - 1 min, Rinse with cold water;
  - 5 min, Wash with alkaline detergent at 131.00°F ±3.6°F (55°C ±2°C)
  - 1 min, Neutralisation with neutralising liquid and deionized water at 89.6°F ±3.6°F (32°C ±2°C);
  - 1 min, Rinse with deionized water at 89.6°F ±3.6°F (32°C ±2°C);
  - 5 min, Thermal disinfection at 199.4°F (93°C) with deionized water.

Thermal disinfection has not been tested experimentally. In compliance with ISO 15883-1, Table B.1 (4) thermal disinfection at a temperature of 194°F (90°C) for 5 min determines a value of A0 3000.

⚠ **CAUTION: After cleaning**

- The use of alkaline detergents removes the lubricant from the torque wrench, causing malfunctions and increasing its wear. Therefore, it is important to lubricate the torque wrench after each cleaning procedure using medical-grade lubricants (see chapter "2- Lubrication"). The dilution of the alkaline detergent required for automatic cleaning is such that it does not damage the silicone parts of the torque wrench (o-ring).
- Use filtered compressed air to flush out any moisture from cavities, grooves, crevices, and other hard-to-reach areas. This prevents the appearance of stains or streaks on the surface and oxidation. Make sure that the device is completely dry even internally before starting the sterilization cycle.
- Check that the device is actually cleaned and repeat the process if the water does not run clear in the final rinse, or if dirt residues are still visible.
- All parts must be inspected before their sterilisation. In general, it is sufficient to conduct a visual inspection in good lighting conditions without the use of magnifying means. All parts must be checked, looking for visible dirt, damage and/or corrosion, signs of wear. Particular attention must be paid to:
  - through-holes,
  - threading
  - internal channels,
  - the edges of the sharp tips
  - edges that may be deformed or largely affected areas: the edges must be continuous.

If surfaces are not visually clean repeat cleaning and inspect again. Discard damaged tools.

2- Lubrication

Before sterilisation, the torque wrench must be lubricated with a commercial medical-grade lubricant.

⚠ **CAUTION:** Do not lubricate the Piezo Tip.

- The lubricant must be applied by spraying it directly onto the peripheral contact surface inside the torque wrench, as indicated in Fig. 2;
- after having applied the lubricant, remove any excess of lubricant using a clean lint-free cloth.

3- Packaging

The Piezo Tips, and the torque wrench must be sterilised using a standard, appropriately sized medical steam sterilisation pouch.

⚠ **CAUTION:** The Piezo Tips, and the torque wrench must be packaged separately. Do not package multiple Piezo Tips and/or torque wrenches in the same pouch.

⚠ **CAUTION:** Make sure the pouch is big enough to hold the instrument without tensioning the seals and without ripping the packaging.

4- Sterilisation

Sterilise only by pre-vacuum steam autoclaving.

The steam sterilisation parameters have been validated by the manufacturer to guarantee a sterility level (Sterility Assurance Level - SAL) of 10<sup>-6</sup>.

⚠ **CAUTION:** If multiple instruments need to be sterilised in the autoclave in a single cycle, do not exceed the maximum allowed load.

⚠ **CAUTION:** Risk of contamination. Do not use gravity autoclaves to sterilise the Piezo Tips. The operating cycle of gravity autoclaves does not guarantee the suitable sterilisation of the internal channel, of the cavities and hard-to-reach points.

⚠ **CAUTION:** Use only a pre-vacuum autoclave to sterilise the Piezo Tips and the torque wrench. Do not use other methods of sterilisation insofar as they may be incompatible with the materials used to produce the instruments.

Do not use the following methods of sterilisation: ethylene oxide sterilisation, hot air sterilisation, flash autoclaving, STERRAD sterilisation, sterilisation with STERIS or similar sterilisation systems. Do not use the following substances/systems to sterilise the instruments: hydrogen peroxide, peracetic acid, Formaldehyde, Gluteraldehyde based systems or other equivalent solutions/systems.

Minimum sterilisation parameters to guarantee a sterility level (Sterility Assurance Level - SAL) of 10<sup>-6</sup>:

⚠ **WARNING:** these sterilisation parameters, validated by an independent laboratory, refer exclusively to a PRE-VACUUM autoclave operating cycle.

- Type of cycle: PRE-VACUUM - (3 pre-vacuum stages).
- Minimum temperature: 270°F (132°C) - (interval between 270°F + 275°F (132°C + 135°C)).
- Minimum exposure time (in package): 4 minutes for the tip and for the torque wrench.
- Minimum drying time: 10 minutes.

⚠ **CAUTION:** Do not exceed a temperature of 275°F (135°C).

All stages of sterilisation must be performed by the operator in compliance with the current revision standards: UNI EN ISO 17665-1, UNI EN ISO 556-1 and ANSI/AAMI ST46.

5- Important information relative to safety of Piezo Tips

⚠ **WARNING: Breakage and wear of tips.**

On rare occasions, the oscillations at high frequency and wear may cause a Piezo Tip to break. Do not bend, reshape or sharpen the Piezo Tips in any way. The alteration of a Piezo Tip may cause it to break. Altered Piezo Tips must never be used. If a Piezo Tip breaks during use, recover and remove all fragments of the tip from the patient, with maximum care. During the treatment, ask the patient to breathe through their nose and/or use a dental dam to avoid the risk of swallowing any fragments generated by the breakage of the tip.

⚠ **WARNING:** The Piezo Tip is consumable. Thoroughly inspect the Piezo Tip before and during each treatment for damage and/or excessive wear (see table "Piezo TipS WEARING OUT EXAMPLES"). Do not use a Piezo Tip if it is scratched or rusty. The Piezo Tip may break during use. If you notice any damage or loss in performance, replace the Piezo Tip with a new one.

**NOTE:** A deformed Piezo Tip may preclude the complete functionality of the device used, and, depending on the model, a message may be displayed.

- Before treatment, make sure that the chosen Piezo Tip is firmly attached to the handpiece. Properly tighten the Piezo Tip on the handpiece using the torque wrench supplied with the device.
- Do not change the shape of the Piezo Tip in any way. Bending or forcing the Piezo Tip may cause it to fracture. Never use bent Piezo Tip.
- Do not attempt to sharpen used Piezo Tips as they may break during use.
- Use only original PWR Piezo Tips. The use of non-original tips will void the warranty of the device. Never use non original tips and never use tips with incompatible devices as they will damage the device and may cause injury to operators or the patient. The use of non original tips will damage the threaded pin of the handpiece, thus compromising the secure attachment of tips, even if original, to the handpiece. In addition, the device settings are tested and guaranteed for correct operation using only original PWR Piezo Tips.
- It is advisable to avoid the application of excessive force or prolonged contact of the Piezo Tip on the soft tissues to avoid thermal damage and / or injury.
- Let the ultrasonic vibrations work, do not exert excessive pressure on the Piezo Tips during use. Apply a light force on the Piezo Tip to obtain the best efficiency.
- Applying excessive pressure could cause the Piezo Tip to fracture, which in turn can lead to injury.
- The Piezo Tip must always be kept moving. If the Piezo Tip is blocked, it can cause overheating of the treated part. It is recommended to use continuous movement to minimize contact between the tip and the part. Do not block it against the tissue so as not to cause it to overheat. It is advisable to use high levels of irrigation as the power level increases.
- The Piezo Tips vibrate with a longitudinal oscillation, with forward and backward movement. During treatment, always keep the instrument tangential to the tooth surface. Move the handpiece back and forth while applying light lateral pressure. Do not apply pressure to the tip in the axial direction
- Do not aim the instrument directly on the surface of the enamel or implant. Position the tip/ operative part only tangentially to the surface of the tooth or implant.
- When the Piezo Tip is used in the interproximal spaces, do not block the instrument or leverage the operative part. The Piezo Tips must be left free to vibrate.
- Check the threaded parts of the Piezo Tip and those on the handpiece. These parts must be thoroughly cleaned.
- Check the state of wear of the tip and its integrity before and during each use.
- If you experience a loss of performance provide for its replacement.
- Use of a worn tip shall lower performance in terms of cutting power.

Treat the following objects as hospital waste:

- Piezo tips, when worn or broken;
- Torque wrench, when worn or broken.

6- Use of the torque wrench.

Refer to Fig. 1.

⚠ **WARNING:** During operations to secure and remove the Piezo Tip, the user must pay special attention to the sharp and cutting edges of the Piezo Tips.

6.1- Securing the Piezo Tip

For correct use of the HuFriedy/Group torque wrench, operate as follows:

- 1 Screw the pre-selected Piezo Tip all the way down onto the handpiece;
- 2 Insert the Piezo Tip into the torque wrench, as illustrated in Figure 1;
- 3 Firmly grip the body of the handpiece.

⚠ **CAUTION:** Do not grip the handpiece by the end and/or cord, but only by the body (Figure 1 - Ref. C). Do not rotate the body while securing the Piezo Tip. Firmly grip the body of the handpiece and turn the torque wrench only.

4 Turn the torque wrench in a clockwise direction until the friction clicks (the external body of the torque wrench turns with respect to the body of the handpiece, emitting mechanical "CLICK" sounds).

5 The Piezo Tip is now properly locked in.

6.2- Removing the Piezo Tip

Position the torque wrench on the Piezo Tip and screw in an anti-clockwise direction.

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All company and product names are trademarks of Hu-Friedy Mfg Co., LLC. Its affiliates or related companies, unless otherwise noted.
©2024 Hu-Friedy Mfg Co., LLC. All rights reserved
Hu-Friedy Mfg Co., LLC