Piezo Tips

INSTRUCTIONS FOR USE

EN

Intended Purpose: Piezo Tips

The Piezo Tips are the working-end devices that remove calculus deposits from teeth and debris from gum lines during ultrasonic dental cleaning and periodontal therapy.

Indications:

The Piezo Tips generate ultrasonic vibrations intended for ultrasonic procedures:

1. General Supra and Subgingival scaling applications
2. Periodontal debridement for non-surgical treatment of Periodontal diseases.

Contraindications:

The Piezo Tips should not be used in patients with cardiac pacemakers. Devices which show evidence of corrosion should not be used as they may pose an elevated hypersensitivity risk to patients with nickel-allergy. Do not expose to phenols or iodophors. Do not use dry heat sterilization, or heat above 135°C (275°F).

Intended Patient Population:

The intended patient population is adult and pediatric patients under the care of a dental professional.

General Considerations:

First use instructions:

- After receiving a new Piezo Tip, make sure to follow initial cleaning, disinfection and sterilization steps before using for the first time. All non-sterile instruments are to be cleaned, disinfected, and sterilized prior to first use and each subsequent use.

Sterilizer and automatic cleaning equipment best practices:

- Sterilization and cleaning equipment should be maintained and checked per the manufacturer's recommendation.
- Consider the legal provisions valid for your country as well as any hygienic instructions of the doctor's practice or hospital.

Water quality

- The use of completely deionized or distilled water can help avoid most water quality problems.

Installation

 **Only use the Guardian to insert the Piezo Tips.**

- Your Piezo tip comes with its own tip carrier/torque limiting wrench (the Guardian).
- Attach tip to handpiece using the Guardian by aligning the threaded end of the tip to the handpiece.
- Rotate the Guardian clockwise.
- Complete installation of the tip by rotating the Guardian until you hear clicking.
- The Guardian must be used to completely tighten the tip during installation. Incomplete tightening may cause the tip to become loose or lack vibration.
- When using the Guardian without a tip previously inserted into the tip carrier, take the following step: place the Guardian over the tip by slipping tip through the small hole in the bottom, closed end of the tool. This will position the tip within the Guardian for safety and convenience.

Tip Removal

- Hold the handpiece vertically. Insert the tip through the hole in the bottom of the Guardian.
- Depress wrench slightly while rotating to fully secure tip to wrench.
- Turn the Guardian counter-clockwise to loosen and remove tip from handpiece.
- Once removed from handpiece, do not remove the tip from the Guardian carrier.

Helpful Advice

- Color of the Guardian (torque limiting wrench/tip carrier) indicates equipment application mode (if applicable). Refer to your equipment instructions for use for Scaling and Perio application mode.
- Pre-select the appropriate power setting prior to patient treatment. See recommended power setting guide. The use of higher than recommended power setting may cause excessive wear or tip breakage.
- Always remember to properly tighten the tip to handpiece using the Guardian wrench.
- Only apply the lateral edge of the tip to the tooth. The point may be applied only when using rounded tips.
- Apply tip at a 0-15° angle to tooth surface.
- Use light pressure and keep the tip moving at all times when in contact with the tooth surface.
- Activate the tip prior to placement on the tooth to lubricate the tip for gentle placement subgingivally.
- Blade-shaped tips can not be sharpened.
- Scaling efficiency can significantly diminish with worn tips. We recommend checking tips for signs of wear regularly.

Warnings:

The Piezo Tips (i) should be replaced when "2mm wear occurs"; (ii) finger stick/cut, scrapped finger on insert, cut hand/finger, or needle stick/puncture may occur with potential for the transfer of bloodborne pathogens during tip removal from handpiece if care is not exercised; (iii) 'material fragmentation' (e.g. broken tip) may occur if used improperly; (iv) accidental swallowing of inserts/tips/tip fragments may occur; (v) device fracture/break may occur if used improperly;

(vi) foreign body in patient may occur (e.g. tip swallowed); (vii) device embedded in tissue or plaque may occur; (viii) perforation of tooth may occur; (ix) patient discomfort (e.g. pain, burn or bleeding) may occur; (x) exposure to blood/bodily fluids; (xi) Patient infection/Sepsis may occur; (xii) Unintentional electrical shock from the generator (unit) may occur; (xiii) leaking may occur; and (xiv) loose housing may occur.

Precautions:

- Symmetry S-Series Piezoelectric tips are compatible with Symmetry IQ™, Satelec and NSK Ultrasonic dental generators only.
- Symmetry E-Series Piezoelectric tips are compatible with EMS devices.
- Handle and dispose of sharps in appropriate waste receptacle per facility protocols.
- For intended use only by dental professionals.
- Wash hands before and after treating each patient, before putting on gloves, and immediately after removing gloves.
- Use aerosol mitigation protocols (pre-treatment with antibacterial rinses, high-volume evacuator, etc.) and proper PPE (surgical mask covering nose and mouth, gowns/jackets in work area, medical grade single-use gloves, etc.) to protect from infection and burns.
- Inspect, clean and sterilize before each use. Inserts that have broken, bent, or worn tips should be removed from service before and during use.
- Verify tip connection is secure to the handpiece prior to procedure.
- Do not apply excessive force during installation.

Cleaning and Sterilization

- Piezo tips remain in the Guardian Tip carrier during the complete cleaning, disinfection and sterilization reprocessing cycle.
- Clean using ultrasonic cleaner or thermodisinfector.
- Steam sterilize for at least 5 minutes at 270°F/132°C. Recommend 30 minute dry time after sterilization cycle. Do not heat above 275°F/135°C.
- Clean, disinfect, inspect and steam sterilize before each use following the respective equipment manufacturers users instructions.
- Do not use phenols or iodophors.

The fundamental suitability of the instruments for effective automated cleaning and disinfection was demonstrated by an independent accredited test laboratory under the following conditions:

Washer-Disinfector	Miele Professional G 7836 CD
Racks	Mobile injector unit (Miele) E429, Four-level rack (Miele) E 493
Cleaning Cycle	2 minute pre-cleaning with cold tap water, Draining 5 minute cleaning with 55 °C cleaning solution, Draining 3 minute rinsing with cold deionized water, Draining 2 minute rinsing with cold deionized water, Draining
Cleaning Solution	0.5 % cleaning solution Neodisher® Mediclean Dental (Chemische Fabrik Dr. Weigert, Hamburg)
Validation Report	Project Number: 00418-1 23919 Examination of an Automated Cleaning Process using quantitative Detection of Protein and Hemoglobin and the Radionuclide Method

The fundamental suitability of the instruments for an effective sterilization was demonstrated by an independent accredited test laboratory under the following conditions:

Sterilization Method	Pre-vacuum Mode
Sterilizer	W & H Lisa MB 17 Steam Sterilizer
Sterilization Temperature	134 °C (273°F)
Pre-Vacuum Phases	3
Holding (full cycle)	4 minutes
Drying Time	30 minutes
Validation Report	Project Numbers: 25517-1; 25517-2; 24219 Validation of a Sterilization Process using Steam Sterilization in Pre-vacuum Mode Method MD 4.0: Sterilization validation of medical products with moist heat Project Numbers: 10918-1; 10918-2 Determination of Residual Moisture after Sterilization using Steam Sterilization in Pre-vacuum Mode

Compatible Piezo Generators (Units):

Symmetry S-Series Piezoelectric tips are compatible with Symmetry IQ, Satelec and NSK Ultrasonic dental generators only.

TIP	USE	APPLICATION				SATELEC							NSK
		SYMMETRY S-SERIES MODE	SATELEC MODE	NSK MODE	POWER	INTEGRATED MODULES			TABLE TOP			INTEGRATED/ TABLE TOP	
						SP3055 Generation	SP4055 Generation	SP NEUTRON Generation	P-MAX/ P-MAX Lux PROPHY MAX/ PRO- PHY-MAX Lux PROPHY-MAX S RANSC.MAX/ TANSC.MAX Lux	P-5 Booster	P-Max II PU-Max II		
USS10	Remove heavy su- pragingival deposits	SCALING	SCALING	SCALING	MED	1-5	1-5	8-10	1 - 5	7 - 9	8 - 10	Min - 3	
USP10	remove light deposits and plaque	SCALING	SCALING	SCALING	MED	1-5	1-5	8-10	1 - 5	7 - 8	8 - 10	Min - 3	
US100	Remove light deposits in shallow to moderate pockets	SCALING	SCALING	SCALING	MED	1-4	1-4	6-8	3 - 7	4 - 6	6 - 8	Min - 1	
USU	Remove light deposits and plaque	PERIO	PERIO	PERIO	LOW	1-6	1-6	1-6	1 - 7	7 - 8	1 - 6	Min - 4	
US1D	removal of residual debris and designed to enhance access to deep periodontal pockets primarily in anterior areas	PERIO	PERIO	PERIO	LOW	1-3	1-3	1-4	1 - 10	1 - 3	1 - 4	Min - 4	
US2LD	removal of residual debris and designed to enhance access to deep periodontal pockets and furcations	PERIO	PERIO	PERIO	LOW	1-3	1-3	1-4	1 - 10	1 - 3	1 - 4	Min - 4	
US2RD	removal of residual debris and designed to enhance access to deep periodontal pockets and furcations	PERIO	PERIO	PERIO	LOW	1-3	1-3	1-4	1 - 10	1 - 3	1 - 4	Min - 4	
USP3	remove heavy deposits and necrotic cemen- tum on anteriors and premolars	PERIO	PERIO	PERIO	LOW	1-3	1-3	1-4	1 - 10	1 - 3	1 - 4	Min - Max	
US4L	remove calculus on all surfaces of premolars and molars, especially difficult to reach inter- proximal surfaces	PERIO	PERIO	PERIO	LOW	1-3	1-3	1-4	1 - 10	1 - 3	1 - 4	Min - 4	
US4R	remove calculus on all surfaces of premolars and molars, especially difficult to reach inter- proximal surfaces	PERIO	PERIO	PERIO	LOW	1-3	1-3	1-4	1 - 10	1 - 3	1 - 4	Min - 4	

TIP	USE	APPLICATION		EMS
		SYMMETRY E-SERIES MODE	POWER	Air Flow Piezon Master
UES10	Remove heavy supragingival deposits	SCALING	MED	(1-6)
UEP	remove light deposits and plaque	SCALING	MED	(1-6)
UE100	Remove light deposits in shallow to moderate pockets	PERIO	LOW	(1-3)
E100M	Remove light deposits and plaque	PERIO	LOW	(1-3)
UE1D	removal of residual debris and designed to enhance access to deep periodontal pockets primarily in anterior areas	PERIO	LOW	(1-3)
UE2LD	removal of residual debris and designed to enhance access to deep periodontal pockets and furcations	PERIO	LOW	(1-3)
UE2RD	removal of residual debris and designed to enhance access to deep periodontal pockets and furcations	PERIO	LOW	(1-3)
UEP3	remove heavy deposits and necrotic cementum on anteriors and premolars	PERIO	LOW	(1-3)
UEAFL	remove calculus on all surfaces of premolars and molars, especially difficult to reach interproximal surfaces	PERIO	LOW	(1-3)
UEAFR	remove calculus on all surfaces of premolars and molars, especially difficult to reach interproximal surfaces	PERIO	LOW	(1-3)
UEAFS	remove calculus on all surfaces of premolars and molars, especially difficult to reach interproximal surfaces	PERIO	LOW	(1-3)

In case of a serious incident that has occurred in relation to Hu-Friedy instruments within the European Union, please report by either calling Hu-Friedy for assistance Tel 00800 4837 4339 or send an e-mail to info@hu-friedy.eu. In addition, please inform your national competent authority.

Made in USA

Tips manufactured by and
wrenches manufactured for:

1-800-Hu-Friedy | HuFriedyGroup.com

Hu-Friedy Europe LLC & Co. KG
Bogenstrasse 6 - 8
78576 Emmingen - Liptingen
Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland


MDSS-UK RP Ltd.
Parkway House, Palatine Rd,
Northenden, Wythenshawe,
Manchester M22 4DB
United Kingdom



All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
Satelec and NSK are trademarks of Satelec and Nakanishi Inc., respectively.
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved.

K2044 / REV 9 - 2601


Piezo Tip Wear Guide Instructions

1. Locate the diagram of the tip to be measured.
2. Position point A to point B as shown and/or align the curvature of the tip to match the guide.
3. Rotate until tip is flat against the card.
4. Evaluate tip condition, comparing length to red/blue bars.

Try our probe sized tips which are designed to enhance access into periodontal pockets and improve adaptation in interproximal areas, allowing you to more thoroughly debride difficult to reach areas to achieve better clinical results.

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EU REP Hu-Friedy Europe LLC & Co. KG
CH REP MDSS CH GmbH
UK REP MDSS-UK RP Ltd.
©2026 Hu-Friedy Mfg. Co., LLC. Alla rättigheter förbehållna.