Nº16 - DICIEMBRE 2017

EL EXPLORADOR

DE HU-FRIEDY

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Orthodontic Practice
Clear Collection instruments for clear aligner treatments

In part 1 of a series, Dr. S. Jay Bowman explores instruments that help increase the utility of aligners and expand the scope of appropriate applications.

Let's be clear: enhancing aligner treatments
As clear aligner treatments have evolved and been progressively refined in the past 10 years, more advanced applications and increasingly complex malocclusions have been addressed.^{1,2} In addition, there has been more interest in tracking some of the specific limitations of moving teeth with plastic.^{3,4} Clearly, the progress of orthodontic treatment with traditional wires and braces has been advanced throughout the past century, but the techniques associated with a sequence of aligners are relatively recent phenomena.

Orthodontists with the most experience and attention in using aligners have worked diligently to improve the options through better understanding of the limitations involved in non-scheduled treatment plans, especially when using indirect to enhance the associated effectiveness.^{5,6,7} The Clear Collection instruments from Hu-Friedy (Chicago, Illinois), all of which have been developed with the aim in mind, namely, to help to increase the utility of aligners and expand the scope of applicable applications (Figure 1).⁸

Adding infra- and intermediate forces
The Clear Collection instruments have been designed and standardized to provide the orthodontist with already been providing, while opening up more options for other methods to accomplish or refine treatment in the past, attempting to add intermediate forces for correction of Class II, III, or midline was problematic.

"Innovation is taking two things that already exist and putting them together in a new way."
— Tom Freston

Some were cutting teeth using conventional aligners to achieve better fit for the aligners.

The Tear Drop
There are occasions where it is necessary to add elastic buttons or different localized force on either the buccal or lingual side of the tooth to correct a space of a few millimeters that cannot be developed and that only require a system of use of clear aligners. The Tear Drop instrument, created for adding a tension or force on clear aligners, not only cut in a narrow, tapered "hole" also in place on the aligner.

The Hole Punch
The Hole Punch instrument, created for adding a tension or force on clear aligners, not only cut in a narrow, tapered "hole" also in place on the aligner.

Clear Collection instruments for clear aligner treatments: part 2

In part 2 of a series, Dr. S. Jay Bowman continues his look at instruments that help increase the utility of aligners and expand the scope of appropriate applications.

Contact points to accent aligners
The two second parts in the Clear Collection (i.e., The Horizontal and The Vertical) are designed to enhance desired tooth movements by controlling "contact" points. Although orthodontics is a contact-based that is unique to aligner treatment planning, there are occasions when the virtual aligner does not produce the desired result in fact, sometimes there reported that a percentage of tooth movement prescribed for a setup is simply not translated from objects into the dental results.^{9,10} The flexibility of plastic, the contact areas transmitted from the buccal side (UP) or lingual "expressions" and contact of models, impressions in the vacuum process of fabricating aligners, and the fact that all teeth do not move to the same degree when exposed to forces, created by the force can all lead to incomplete correction.

To improve the predictability of aligned tooth movement, The Vertical and The Horizontal aligners were designed to produce

"We will control the horizontal. We will control the vertical. You are about to participate in a great adventure. You are about to experience the awe and mystery which reaches from the inner mind to —"
— "The Outer Limits" TV Series (1963)

Figure 1: An instrument designed specifically for enhancing the correction of misaligned teeth with clear aligners or even during minor tooth movement when using single clear retainers. Pushing upper aligner and supports to often perforations—especially since aligners do not have a large surface area contact on buccal. Consequently, these incisors can get left behind, resulting

in another form of "lag" or loss of tracking. The Vertical is used to produce an indentation at the mesial or distal of a specific tooth in the buccal and/or lingual aspects of the aligner plastic. These indentations are made without heating the plastic and at a very shallow depth to avoid compression of the integrity of the plastic (Figure 1). The start is to add contact points to accent

Figure 2: The Clear Collection instruments are designed to provide the orthodontist with already been providing, while opening up more options for other methods to accomplish or refine treatment in the past, attempting to add intermediate forces for correction of Class II, III, or midline was problematic.

Figure 3: The Tear Drop instrument, created for adding a tension or force on clear aligners, not only cut in a narrow, tapered "hole" also in place on the aligner.

Figure 4: The Hole Punch instrument, created for adding a tension or force on clear aligners, not only cut in a narrow, tapered "hole" also in place on the aligner.

Bioscience mechanics
If certain teeth are not "tracking" or are lagging" behind (i.e., not fitting into the tray; Figure 7) Aligner ChewiesSM are employed (ChewiesSM Aligner Tray Seaters, Dentistry Resource, Inc., York, Pennsylvania). Patients are asked to hold the Chewie between the teeth in question and squeeze 10-15 seconds, release, and repeat 5-6 minutes, 2-3 times per day.¹¹ (Figure 8). In addition, Chewies are routinely used for the first few days when patients switch to a new pair of aligners to help them seat more completely.

Another option is to add elastics to forcibly erupt the teeth into the aligner tray.^{12,13} There are several methods to generate so-called "boosting" mechanics" (Figure 9).¹⁴

1. An orthodontic elastic is stretched over the plastic aligner to connect to a combination of bonded buttons on both lingual and buccal of the tooth.
2. A bonded button is placed on the lingual of the tooth (plastic closed with the Hole Punch). The Tear Drop notches are cut into the tray

at the gingival embrasure spaces in the buccal plastic on both the mesial and distal of the problematic tooth, and an elastic hooked from the buccal notches and stretched over the aligner tray to the button to produce an eruptive force.

The Tear Drop and/or Hole Punch cuts are made in each tray in the series

of aligners prior to their being delivered to the patient. This concept is also useful when the extrusion of a tooth is needed for restorative purpose.

Other biomechanics that benefit from the application of elastic forces to aligner trays include the correction of severe posterior or anterior intrusion mechanics for open bites,^{15,16} and forces to control

Figure 5: The Clear Collection instruments are designed to provide the orthodontist with already been providing, while opening up more options for other methods to accomplish or refine treatment in the past, attempting to add intermediate forces for correction of Class II, III, or midline was problematic.

Figure 6: The Tear Drop instrument, created for adding a tension or force on clear aligners, not only cut in a narrow, tapered "hole" also in place on the aligner.

Figure 7: "Aligner" or "lag" tracking is most often characterized as an "air gap" between the buccal or lingual of tooth and the plastic, indicating teeth are not following the prescribed tooth movement.

Figure 8: Aligner Chewies are held tightly between specific "lagging" teeth for 10-15 seconds. This process is repeated for 5 minutes, 2-3 times daily, especially when changing to a new pair of aligner trays.

Figure 9: Buccal mechanics to forcibly erupt a "lagging" tooth using orthodontic elastics (SA 08). The Hole Punch is employed to clear aligner plastic to permit the addition of bonded buttons. The Tear Drop is used to cut narrow mesial and distal embrasures (SA 09).

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