HU-FRIEDY PEDO & PERMANENT CROWNS

STAINLESS STEEL PRIMARY AND PERMANENT MOLAR CROWNS

Instructions for Use

Hu-Friedy's Stainless Steel Crowns are intended to provide full coverage, stronger, provisional treatment for primary molars and fill the need for long term provisional coverage on permanent molars. To help minimize contouring and trimming, these stainless steel primary and permanent molar crowns are precrimped, prebelled and prefestooned.

Crown Selection and Preparation

1. Select the proper crown size. Measure the mesio-distal dimension of the tooth in the patient's mouth. Next, measure the crowns in the kit to determine the corresponding size.
2. Once the appropriate size has been identified, test for accurate fit in the mouth. Note the occlusal relationship between the opposing teeth and the relationship to the adjacent dentition.
3. Using an explorer, probe for any open margins between the preparation and the crown. If necessary, reshape the crown to achieve a better marginal fit.
4. If there is excessive blanching of the gingival tissue when the crown is seated, the crown will need to be trimmed. Using an indelible pen or wax pencil, mark the crown following the contour of the gingiva. Use crown scissors to trim 1 mm below the mark. Try to cut in a smooth motion to avoid producing sharp metal burs that might irritate the patient's tissue. If sharp areas are produced, use a heatless stone to remove the burs and polish with a rubber wheel. Use contouring pliers to reshape and crimp the crown at the trimmed margins. Complete another test fit to ensure accuracy, proper contacts and occlusion.
5. Complete another try-in fit to ensure accuracy, proper contact and check the occlusion.
6. When checking the occlusion, carefully notice any differences in the maxillary/mandibular relationship. If the crown is fully seated and the crown prevents proper occlusion, further reduction of the preparation may be necessary. Slightly reduce the tooth further and re-try the crown once again to confirm proper occlusion.

Seating the Crown

1. Isolate the prepared tooth by placing rubber dam or by using cotton rolls.
2. Wash and thoroughly dry the preparation.
3. Mix the cement and fill the selected crown, following cement manufacturer's instructions.
4. Seat the crown. Recommendation: for mandibular crowns start at the lingual surface, pressing across to the buccal surface; for maxillary crowns start at the buccal surface, pressing across to the lingual. This will assist in the residual cement being squeezed out. Do not use excessive force during final seating.
5. During the initial seating of the crown, verify proper alignment by checking the relationship with adjacent teeth and opposing teeth. When proper alignment is confirmed, instruct the patient to bite on gauze or an occlusal bite stick; this will ensure the crown is fully seated into the proper occlusion. If using a device to assist in crown seating during cementation, please use care to obtain proper occlusion. Do not use excessive force in order to avoid seating the crown below the level of proper occlusion.
6. Following individual cement manufacturer's instructions, remove the excess cement at the allocated time. Cement removal can be done with an explorer, scaler, or excavator. Use dental floss to remove cement from the interproximal surfaces and subgingival areas.

Hu-Friedy Stainless Steel Crowns are made with a stainless steel alloy that contains elemental nickel and chromium. A small percentage of the population is known to be allergic to these metals. If an allergic reaction occurs, direct the patient to consult a physician.

This product is intended for use only by dental professionals.

This product is single use (disposable); reuse is prohibited. Potential risks if the product was reused include:

- Introduction of microbes into the body due to reuse of item which may have been contaminated during its first use, which can result in subsequent health deterioration of patient
- Soft tissue injury from any rough edges if the crown was previously trimmed

Keep product at room temperature and out of direct sunlight.

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Hu-Friedy Mfg. Co., LLC: