English

INSTRUCTION FOR USE

PWR Piezo Implant Cleaning Set consists of two parts: a Tip Holder (UPICTH, Fig. 1 - B) and a Tip (UPIC1, Fig. 1 – A). This set has been designed for use in combination with PWR Pair Combo and PWR Piezo ultrasonic scalers distributed by HuFriedyGroup and in accordance with the intended use.

PWR Piezo Implant Cleaning Set can be used to remove plaque or biofilm from the surfaces of implants, metal and ceramic dentures and natural teeth.

The torque wrench supplied with the Tip Holder, or sold separately, must only be used to secure and remove the Tip Holder from PWR Pair Combo and PWR Piezo devices distribuited by HuFriedyGroup.

⚠ **WARNING:** The device must be used only by specialised and properly trained personnel, such as Physician, Dentist and/or Dental Hygienist, able-bodied adults of any weight, age, height, gender and nationality.

NOTE: Regarding the intended patient group, patient selection criteria, indications for use, contraindications and warnings, refer to the PWR Pair Combo and PWR Piezo devices' "USE & MAINTENANCE MANUAL".

IMPORTANT: The terms "**WARNING**" and "**CAUTION**" and "**NOTE**" identify conditions that require particular caution and must, therefore, be examined carefully. In particular:
⚠ **WARNING:** Identifies conditions and practices that present a risk of death or serious injury for the patient and/or user.
⚠ **CAUTION:** : Identifies conditions and practices that may cause minor injuries and/or damage to the devices.
NOTE: Warning not associated with a hazardous or imminently dangerous condition.

⚠ **CAUTION:** The information contained in this document is limited to that strictly necessary for the first use of the product and must be consulted as a supplement to the PWR Pair Combo and PWR Piezo devices' "USE & MAINTENANCE MANUAL". Read the PWR Pair Combo and PWR Piezo's "USE & MAINTENANCE MANUAL" paying particular attention to the "Safety Requirements" chapter, before performing any type of operation on the system and commencing any type of treatment.

⚠ **WARNING:** If a serious accident attributable to the device occurs during correct and intended use, it is recommended to report it to the Competent Authority and to the company indicated on the product label.

⚠ **WARNING: SINGLE USE ONLY**

The polymer Tip UPIC1 is intended to be used on an individual patient for single use only.
- Discard it after each use.
- Do not attempt to clean or sterilize it after use.

⚠ **WARNING: Before first use**

The Tip Holder, the Tip and the torque wrench are supplied "non-sterile", therefore, before first use, they must be reprocessed in accordance with the procedures summarily described below and reported in detail in the "Cleaning and Sterilization" section of the PWR Pair Combo and PWR Piezo devices' "USE & MAINTENANCE MANUAL".

After each use, the torque wrench and the Tip Holder must be immediately cleaned and sterilized.

1- Cleaning procedures

IMPORTANT: Submerge, the Tip Holder and/or the torque wrench in demineralised water or in an enzymatic cleaning solution immediately after use. Do not leave residue or blood deposits on the Tip Holder and torque wrench, eliminate larger impurities with a disposable cloth or paper towel.

The cleaning procedure can be performed using two different methods: manual cleaning and automatic cleaning. These two methods are alternative and only one either the other must be performed. Repeated reconditioning has a minimal effect on these devices and instruments. The end of service life is generally determined by wear or damage due to use.

⚠ **CAUTION:** Do not use hydrogen peroxide to clean the device.

⚠ **WARNING:** To avoid injury to operators, treat the Tip Holders separately. Particular attention must be paid while handling sharp and potentially infected instruments to avoid the risk of wounds and infections. Health personnel in charge of use and cleaning of the devices must use suitable protection and always use care when handling instruments with sharp edges and pointed tips.

1.1- Manual Cleaning

NOTE: The Tip Holder, the Tip, and torque wrench must be processed separately.
- Do not leave the Tip Holder inserted into the torque wrench.
- Do not leave the Tip installed on the Tip Holder.
- Do not place the Tip Holder, Tip or torque wrench in the same container.
- Do not use the same toothbrush for the Tip Holder, Tip, and torque wrench.

NOTE: Process independently validated with Enzymec enzymatic detergent and 0.8% dilution.

1 Prepare a pH-neutral (6-9) enzymatic detergent solution, following the manufacturer's instructions.

⚠ **CAUTION:** Once used, dispose of the enzymatic cleaner solution properly, do not recycle.

2 Place the Tip Holder, the Tip, or torque wrench in a clean container, in a horizontal position. Add enough of the prepared solution to completely cover the device to be cleaned.

3 Soak the device in the enzymatic detergent solution for 10 minutes at room temperature. This process reduces the amount of blood, protein, and mucus on the instrument.

4 After 10 minutes of immersion in the enzymatic solution, gently brush all surfaces until any visible dirt is removed.
For the Tip Holder and Tip: use a clean toothbrush with soft nylon bristles for the outer surfaces, a clean brush with soft nylon bristles for the inner cavities and crevices. Thoroughly clean hard-to-clean areas such as sharp edges.
For the torque wrench: Use a clean soft-bristled nylon toothbrush. Brush thoroughly, for about 20 seconds, all of the following parts of the torque wrench:
- through holes and internal channels;
- external metal ring;
- internal cavities, grooves and fissures.

5 Remove the device being cleaned from the enzymatic detergent solution.

6 Thoroughly rinse and brush all the surfaces of the device (see point 4) under running water for:
- at least 10 minutes for the wrench;
- at least 1 minute for the Tip and Tip Holder.

7 For the Tip Holder and Tip add the following steps: use a 20 ml disposable syringe to flush the enzymatic detergent solution into hard-to-reach areas (such as inner channel, through-holes or /cannulas). Repeat this step three times to ensure the effective removal of dirt from the internal surfaces. Each time use freshly prepared solution.

8 Place the device to be cleaned in a stainless-steel instrument tray with a mesh bottom. Place the tray in the ultrasonic cleaning tub and make sure that in the bath, at 75.2°F ±3.6°F (24°C ±2°C), the instruments are submerged in fresh enzymatic solution (see step 1).

9 Ultrasonicate the torque wrench for at least 20 minutes and the Tip Holder and Tip for 10 minutes, or as instructed by the manufacturer of the enzymatic cleaner and ultrasonic tank.

10 Remove the device from the ultrasonic cleaning tank and repeat step 6.

11 **For the Tip Holder and Tip:** flush three times the inner channel of the device to be cleaned with injected demineralized water using a 20 ml syringe in order to remove any dirt residues. Soak the device to be cleaned in demineralized water for 1 minute.
For the torque wrench: Soak the device in demineralized water for 5 minutes.

1.2- Automatic Cleaning

NOTE: The Tip Holder, the Tip, and torque wrench must be processed separately.
- Do not leave the Tip Holder inserted into the torque wrench.
- Do not leave the Tip installed on the Tip Holder.
- Do not place the Tip Holder, Tip or torque wrench in the same container.
- Do not use the same toothbrush for the Tip Holder, tip, and torque wrench.

Automatic cleaning involves the use of a thermodisinfector and the following materials:
- Alkaline detergent. For example, neodisher® FA (0.2 % v/v);
- Neutralizing liquid. For example, neodisher® Z (0.1 % v/v);
- Metal basket;
- Adaptors.

NOTE: Make sure that the accessories are appropriately blocked in the basket and cannot move during washing. Any blows could damage them. Position the instruments in a way that the water can flow through all the surfaces, even internal.

⚠ **WARNING:** Avoid overloading the thermal disinfector as this could compromise cleaning effectiveness.

1 Place the accessories in a metal basket using the special adaptors (supplied as optionals).

2 Set the sequence and parameters for the cycle in the thermal disinfector:
- 1 min, Rinse with cold water;
- 5 min, Wash with alkaline detergent at 131°F ± 3.6°F (55°C ±2°C);
- 1 min, Neutralization with neutralizing liquid and deionized water at 89.6°F ± 3.6°F (32°C ±2°C);
- 1 min, Rinse with deionized water at 89.6°F ± 3.6°F (32°C ±2°C);
- 5 min, Thermal disinfection at 199.4°F (93°C) with deionized water.

Thermal disinfection has not been tested experimentally. In compliance with ISO 15883-1, Table B.1 [4] thermal disinfection at a temperature of 194°F (90°C) for 5 min determines a value of A0 3000.

⚠ **CAUTION: After cleaning**
- The use of alkaline detergents removes the lubricant from the torque wrench, causing malfunctions and increasing its wear. Therefore, it is important to lubricate the torque wrench after each cleaning procedure using medical-grade lubricants (see chapter "2- Lubrication"). The dilution of the alkaline detergent required for automatic cleaning is such that it does not damage the silicone parts of the torque wrench (O-ring).
- Use filtered compressed air to flush out any moisture from cavities, grooves, crevices, and other hard-to-reach areas. This prevents the appearance of stains or streaks on the surface and oxidation. Make sure that the device is completely dry even internally before starting the sterilization cycle.
- Check that the device is actually clean and repeat the process if the water does not run clear in the final rinse, or if dirt residues are still visible.
- All parts must be inspected before their sterilization. In general, it is sufficient to conduct a visual inspection in good lighting conditions without the use of magnifying means. All parts must be checked, looking for visible dirt, damage and/or corrosion, signs of wear. Particular attention must be paid to:
 - through-holes,
 - threading
 - internal channels,
 - the edges of the sharp tips
 - edges that may be deformed or largely affected areas: the edges must be continuous.

If surfaces are not visually clean repeat cleaning and inspect again. Discard damaged tools.

2- Lubrication

Before sterilization, the torque wrench must be lubricated with a commercial medical-grade lubricant.

⚠ **CAUTION:** Do not lubricate the Tip Holder or the Tip.

The lubricant must be applied by spraying it directly onto the peripheral contact surface inside the torque wrench, as indicated in Fig. 3;
after having applied the lubricant, remove any excess of lubricant using a clean lint-free cloth.

3- Packaging

The Tip Holder, the Tip and the torque wrench must be sterilized using a standard, appropriately sized medical steam sterilization pouch.

⚠ **CAUTION:** The Tip Holder, the Tip and the torque wrench must be packaged separately. Do not package multiple tips and/or torque wrenches in the same pouch.

⚠ **CAUTION:** Make sure the pouch is big enough to hold the instrument without tensioning the seals and without ripping the packaging.

4- Sterilization

Sterilize only by pre-vacuum steam autoclaving.

The steam sterilization parameters have been validated by the manufacturer to guarantee a sterility level (Sterility Assurance Level - SAL) of 10^{-6}.

⚠ **CAUTION:** If multiple instruments need to be sterilized in the autoclave in a single cycle, do not exceed the maximum allowed load.

⚠ **CAUTION:** Risk of contamination. Do not use gravity autoclaves to sterilize the tips. The operating cycle of gravity autoclaves does not guarantee the suitable sterilization of the internal channel, of the cavities and hard-to-reach points.

⚠ **CAUTION:** Use only a pre-vacuum autoclave to sterilize the tips and the torque wrench. Do not use other methods of sterilization insofar as they may be incompatible with the materials used to produce the instruments.

Do not use the following methods of sterilization: ethylene oxide sterilization, hot air sterilization, flash autoclaving, STERRAD sterilization, sterilization with STERIS or similar sterilization systems. Do not use the following substances/systems to sterilize the instruments: hydrogen peroxide, peracetic acid, Formaldehyde, Glutaraldehyde based systems or other equivalent solutions/systems.

Minimum sterilization parameters to guarantee a sterility level (Sterility Assurance Level - SAL) of 10^{-6}:

⚠ **WARNING:** these sterilization parameters, validated by an independent laboratory, refer exclusively to a PRE-VACUUM autoclave operating cycle.
- Type of cycle: PRE-VACUUM - (3 pre-vacuum stages).
- Minimum temperature: 270 °F (132 °C) (interval between 270°F ÷ 275°F (132°C ÷ 135°C)).
- Minimum exposure time (in package): 4 minutes for the tips/tip holders and for the torque wrench.
- Minimum drying time: 10 minutes.

⚠ **CAUTION:** Do not exceed a temperature of 275°F (135 °C).

All stages of sterilization must be performed by the operator in compliance with the current revision standards: UNI EN ISO 17665-1, UNI EN ISO 556-1 and ANSI/AAMI ST:46.

5- Important information relative to safety of tips: UPIC1, UPICTH

⚠ **WARNING: Breakage and wear of tips.**

On rare occasions, the oscillations at high frequency and wear may cause a tip to break. Do not bend, reshape or sharpen the tips in any way. The alteration of a tip may cause it to break. Altered tips must never be used. If a tip breaks during use, recover and remove all fragments of the tip from the patient, with maximum care. During the treatment, ask the patient to breathe through their nose and/or use a dental dam to avoid the risk of swallowing any fragments generated by the breakage of the tip.

⚠ **WARNING:** The tips are consumable. Thoroughly inspect the tips before and during every treatment, looking for damage and/or excessive wear. Do not use a tip if it is scratched or rusted. The tip may break during use. If damage or deterioration in performance is noted, replace the tip with a new one.

NOTE: A deformed tip may preclude the complete functionality of the device used, and, depending on the model, a message may be displayed.
- Before the treatment, check that the Tip is firmly secured to the Tip Holder. Properly secure the Tip Holder onto the handpiece using the torque wrench. And check that the Tip is firmly secured on the Tip Holder.
- Do not change the shape of the tip in any way. Bending or forcing the tip may cause it to fracture. Never use bent tips.
- Do not attempt to sharpen used tips as they may break during use.
- Replace the tip/s only with original HuFriedyGroup spare parts. Use of tips other than original HuFriedyGroup tips shall void the warranty of the device. Never use tips other than original HuFriedyGroup tips insofar as they will damage the device and may cause injury to operators or the patient. The use of Tip Holder other than original HuFriedyGroup ones shall damage the threaded pin of the handpiece, thus compromising the secure attachment of tips, even if original, to the handpiece.
- It is advisable to avoid the application of excessive force or prolonged contact of the tip on the soft tissues to avoid thermal damage and / or injury.
- Let the ultrasonic vibrations work, do not exert excessive pressure on the tips during use. Apply a light force on the tip to obtain the best efficiency.
- The application of excessive pressure may cause the Tip to fracture, which may in turn cause injuries. The use of tips other than the original ones will invalidate the device warranty.
- The tip must always be kept moving. If the tip is blocked, it can cause overheating of the treated part. It is recommended to use continuous movement to minimize contact between the tip and the part. Do not block it against the tissue so as not to cause it to overheat. It is advisable to use high levels of irrigation as the power level increases.
- The tips vibrate with a longitudinal oscillation, with forward and backward movement. During treatment, always keep the instrument tangential to the tooth surface. Move the handpiece back and forth while applying light lateral pressure. Do not apply pressure to the tip in the axial direction
- Do not aim the instrument directly on the surface of the enamel or implant. Position the tip/operative part only tangentially to the surface of the tooth or implant.
- When the tip is used in the interproximal spaces, do not block the instrument or leverage the operative part. The tips must be left free to vibrate.
- Check the threaded parts of the Tip Holder and those on the handpiece. These parts must be thoroughly cleaned.
- The tips must be replaced as soon as there is visible wear.
- Use of a worn Tip shall lower its performance.

⚠ **WARNING: Hospital waste**

Treat the following objects as hospital waste:
Tips and/or Tip Holders, when worn or broken.
Torque wrench, when worn or broken.

6- Use of the torque wrench

Refer to Fig. 2

6.1- Securing the Tips

1 Tighten the Tip Holder UPICTH (Fig. 1 - B) onto the handpiece using the torque wrench (see Figure 3).

2 Select the power and irrigation levels depending on the HuFriedyGroup device to be used (see Table 1);

⚠ **WARNING:** Tip holder and Tip (UPICTH+UPIC1) must be used together respecting the "Power level" and "Irrigation level" as shown in Table 1. Irrigation must always be present and adequate for the treatment.

3 Run the FLUSH function to load the cord/handpiece/Tip system with the irrigation liquid. Or press the footswitch until the liquid starts spilling out of the Tip Holder.

⚠ **WARNING:** If the liquid is not loaded beforehand into the cord/handpiece/Tip system, there will be a risk of breakage of the UPIC1 Tip.

4 Screw the UPIC1 Tip onto the Tip Holder, position it against the surface and tighten it hard enough as allowed by the strength of your fingers. (see Fig. 1 - C).

⚠ **WARNING:** Care is required in the initial stage of screwing, when the threads must match precisely. If this is not done correctly, the thread of the UPIC1 Tip could get damaged, causing it to come unscrewed during treatment.

5 Press the footswitch to start the treatment. Do not apply excessive pressure when using the UPIC1 tip. A slight pressure will ensure greater care for the surface of the implant.

⚠ **WARNING:** During the operation, check that the UPIC1 Tip is always tightened against the Tip Holder; if the UPIC1 Tip is no longer correctly screwed into place, tighten it again hard onto the Tip Holder.

⚠ **CAUTION:** The UPIC1 Tip is prone to deterioration in the form of wear and deformation. Wear will lead to a reduction in the length of the Tip. Deformation, normally localised on the end of the tip, will make the shape unsuitable for use. In both cases, the damage will lead to a progressive drop in functional performance. When it becomes visible, it will be necessary to replace the UPIC1 Tip.

⚠ **CAUTION:** While securing the Tip Holder onto the handpiece make sure not to grip the handpiece by the end and/or cord part, but only by the body (Fig. 2 - C). Do not rotate but firmly grip the body of the handpiece, while rotating the torque wrench only. (Fig. 4)

6.2- Removing the Tip holder

Position the torque wrench on the Tip Holder and screw in an anti-clockwise direction.

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All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
© 2026 Hu-Friedy Mfg. Co., LLC. All rights reserved
Hu-Friedy Mfg. Co., LLC