

An IMS® Product

COLOR CODE TABS

Steam sterilize for at least
4 minutes at 270°F/132°C or
30 minutes at 250°F/121°C.
Do not exceed 350°F/177°C.
Recommend 30-minute dry
time after sterilization cycle.

Manufactured for:
Hu-Friedy Mfg. Co., LLC
3232 N. Rockwell St.
Chicago, IL 60618 | USA
Hu-Friedy.com

IMS-12911N | 5 TABS/OCEAN BLUE

N1232/0318



+ D6719198961P