


# FASTER, SAFER, MORE EFFICIENT—

## ALL-IN-ONE INTEGRATED CASSETTE SYSTEM



It's simple—standardizing procedural instrument setups with sterilization cassettes saves time and money. A systematic instrument management workflow, from cleaning to chairside, enhances productivity, generates revenue, and extends the life of instruments.

Infinity Series cassettes are Hu-Friedy's latest advancement to simplify workflow and enhance sterilization effectiveness. The proprietary hole pattern and rail design increases exposure to instruments during reprocessing, creating the optimal environment for cleaning and sterilization. Compatible with all cleaning equipment, including automated washers and ultrasonic units, Infinity Series cassettes are constructed of electropolished stainless steel for long-term durability and resistance to breakage. Cassettes are available in 17 configurations and 11 options for color coding.

"Hu-Friedy is dedicated to educating dental practitioners on infection prevention protocols, conducting in-practice

audits of instrument management procedures, and providing resources for infection prevention," said Hu-Friedy president Ken Serota. "Helping provide the safest environment for patients and staff is a priority for us."



**HU-FRIEDY MFG CO, INC**  
FREE Product Information  
CALL 800.420.8990 ext. 24153  
CLICK [www.dps.li/a/4QH-153](http://www.dps.li/a/4QH-153)  
Circle Reader Service 153