

# Ultrasonic Scaling Inserts



## Popular Styles For Your Power Scaling Needs

	#10 Universal	#100 Thin	#1000 Triple Bend
<b>Tip</b>	Medium length straight shank extending to a tapered tip. <b>Tip Cross Section:</b> Round	Long thin straight shank extending to a thin tapered tip. <b>Tip Cross Section:</b> Round	Unique triple bend shank design allows increased access around line angles and interproximally. Beveled edges on the tip enhance removal of heavy stain and calculus. <b>Tip Cross Section:</b> Diamond Shape
<b>Area of Use</b>	Supragingival; Subgingival only if tissue allows	Subgingival — designed specifically to reach into periodontal pockets; Also can be used supragingivally	Supragingival; subgingival only if tissue allows
<b>Deposits</b>	Light to heavy, tenacious calculus, orthodontic cement	Light to Moderate	Moderate to heavy, tenacious calculus, orthodontic cement
<b>Power Setting</b>	Low to High	Low to Medium	Low to High

	XT	XT Triple Bend	#3 Beavertail
<b>Tip</b>	Long, ultra thin curved shank tapers to a thin tip. <b>Tip Cross Section:</b> Round	24% thinner than the #1000 Triple Bend, allowing increased access around line angles and interproximally. Beveled tip allows enhanced calculus and stain removal. <b>Tip Cross Section:</b> Diamond Shape	Thick, broad, flat tip <b>Tip Cross Section:</b> Trapezoidal
<b>Area of Use</b>	Subgingival; ideal for narrow, deep pockets, furcations, developmental concavities and tight areas between misaligned teeth	Supra- and subgingival	Supragingival
<b>Deposits</b>	Light to Moderate	Light to Moderate	Heavy deposits, tenacious stain and orthodontic cement
<b>Power Setting</b>	Low to Medium	Low to Medium	Medium to High

	After Five™ Left, Straight & Right	Power Plus Beveled & Conical
<b>Tip</b>	Probe sized insert tip; Long curved shank narrows to a tapered tip. <b>Tip Cross Section:</b> Round	25% longer tip design to reach deep into periodontal pockets. 2mm longer than #10 Insert. <b>Tip Cross Section:</b> Round (Conical) Diamond Shape (Beveled)
<b>Area of Use</b>	Used to remove calculus from periodontal pockets with improved adaptation in interproximal areas, furcations and developmental concavities	Supragingival and subgingival; wide, deep pockets (standard), narrow, deep vertical pockets and crowded areas (thin)
<b>Deposits</b>	Light to Moderate	Moderate to Heavy (standard), Light to moderate (thin)
<b>Power Setting</b>	Low to Medium	Low to High (standard), Low to Medium (thin)

HuFriedyGroup Inserts are compatible with the SWERV™ 3 Scaler as well as all Cavitron® and Cavitron-compatible generators.\*

# SELECT YOUR HANDLE & WATER DELIVERY STYLE

		Swivel Direct Flow™ Handle	Swivel Handle	Streamline™ Plus Handle	Streamline™ Direct Flow Handle	Streamline™ Handle
		<p>With a slight rolling motion of the fingertips, the handle and tip of a Swivel Insert rotates 360°, allowing for smooth and efficient power scaling.</p> <p><b>Water:</b> Direct flow near the insert tip for enhanced visibility.</p>	<p>The same features as Swivel Direct Flow™ Handles</p> <p><b>Water:</b> Deliver from the base of the tip.</p>	<p>The ergonomic, wide grip results in less hand fatigue and enhanced tactile sensitivity. The tapered cone provides a better view of the treatment area for the clinician.</p> <p><b>Water:</b> Direct flow near the insert tip for enhanced visibility</p>	<p>A lightly textured, wide diameter resin handle designed to reduce finger pinching and enhance clinician comfort.</p> <p><b>Water:</b> Direct flow near the insert tip for enhanced visibility.</p>	<p>The original lightly textured Streamline handle. For those who prefer a slightly smaller diameter handle the Streamline handle is an ideal choice.</p> <p><b>Water:</b> Delivery from the base of the tip.</p>
						  25kHz    30kHz
<b>#10 Universal</b> (Lavender)	30kHz 25kHz	<b>UI30SD10</b> <b>UI25SD10</b>	<b>UI30SS10</b> <b>UI25SS10</b>	<b>UI30SP10</b>	<b>UI30SF10</b> <b>UI25SF10</b>	<b>UI1030KX</b> <b>UI1025K</b>
<b>#100 Thin</b> (Black)	30kHz 25kHz	<b>UI30SD100</b> <b>UI25SD100</b>	<b>UI30SS100</b> <b>UI25SS100</b>	<b>UI30SP100</b>	<b>UI30SF100</b> <b>UI25SF100</b>	<b>UI30K100SX</b> <b>UI10025K</b>
<b>#1000 Triple Bend</b> (Orange)	30kHz 25kHz	<b>UI30SD1000</b> <b>UI25SD1000</b>	<b>UI30SS1000</b> <b>UI25SS1000</b>	<b>UI30SP1000</b>	<b>UI30SF1000</b> <b>UI25SF1000</b>	<b>UI100030KX</b> <b>UI100025K</b>
<b>XT</b> (Purple)	30kHz 25kHz		<b>UI30SSXT</b> <b>UI25SSXT</b>		<b>UI30SFXT</b>	
<b>XT Triple Bend</b> (Purple)	30kHz	<b>UI30SDXTTB</b>		<b>UI30SPXTTB</b>	<b>UI30SFXTTB</b>	
<b>After Five™ Left</b> (Teal)	30kHz 25kHz	<b>UI30SDL</b> <b>UI25SDL</b>		<b>UI30SPL</b>	<b>UI30SFL</b> <b>UI25SFL</b>	
<b>After Five™ Straight</b> (Blue)	30kHz 25kHz	<b>UI30SDS</b> <b>UI25SDS</b>		<b>UI30SPS</b>	<b>UI30SFS</b> <b>UI25SFS</b>	
<b>After Five™ Right</b> (Red)	30kHz 25kHz	<b>UI30SDR</b> <b>UI25SDR</b>		<b>UI30SPR</b>	<b>UI30SFR</b> <b>UI25SFR</b>	
<b>Beavertail</b> (Yellow)	30kHz 25kHz					<b>UI330KX</b> <b>UI325K</b>
<b>Power Plus</b>			<b>UI30PPSC</b> <b>UI30PPTC</b> <b>UI30PPSB</b> <b>UI30PPTB</b>			

Visit us at [HuFriedyGroup.com/Inserts](https://www.HuFriedyGroup.com/Inserts)

\*Cavitron is a registered trademark of DENTSPLY International, INC.  
 Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com  
 All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.  
 ©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-506/0823