



SPORECHECK 24 HOUR IN-OFFICE TEST®


ACCURATE. EASY TO USE. EFFICIENT.



Hu-Friedy's **biological monitoring system** offers a quick and accurate approach to in-office steam sterilizer monitoring. The SporeCheck system provides reliable **test results in 24 hours** while being in compliance with local, state and federal guidelines.

SporeCheck is part of the Sterilization step in Hu-Friedy's Instrument Management System (IMS®)

This 4-step process is designed to keep instruments organized and intact from chairside, through cleaning and sterilization, and back chairside. IMS provides optimal infection control and prevention, resulting in improved efficiency and safety.



How the best perform



POINTS OF PERFORMANCE

- SporeCheck helps you follow CDC guidelines, which recommend that steam sterilizers be monitored at least weekly, preferably daily, as well as with every load that contains implantable devices.
- In-office biological monitoring allows you to take control of your steam sterilizer testing and provides you with the peace of mind that your sterilizer is functioning properly with results in 24 hours.
- The SporeCheck system is easy to use, and the biological indicator (BI) vials work for ALL types of steam sterilizers, including gravity, pre-vacuum and flash sterilization.

SPORECHECK STARTER KIT

Contains: Dry Block Incubator with built-in vial crusher, record notebook and 1 box SporeCheck Biological Indicators, 25 count

IMS-1376

Starter Kit pictured below

INDIVIDUAL PRODUCTS

SporeCheck Biological Indicator Vials, 25 count IMS-1373

SporeCheck Biological Indicator Vials, 100 count IMS-1373H

55-60° Dry Block Incubator IMS-1374


Sterilizer Record Notebook IMS-1375

Replacement Vial Crusher IMS-1377

RECORD NOTEBOOK



SELF-CONTAINED
STEAM BIOLOGICAL
INDICATORS (25 VIALS)



TEST RESULTS ARE EASY
TO INTERPRET WITH A
VISIBLE COLOR CHANGE FOR
POSITIVE RESULTS, LEAVING
NO ROOM FOR ERROR



CONVENIENT, BUILT-IN
VIAL CRUSHER



DRY BLOCK INCUBATOR WITH
BUILT-IN VIAL CRUSHER

Visit us online at Hu-Friedy.com

Call 1-800-Hu-Friedy or contact your authorized Hu-Friedy representative for more information.
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. IMS-355/0212.

