




IMS[®] INFINITY SERIES[™] CASSETTES

Function and fashion combine to create Hu-Friedy's IMS Infinity Series Cassettes. The innovative hole pattern and rail design significantly increase exposure to instruments during reprocessing, creating an optimal environment for cleaning and sterilization.

The sophisticated design of the Infinity Series Cassettes brings a practice to the next level in professionalism and patient care.



POINTS OF PERFORMANCE

Ergonomic latch with positive locking mechanism is safe and easy-to-use, allows one-handed opening and visual confirmation if the cassette is locked / unlocked.

Innovative color-coded silicone rail system, available in 11 color options, significantly reduces instrument contact by up to 59%, allowing for more water flow while protecting the instruments during reprocessing.

Revolutionary new hole pattern design creates 30% additional open area for more access during cleaning and sterilization, and is compatible with current cleaning equipment, including automated washers.

High quality, electropolished stainless steel protects against corrosion while providing a professional appearance.

Smooth round corners and slotted edges increase drainage and reduce drying time.

New beveled edge allows for a more ergonomic grasp.

How the best perform


Hu-Friedy

REFERENCE SECTION

WHY IMS?

Hu-Friedy's Instrument Management System (IMS) standardizes procedures to keep instruments organized and intact from cleaning to chairside, providing opportunity for:

- 5-10 minutes saved per procedure
- Increased staff safety
- Increased instrument life
- Focus on revenue-generating activities
- Reduction in overtime
- Potential to see 1-2 additional patients per day



CASSETTE COLOR CODING

	0 Red		4 Purple (new)		8 Blue
	1 Gray		5 Yellow		9 Green
	2 White		6 Lavender		11 Ocean Blue (new)
	3 Orange		7 Black (new)		

INFINITY SERIES CASSETTE PART CODES

DESCRIPTION	DIMENSIONS	ACCESSORY AREA	PART CODE
5 Instrument Exam Cassette	3" x 8" x 1.25"		IMN605X
8 Instrument Cassette	5.5" x 8" x 1.25"	Yes	IMN508X
10 Instrument Cassette	5.5" x 8" x 1.25"		IMN510X
12 Inst Signa Stat Cassette	6.5" x 10.5" x 1.25"	Yes	IMN812X
12 Instrument Cassette	8" x 11" x 1.25"	Yes	IMN412X
16 Instrument Cassette	8" x 11" x 1.25"	Yes	IMN416X
20 Instrument Cassette	8" x 11" x 1.25"		IMN420X
8 Instrument Double-Decker Cassette	3" x 8" x 1.5"		IMN908X
14 Instrument Double-Decker Cassette	4.5" x 8" x 1.5"		IMN914X
20 Instrument Double-Decker Cassette	6" x 9.7" x 1.5"	Yes	IMN920X
Oral Surgery Cassette	8" x 14.5" x 1.25"	Yes	IMN0SX
Small Ultrasonic Insert/Handpiece Cassette	1.5" x 8" x 1.25"		IMN30HUX
Large Ultrasonic Insert/Handpiece Cassette	3" x 8" x 1.25"		IMN60HX
Radiography Cassette	5.5" x 8" x 1.25"		IMN7000
Large Thin Orthodontic Cassette	8" x 11" x 0.75"		IMNORTHL8
Small Thin Orthodontic Cassette	7" x 8" x 0.75"		IMNORTHS8
Small Double-Decker Ortho Cassette	4.5" x 8" x 1.5"		IMN914X-OR

Visit us online at [Hu-Friedy.com](https://www.hu-friedy.com)

Call 1-800-Hu-Friedy or contact your authorized Hu-Friedy representative for more information.
© 2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. IMS-101/1114.

