



TEMPORARY CROWN REMOVER

FACT


Temporary crowns are necessary in dentistry in order to preserve a site as it awaits its final restoration. This crown plays a vital role in sustaining proper tooth functionality and maintaining a patient's appearance.

CHALLENGE

Separating a temporary crown from its provisional cement can present a difficult task, especially given the environment of the surrounding oral cavity.

SOLUTION

Hu-Friedy's slender Temporary Crown Remover provides an improved grip, enhanced visibility and greater patient comfort.



CRGR - TEMPORARY CROWN REMOVER



POINTS OF PERFORMANCE

- **Carbide inserts** for enhanced grip
- **60 degree curved beak** offers improved intraoral visibility
- **Slender design** provides increased patient comfort
- **Gold-plated handle** for quick and easy identification



Visit us online at [Hu-Friedy.com](https://www.hu-friedy.com)

Call 1-800-Hu-Friedy or contact your authorized Hu-Friedy representative for more information.
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. M7540/0613

How the best perform

