



## SWERV<sup>3</sup> MAGNETOSTRICTIVE ULTRASONIC SCALER

COMFORT FOR YOU AND YOUR PATIENTS



### FACT

For optimal clinical outcomes, clinicians use both hand and ultrasonic scaling instrumentation in order to effectively remove calculus, plaque and debris.

### CHALLENGE

Ultrasonic scalers need to effectively remove calculus during each procedure, plus be comfortable for both the practitioner and the patient.

### SOLUTION

Hu-Friedy's magnetostrictive ultrasonic scaler combines power and handling in a unique ergonomic design that ensures the comfort of dental patients and professionals. It integrates finely tuned electronics that improve control and scaling efficiency in all clinical indications.



“ I LOVE HU-FRIEDY'S SWERV<sup>3</sup>. This new technology creates a smooth ride for the patient and the clinician. Patients will notice a huge improvement in comfort and sound. Plus, Swivel™ Insert technology is now a luxury I could never go without. It makes adapting the insert so much easier and more efficient. ”

– Leah MacPherson, RDH, BS, MHP  
Professor, Middlesex Community College

How the best perform




## POINTS OF PERFORMANCE

- Finely tuned electronics provide superior patient and clinician comfort
- Switching power settings is effortless with the color-coded dual power modes
- The smooth surface makes disinfection simple and effective
- The color-coded, lighted display is easy to identify and adjust with the touch of a button
- The auto-purge function helps remove air from waterlines, saving the clinician time between patients

## SWERV<sup>3</sup> PART-CODES

SWERV<sup>3</sup> Series, 30K | UM3-30

SWERV<sup>3</sup> Series, 25K | UM3-25



LIGHTWEIGHT,  
ERGONOMIC  
HANDPIECE

## HU-FRIEDY ULTRASONIC INSERTS

QUALITY, PERFORMANCE & COMFORT

HU-FRIEDY OFFERS A FULL LINE OF TIP DESIGNS FOR SUPRA- AND SUBGINGIVAL POWER SCALING NEEDS, FEATURING COLOR CODED HANDLES FOR EASY IDENTIFICATION IN THE OPERATORY.

### SWIVEL AND SWIVEL DIRECT FLOW

With a slight rolling motion of the fingertips, the handle and tip of a SWIVEL insert rotates 360°, allowing for smooth and efficient power scaling.

Swivel's lightly textured, wide diameter silicone grip also provides a secure, comfortable grasp.



### STREAMLINE DIRECT FLOW

For clinicians who prefer a more traditional handle, the lightly textured wide diameter resin handle is designed to reduce finger pinching and enhance clinician comfort.



Find out more about Hu-Friedy's quality line of inserts at [Hu-Friedy.com](https://www.hu-friedy.com)

Call 1-800-Hu-Friedy or contact your authorized Hu-Friedy representative for more information.  
© 2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-422/0615

