



STREAMLINE DIRECT FLOW™ INSERTS

NEW FOR MAGNETOSTRICTIVE ULTRASONIC SCALING

POINTS OF PERFORMANCE

- A **comfortable, wide diameter handle** reduces finger pinching to ensure clinician comfort.
- **Through-tip water delivery** offers a targeted water flow, reduces excess spray and increases visibility of the treatment site.
- **Proprietary stainless steel alloy** ensures the highest level of scaling efficiency – exactly what you’ve come to expect from Hu-Friedy.
- **Variety of tips** to choose from to help optimize clinical care.



“ The Streamline Direct Flow XT is my favorite insert. I like it most for use in routine periodontal maintenance patients and for accessing deep chronic periodontal pockets and furcations. **THE XT HAS VERY GOOD SCALING EFFECTIVENESS AND IS COMFORTABLE FOR THE CLINICIAN.** Water delivery and adaptability are also very good. ”

– Deb Beres, RDH

“ As a dental hygienist with more than 30 years of experience, I love to try new products. **I REALLY LOVE** the new #10 Universal insert from Hu-Friedy! ”

– Amy Hazelwood, RDH, BS

“ Not only is the Streamline Direct Flow XT insert from Hu-Friedy extremely thin, it also **EFFICIENTLY REMOVES CALCULUS DEPOSITS.** ”

– Kimberly Long, RDH

How the best perform



STREAMLINE DIRECT FLOW INSERTS



25 kHz Part Code:
30 kHz Part Code:

#10 Universal
| UI25SF10
| UI30SF10

#100 Thin
| UI25SF100
| UI30SF100

#1000 Triple Bend
| UI25SF1000
| UI30SF1000

XT
| UI30SFXT

XT Triple Bend
| UI30SFXTTB

POWER LEVEL

LOW—HIGH

LOW—MEDIUM

LOW—HIGH

LOW—MEDIUM

LOW—MEDIUM

- Efficiently removes moderate to heavy deposits
- A great general-purpose insert

- Great for light to moderate deposits
- Thin tip design allows for subgingival access

- Designed to remove tenacious, heavy calculus
- Good for stains
- Beveled-edge design helps remove calculus quickly


- Tough enough to remove tenacious calculus
- Thin enough to scale subgingivally and interproximally

- Light to moderate deposits even in periodontal pockets
- Beveled edges for efficient calculus removal
- 24% thinner than #1000 Triple Bend for easier access

AFTER FIVE STREAMLINE DIRECT FLOW INSERTS

Ideal for:

- Light to moderate deposits
- Subgingival access
- Furcation areas
- Developmental concavities



25 kHz Part Code:
30 kHz Part Code:

Left
| UI25SFSL
| UI30SFSL

Straight
| UI25SFSS
| UI30SFSS

Right
| UI25SFR
| UI30SFR

Set of Three
| UI25SFSET
| UI30SFSET

POWER LEVEL

LOW—MEDIUM

LOW—MEDIUM

LOW—MEDIUM

LOW—MEDIUM

Visit us online at Hu-Friedy.com

Call 1-800-Hu-Friedy or contact your authorized Hu-Friedy representative for more information.
© 2020 Hu-Friedy Mfg. Co., LLC. All rights reserved. *Cavitron is a registered trademark of Dentsply International Corporation.
HF-509/0315

