

# HU-FRIEDY DIAGNOSTIC INSTRUMENTS

| It All Starts Here



**PROBES:** Qulix™ probes are color-coded with accurate, smoothly finished black markings, which will not chip, flake or fade.

**EXPLORERS:** Hu-Friedy explorers are hand crafted from Duraspond™ steel, a specially blended steel alloy which provides flexibility and resilience for even the most precise diagnostic procedures.

**EXPROS:** A double-ended instrument with an explorer on one end and a probe on the other, Expros are excellent for dental examinations and charting.

**MIRRORS:** All Cone Socket Mouth Mirrors feature rhodium coating to resist scratching, front surface for sharp and distortion free images and smooth finish for patient comfort.

**DRESSING PLIERS:** Used to grasp or transfer materials in and out of the oral cavity, all Hu-Friedy dressing pliers are made of heavy-gauge, heat-treated stainless steel to prevent tip bending and misalignment.

## 5 + 1 DIAGNOSTIC

Buy any 5 Diagnostic Instruments, **RECEIVE 1 FREE\*\***

Promo code: 2063

Visit us online at [Hu-Friedy.com](http://Hu-Friedy.com)

How the best perform



©2026 Hu-Friedy Mfg. Co., LLC. \*Offer valid in the 50 United States and District of Columbia only. Does not apply to school, government or institution offers. Offer valid through December 31, 2013.