



DURASPOND™ EXPLORERS

PRECISION WITH A TWIST



THIS EXPLORER WAS TIED IN A KNOT TO DEMONSTRATE THE SUPERIOR QUALITIES OF DURASPOND STEEL.

FACT

Explorers are used for caries and calculus detection, exploration of pocket characteristics, furcations and restorations, making them crucial to an instrument set-up.

CHALLENGE

In order for a clinician to accurately chart and diagnose a patient during their regular appointment, they need to rely on the resilience and precision of their explorer to deliver optimal results.

SOLUTION

Hu-Friedy explorers are hand crafted from Duraspond steel, a specially blended steel alloy which provides flexibility, resilience and enhanced tactile sensitivity for even the most precise diagnostic procedures. Available in over 30 styles with several handle options, there is a Hu-Friedy explorer for every clinician and every patient.

SATIN STEEL XTS® EXPLORERS

Non-stick black XTS coating makes instrument perfect for contouring anatomy during composite restorations and reduces glare



| EXD5XTS

5 XTS Explorer (shown)
Handle: #6 Satin Steel

| EXD3CHXTS

3CH Cowhorn XTS Explorer
Handle: #6 Satin Steel

| EXD2XTS

2 Pigtail XTS Explorer
Handle: #6 Satin Steel

| EXTU1723XT

TU 17/23 XTS Explorer
Handle: #6 Satin Steel

| XP23/12XTS

23/CP-12 XTS Color-Coded Expro
Handle: #6 Satin Steel

AFTER FIVE® EXPLORERS

Longer terminal shank allows for easier access to deeper pockets



| EXD11/12A6

Terminal Shank is 3 mm longer than popular ODU11/12 Explorer
Handle: Available in all Explorer handle options

How the best perform

Hu-Friedy®

POPULAR EXPLORERS



| EXD5

- #17 end perfect for subgingival calculus detection
- #23 end ideal for caries detection
- Also available in Wilkins/Tufts design



| EXS23

- Ideal for caries detection



| EXD11/12

- Ideal for posterior calculus detection in deep pockets especially on proximal surfaces
- Designed to explore entire dentition with a single instrument



| EXD3CH

- Paired curled working ends for easy access to interproximal areas

EXPLORER HANDLE OPTIONS



#8 ResinEight



#7 Satin Steel Colours



#6 Satin Steel



#31 Standard Diagnostic

Visit us online at [Hu-Friedy.com](https://www.hu-friedy.com)

Call 1-800-Hu-Friedy or contact your authorized Hu-Friedy representative for more information.
© 2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-247/1213

