



Hu-Friedy
**INFECTION CONTROL
STARTER PACKAGES**

Hu-Friedy offers complete Infection Control Starter Packages that include everything you need to ensure a safer and smarter practice for you and your patients from day one. Simply choose the appropriate package based on whether or not you are a cassette user.

With each package, you will receive some of the most popular Hu-Friedy sterilization, cleaning and care, and hand hygiene products, all for a discounted price!

How the best perform



The Infection Control Starter Packages feature the following products:



STERILIZATION PRODUCTS

SporeCheck 24 Hour In-Office Test®

The easy-to-use in-office biological monitoring system provides you with the peace of mind that your sterilizer is functioning properly with results in 24 hours. It can be used with all types of steam sterilizers including gravity, pre-vacuum and Flash Sterilizers. Kit includes: self-contained, steam biological indicators (box of 25), record notebook and dry block incubator with vial crusher.

Available on ICSP, ICSPCONSBL and ICSPCASS



Universal Cassette Sterilization Wrap

High performance, fiber free wrap at a lower cost for complete sterile protection during the sterilization cycle and during storage.

Available on ICSP and ICSPCASS



Bagette® Sterilization Pouches

Constructed of high quality materials and thick seals to help prevent instrument and cassette tear through with built in internal and external process indicators to provide time and temperature verification.

Available on ICSPCONSBL



Steam Sterilization Integrators

Integrators are a Type 5 sterilization monitor designed to react to all critical variables in the sterilization cycle: time, temperature and presence of steam.

Available on ICSP, ICSPCONSBL and ICSPCASS



IMS® Autoclave Monitor Tape and Dispenser

Chemical indicator tape for use with steam autoclaves includes 60 yards (55 m) per roll. 3/4" (19 mm) wide. The IMS Tape Dispenser holds up to 9 rolls of IMS Monitor Tape.

Available on ICSP and ICSPCASS



CLEANING & CARE PRODUCTS



Enzymax Earth™ Products

Specially formulated for efficient ultrasonic cleaning or pre-soaking of cassettes and instruments, plus a safer product for the environment. Contains Steelgard™ agents and corrosion inhibitors to protect instruments and cassettes from hard water deposits and corrosion.

Available on ICSP, ICSPCONSBL and ICSPCASS



AdvantaClear® Surface Disinfectant

Designed to clean and disinfect dental surfaces while working extremely fast and effectively to kill a broad spectrum of disease-causing microorganisms. Safe on hard, non-porous healthcare surfaces, including glass, plastics, stainless steel, and vinyl.

Available on ICSP, ICSPCONSBL and ICSPCASS



Hu-Friedy Waterline Filter

The Hu-Friedy Waterline Filter continuously works to deliver clean water to patients and staff for up to one year. The filter features a convenient, ready-to-use cartridge that is easy to install and inhibits microbial contamination in the waterline. With a universal design, it can be used with either a municipal waterline or an independent water reservoir. The filter is non-corrosive and does not contain silver or other hazardous materials, allowing for quick and easy disposal.

Available on ICSP, ICSPCONSBL and ICSPCASS



Enzymax® Spray Gel

Ready-to-use instrument pre-cleaner maintains moisture to help prevent debris from drying on instruments before cleaning. Can also be used as a linen spot remover.

Available on ICSP, ICSPCONSBL and ICSPCASS



Cleaning Monitor and Holder

Monitors the efficiency of the cleaning process in ultrasonic cleaners through a non-toxic, synthetic test soil that mimics the presence of blood and tissue. Assesses variables including cavitation, time, temperature and detergent. Test strips should be used in conjunction with the holder to provide an accurate result.

Available on ICSP

HAND HYGIENE PRODUCTS



Hu-Friedy Hand Essentials™ Skin Repair Cream

Scientifically proven to help restore and maintain skin's natural moisture balance while delivering vital nutrients. Contains Olivamine®, a patented blend of ingredients to aid in repairing the skin. Will not interfere with the antimicrobial properties of CHG (Chlorhexidine Gluconate) and is latex-friendly.

Available on ICSP, ICSPCONSBL and ICSPCASS



Hand Essentials Instant Hand Antiseptic

Hand Essentials Instant Hand Antiseptic is an essential germ-killing aid, containing 70% (v/v) ethyl alcohol to disinfect hands of most common disease-causing bacteria. The fragrance-free formula contains aloe vera and vitamin E to moisturize and soothe the skin.

Available on ICSP, ICSPCONSBL and ICSPCASS



Sterillium® Comfort Gel™ Hand Antiseptic

Scientifically proven, moisturizing, easy-to-use gel that not only meets all test and key use criteria - it also far exceeds them. Kills ≥99.999% of a broad range of nosocomial pathogens within 15 seconds.

Available on ICSP, ICSPCONSBL and ICSPCASS

Hu-Friedy's Starter Packages offer you three choices for your everyday instrument processing and infection control products. Simply choose the package that best fits the needs of your practice.

CHOOSE THE APPROPRIATE PACKAGE TO FIT YOUR NEEDS

Cassette User | ICSPCASS

Cassette Infection Control Starter Package

| Description | Part Code |
|---|----------------------------|
| IMS Universal Wrap 20"x20", 100 Sheets | IMS-1216H |
| IMS Universal Wrap 15"x15", 100 Sheets | IMS-1215H |
| IMS Lead/Latex Free Autoclave Mon. Tape | IMS-1460 |
| IMS Integrator Strips, 100/Pkg | IMS-1241 |
| SporeCheck Starter Kit | IMS-1376 |
| Enzymax Earth, 1 Gallon Bottle | IMS-1336 |
| Pump for Enzymax Gallon Bottle | IMS-1226P |
| Enzymax Earth Powder Sample Pack (Qty 2) | IMS-1332SP |
| Enzymax Spray Gel, 24 oz | IMS-1229 |
| AdvantaClear Surface Disinfectant Wipes, 160/Cstr (Qty 3) | IMS-2160 |
| Hu-Friedy Waterline Filter | IMS-1440 |
| Hand Essentials Skin Repair Cream, 2 oz | IMS-1502 |
| Sterillium Comfort Gel, 475 mL | IMS-1504 |
| Hand Essentials Instant Hand Antiseptic, 16 oz | IMS-1505 |

Non-Cassette User | ICSPCONSBL

Infection Control Starter Package

| Description | Part Code |
|---|----------------------------|
| IMS Bagette 3-1/2"x 9", 200 Pouches | IMS-1236 |
| IMS Integrator Strips, 100/Pkg | IMS-1241 |
| SporeCheck Starter Kit | IMS-1376 |
| Enzymax Earth, 1 Gallon Bottle | IMS-1336 |
| Pump for Enzymax Gallon Bottle | IMS-1226P |
| Enzymax Earth Powder Sample Pack (Qty 2) | IMS-1332SP |
| Enzymax Spray Gel, 24 oz | IMS-1229 |
| AdvantaClear Surface Disinfectant Wipes, 160/Cstr (Qty 3) | IMS-2160 |
| Hu-Friedy Waterline Filter | IMS-1440 |
| Hand Essentials Skin Repair Cream, 2 oz | IMS-1502 |
| Sterillium Comfort Gel, 475 mL | IMS-1504 |
| Hand Essentials Instant Hand Antiseptic, 16 oz | IMS-1505 |

Cassette User | ICSP

Cassette Infection Control Starter Package, Expanded

| Description | Part Code |
|---|----------------------------|
| IMS Universal Wrap 20"x20", 500 Sheets | IMS-1216 |
| IMS Universal Wrap 15"x15", 1000 Sheets | IMS-1215 |
| IMS Lead/Latex Free Autoclave Mon. Tape (Qty 2) | IMS-1460 |
| IMS Monitor Tape - Exam | IMS-1256 |
| IMS Monitor Tape - Prophy | IMS-1259 |
| IMS Monitor Tape - Restorative | IMS-1247 |
| IMS Monitor Tape - Endo | IMS-1255 |
| IMS Multiple Tape Dispenser, 9 Roll Cap. | IMS-1270 |
| IMS Integrator Strips, 100/Pkg | IMS-1241 |
| SporeCheck Starter Kit | IMS-1376 |
| Enzymax Earth, 1 Gallon | IMS-1336 |
| Pump for Enzymax Gallon Bottle | IMS-1226P |
| Enzymax Earth Powder Sample Pack | IMS-1332SP |
| Enzymax Spray Gel, 24 oz | IMS-1229 |
| AdvantaClear Surface Disinfectant Wipes, 160/Cstr (Qty 6) | IMS-2160 |
| Hu-Friedy Waterline Filter | IMS-1440 |
| Cleaning Monitor Holder | IMS-1200H |
| Cleaning Monitors - Ultrasonic, 50/Pkg | IMS-1200U |
| Hand Essentials Skin Repair Cream, 4 oz | IMS-1501 |
| Sterillium Comfort Gel, 16.06 oz | IMS-1504 |
| Hand Essentials Instant Hand Antiseptic, 16 oz | IMS-1505 |

Please note contents of starter packages are subject to change. To verify the most current configuration or for more information about Hu-Friedy, please visit

[Hu-Friedy.com](https://www.hu-friedy.com)

Hu-Friedy Mfg. Co., LLC: 3232 N. Rockwell St. | Chicago, IL 60618 | USA
Telephone: 1-800-Hu-Friedy or 1-773-975-6100 | Website: [Hu-Friedy.com](https://www.hu-friedy.com)

