

IMS™ Signature Series™ Air/Water Syringe Double-Decker™ Cassette Clip

REF IM1006

QTY 1

Stainless Steel / Acier Inoxydable

STERILIZE BEFORE USE / STÉRILISER AVANT UTILISATION

Recommended for use with IMS™ Signature Series™ Cassettes only.

Est recommandé pour une utilisation uniquement avec les cassettes IMS Signature Series.




Sterilizable in steam sterilizer
at 134 °C / 273 °F for 4 min



Washer-disinfector for thermal
disinfection possible

∴HuFriedyGroup

The Best In Practice



www.hufriedygroup.com/Reprocessing/
1-800-Hu-Friedy; info@hu-friedy.com

INSTRUCTIONS

Use clip on side of cassette. Insert air/water syringe tip into clip.

Utiliser le clip sur le côté de la cassette. Insérer l'embout de la seringue air/eau dans le clip.



Made in USA

Hu-Friedy



Hu-Friedy Mfg. Co., LLC
2501 Sanders Rd.
Northbrook, IL 60062 | USA



All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC
its affiliates or related companies, unless otherwise noted.

©2025 Hu-Friedy Mfg. Co., LLC. All rights reserved.

B7511 / REV 7 - 2508