

IMS™ Hinged Instrument Clip Clip d'instrument à charnière

REF IM1000
Stainless Steel
acier inoxydable

QTY 2



Sterilizable in steam sterilizer
at 134 °C / 273 °F for 4 min



Washer-disinfector for thermal
disinfection possible

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC
its affiliates or related companies, unless otherwise noted.
©2025 Hu-Friedy Mfg. Co., LLC. All rights reserved.

∴HuFriedyGroup
The Best In Practice


INSTRUCTIONS

Use clip on lid of cassette. Utiliser le clip dans le couvercle de la cassette
Clip is recommended for use with IMS™ Infinity & IMS™ Signature Series™ Cassettes only.
L'utilisation de ce clip est recommandée uniquement avec les cassettes IMS Signature Series et Infinity

Hu-Friedy



Hu-Friedy Mfg. Co., LLC
2501 Sanders Rd.
Northbrook, IL 60062 | USA




Made in USA



B2932
REV 13 - 2508



www.hufriedygroup.com/Reprocessing/
1-800-Hu-Friedy; info@hu-friedy.com