

QTY 1

∴HuFriedyGroup
The Best In Practice

Hu-Friedy

Hu-Friedy Mfg. Co., LLC
2501 Sanders Rd.
Northbrook, IL 60062 | USA



Sterilizable in steam sterilizer
at 134° C / 273° F for 4 min



Washer-disinfector for thermal
disinfection possible



www.hufriedygroup.com/Reprocessing/
1-800-Hu-Friedy; info@hu-friedy.com



All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC
its affiliates or related companies, unless otherwise noted.
©2025 Hu-Friedy Mfg. Co., LLC. All rights reserved.

C1091 / REV 12 - 2508