

Air/Water Syringe Clip Clips de seringue air/eau

for IMS™ Resin Cassettes
pour cassettes de résine IMS

REF **IMS-1278**
Stainless Steel
Acier Inoxydable

QTY 1
Clean and sterilize before use.
Nettoyer et stériliser avant utilisation.




Sterilizable in steam sterilizer
at 134 °C / 273 °F for 4 min



Washer-disinfector for thermal
disinfection possible

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC
its affiliates or related companies, unless otherwise noted.
©2025 Hu-Friedy Mfg. Co., LLC. All rights reserved.

HuFriedyGroup
The Best In Practice



www.hufriedygroup.com/Reprocessing/
1-800-Hu-Friedy; info@hu-friedy.com

INSTRUCTIONS

Place the flat side on the
outside of the cassette.

Placez le côté plat à
l'extérieur de la cassette



Hu-Friedy



Hu-Friedy Mfg. Co., LLC
2501 Sanders Rd.
Northbrook, IL 60062 | USA

A4435
REV 6 - 2508



Made in USA

