

Air/Water Syringe Clip

clip de seringue air/eau

for IMS™ Resin Cassettes | pour cassettes de résine

REF IMS-1279
Stainless Steel | acier inoxydable

QTY 1



Sterilizable in steam sterilizer
at 134° C / 273° F for 4 min




Washer-disinfector for thermal
disinfection possible

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC
its affiliates or related companies, unless otherwise noted.

©2025 Hu-Friedy Mfg. Co., LLC. All rights reserved.

..HuFriedyGroup
The Best In Practice



www.hufriedygroup.com/Reprocessing/
1-800-Hu-Friedy; info@hu-friedy.com

INSTRUCTIONS

Place the flat side on the outer of
the resin cassette.

placez le côté plat à l'extérieur de
la cassette de résine

Hu-Friedy



Hu-Friedy Mfg. Co., LLC
2501 Sanders Rd.
Northbrook, IL 60062 | USA



Made in USA



B7847 / REV 5 - 2508