

IMS™ Needle Cap Holder

REF IM1001

Stainless Steel
Acier Inoxydable

QTY 2

Clean and sterilize before use.
Nettoyer et stériliser avant utilisation.

Recommended for use with cassettes only.

Est recommandé pour une utilisation avec des cassettes uniquement.



www.hufriedygroup.com/Reprocessing/
1-800-Hu-Friedy; info@hu-friedy.com



Sterilizable in steam sterilizer
at 134 °C / 273 °F for 4 min



Washer-disinfector for thermal
disinfection possible

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC
its affiliates or related companies, unless otherwise noted.
©2025 Hu-Friedy Mfg. Co., LLC. All rights reserved.

..:HuFriedyGroup

The Best In Practice

For How-To Videos, visit
www.HuFriedyGroup.com/IMSMaintenance

Hu-Friedy



Hu-Friedy Mfg. Co., LLC
2501 Sanders Rd.
Northbrook, IL 60062 | USA



EU | REP Hu-Friedy Europe LLC & Co. KG
Bogenstrasse 6 - 8
78576 Emmingen - Liptingen
Germany

CH | REP MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

UK | REP MDSS-UK RP Ltd.
Parkway House, Palatine Rd,
Northenden, Wythenshawe,
Manchester M22 4DB
United Kingdom



B2935 / REV 5 -2508

