

## IMS™ Needle Cap Holder

**REF** IM1003

Stainless Steel  
Acier Inoxydable

**QTY** 6

Clean and sterilize before use.  
Nettoyer et stériliser avant utilisation.

Recommended for use with cassettes only.

Est recommandé pour une utilisation avec des cassettes uniquement.



[www.hufriedygroup.com/Reprocessing/](http://www.hufriedygroup.com/Reprocessing/)  
1-800-Hu-Friedy; [info@hu-friedy.com](mailto:info@hu-friedy.com)



Sterilizable in steam sterilizer  
at 134 °C / 273 °F for 4 min



Washer-disinfector for thermal  
disinfection possible

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC  
its affiliates or related companies, unless otherwise noted.

©2025 Hu-Friedy Mfg. Co., LLC. All rights reserved.

**∴HuFriedyGroup**  
The Best In Practice

For How-To Videos, visit  
[www.HuFriedyGroup.com/IMSMaintenance](http://www.HuFriedyGroup.com/IMSMaintenance)

**Hu-Friedy**



**Hu-Friedy Mfg. Co., LLC**  
2501 Sanders Rd.  
Northbrook, IL 60062 | USA

B2944 / REV 9 -2508



Made in USA



**∴HuFriedyGroup**  
The Best In Practice