


# PRACTICE SMARTER WITH NEW ULTRASONIC INSERT KITS

*Safely and effectively manage your inserts*

ULTRASONIC  
SCALING


INSTRUMENT  
MANAGEMENT



Ultrasonic inserts remove even the most tenacious debris, but they can be delicate. The stack should be straight for optimal vibration, and tips need to keep their original shape for superior effectiveness.

Get protection for your inserts and peace of mind that your inserts will function optimally each and every time with **new Hu-Friedy Ultrasonic Kits featuring Swivel Direct Flow Inserts.**

- **Superior Protection**  
from weight of other instruments during cleaning and care process
- **Extended Insert Lifespan**  
from minimized potential damage to insert tip and stack
- **Safer Protocol** as less handling of inserts reduces chances for sharps injuries
- **Professional Appearance** on procedural tray
- **Color-coded** for at-a-glance identification



Left/Right/Straight Swivel Direct Flow Set, 30k  
+ Signature Series Cassette with Blue Rails

| IMINSRTSET

#10 Swivel Direct Flow Insert, 30K  
+ Signature Series Cassette with Lavender Rails

| IMINSRT10

#100 Swivel Direct Flow Insert, 30k  
+ Signature Series Cassette with Black Rails

| IMINSRT100

#1000 Triple Bend Swivel Direct Flow Insert, 30K  
+ Signature Series Cassette with Orange Rails

| IMINSRT1000

XT Triple Bend Swivel Direct Flow Insert, 30K  
+ Signature Series Cassette with Lavender Rails

| IMINSRTXTB

For more on Hu-Friedy's ultrasonic insert line,  
visit [Hu-Friedy.com/Inserts](https://www.hu-friedy.com/Inserts)

Call 1-800-Hu-Friedy and request HF-512, the Ultrasonic Insert Reference Guide.

How the best perform

