



SURGICAL



ATLAS FORCEPS: INNOVATION IN THE PALM OF YOUR HAND

Atlas Forceps are uniquely designed to enhance performance and comfort by providing strength, precision and flexibility throughout the surgical procedure.



How the best perform



ATLAS FORCEPS COLLECTION

It takes years of hard work, education and training to determine the proper balance of the firm yet delicate technique needed during an extraction procedure to ensure patient comfort. Engineered to provide simplicity and precision, Atlas Forceps have the power to effectively perform even the most difficult extractions. The innovative internal spring helps reduce hand manipulation, increases comfort and control, and is available in apical designs for an atraumatic extraction. Atlas Forceps from Hu-Friedy help you perform at your best—every time.

Internal springs maintain the forceps in the open position, **helping reduce manipulation during surgery, enhancing comfort and control**

Apical designs provide a **firm, secure grasp of the crown and root**

Thin, sharp tapered beaks for **ease in penetrating the PDL and apical progression***

Innovative spring hinge ensures forceps are in the **open position during reprocessing**

Laser marked hinge for **easier identification when bagged during cleaning and sterilization**


Meticulously hand crafted from the **finest grade surgical stainless steel**

*Apical designs with tapered beaks during subgingival apical progression in an extraction

Internal springs maintain open position, enhancing comfort and control

Improved beak geometries for a firm, secure grasp of the crown and root

Laser marked hinge for easy identification





| F150IS

150 Upper Universal, Atlas Forceps



| F151IS

151 Lower Universal, Atlas Forceps



| F11S

1 Upper Anterior, Atlas Forceps



| F23IS

23 Cowhorn, Lower Molars, Atlas Forceps



| F222IS

222 Lower Molars, Atlas Forceps



| F69IS

69 Tomes Upper/Lower Roots, Atlas Forceps



| F88LIS

88L Nevius Upper Molar, Left, Atlas Forceps



| F88RIS

88R Nevius Upper Molar, Right, Atlas Forceps



| F53LIS

53L Upper Molar, Left, Atlas Forceps



| F53RIS

53R Upper Molar, Right, Atlas Forceps



| FAF1IS

1 Upper Anterior, Apical, Atlas Forceps



| FAF150IS

150 Upper Universal, Apical, Atlas Forceps



| FAF151IS

151 Lower Universal, Apical, Atlas Forceps



| FAF222IS

222 Lower Molars, Apical, Atlas Forceps


PREVENTIVE




RESTORATIVE



ORTHODONTIC




INSTRUMENT
MANAGEMENT



ENDODONTIC



SURGICAL



POWER



FOR MORE INFORMATION

visit [Hu-Friedy.com](https://www.hu-friedy.com)
or call 1-800-Hu-Friedy (1-773-975-6100)

© 2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-293/0819

