

# ADA CDT Coding Collection Series



Infection Prevention and  
Instrument Management

Start pairing patients and insurance codes more effectively with procedural instrument kits that follow American Dental Association (ADA)—Current Dental Terminology (CDT) Coding.

Begin each preventative appointment with an Evaluation Kit and add the supplemental treatment kit based on the patient assessment. Each kit is uniquely color-coded with a compact cassette for ensured equipment compatibility and safety.



IM-D0120/80



## EVALUATION KIT

Exam 5 Instrument Cassette, Blue	<b>IM6058</b>
Mirror Handle, SE Satin Steel #6 Handle	<b>MH6</b>
5 Mirror Head, Double Sided High Definition	<b>MIR5DSHD</b>
#5 Explorer, DE Satin Steel #6 Handle	<b>EXD56</b>
UNC Probe, SE Satin Steel #6 Handle	<b>PCPUNC126</b>
#11/12 ODU DE After 5™ Explorer, DE Satin Steel #6 Handle	<b>EXD11/12A6</b>
Nabors Probe 3-6-9-12, DE Satin Steel #6 Handle	<b>PQ2N6</b>

IM-D1110



## PROPHYLAXIS KIT

Exam 5 Instrument Cassette, Green	<b>IM6059</b>
Plastic Test Stick — Individual	<b>PTS-IND</b>
Nevi™ 2 Sickle Scaler, DE EverEdge™ 2.0 #9 Handle	<b>SCNEVI29E2</b>
Nevi™ 4 Sickle Scaler, DE EverEdge™ 2.0 #9 Handle	<b>SCNEVI49E2</b>
5/6 Barnhart Universal Curette, DE EverEdge™ 2.0 #9 Handle	<b>SBH5/69E2</b>
15/16 Gracey Curette, DE EverEdge™ 2.0 #9 Handle	<b>SG15/169E2</b>
13/14 Gracey Curette, DE EverEdge™ 2.0 #9 Handle	<b>SG13/149E2</b>

IM-D4341/42



## SCALING AND ROOT PLANING KIT

Double Decker, 8 Instrument Cassette, Red	<b>IM9080</b>
Plastic Test Stick — Individual	<b>PTS-IND</b>
3/4 Langer After 5™ Universal Curette, DE EverEdge™ 2.0 #9 Handle	<b>SL3/4AF9E2</b>
Nevi™ 2 Sickle Scaler, DE EverEdge™ 2.0 #9 Handle	<b>SCNEVI29E2</b>
Nevi™ 4 Sickle Scaler, DE EverEdge™ 2.0 #9 Handle	<b>SCNEVI49E2</b>
1/2 Mini 5™ Gracey Curette, DE EverEdge™ 2.0 #9 Handle	<b>SAS1/29E2</b>
7/8 After 5™ Gracey Curette, DE EverEdge™ 2.0 #9 Handle	<b>SRPG7/89E2</b>
15/16 After 5™ Gracey Curette, DE EverEdge™ 2.0 #9 Handle	<b>SRPG15/169E2</b>
13/14 After 5™ Gracey Curette, DE EverEdge™ 2.0 #9 Handle	<b>SRPG13/149E2</b>
17/18 Gracey Curette, DE EverEdge™ 2.0 #9 Handle	<b>SG17/189E2</b>

IM-D1120



## PEDO PROPHY KIT

Exam 5 Instrument Cassette, Lavender	<b>IM6056</b>
Plastic Test Stick — Individual	<b>PTS-IND</b>
Mirror Handle, SE Satin Steel #6 Handle	<b>MH6</b>
5 Mirror Head, Double Sided High Definition	<b>MIR5DSHD</b>
#5 DE Explorer, #6 Handle	<b>EXD56</b>
13/14 Columbia Universal Curette, DE EverEdge™ 2.0 #9 Handle	<b>SC13/149E2</b>
Nevi™ 3 Sickle Scaler, DE EverEdge™ 2.0 #9 Handle	<b>SCNEVI39E2</b>
Nevi™ 1 Sickle Scaler, DE EverEdge™ 2.0 #9 Handle	<b>SCNEVI19E2</b>

# ADA CDT Coding Collection Series



Infection Prevention and Instrument Management

IM-D4346



## PRESENCE OF GINGIVAL INFLAMMATION KIT

Exam 5 Instrument Cassette, Yellow	IM6055
Plastic Test Stick — Individual	PTS-IND
Nevi™ 3 Sickle Scaler, DE EverEdge™ 2.0 #9 Handle	SCNEVI39E2
Nevi™ 4 Sickle Scaler, DE EverEdge™ 2.0 #9 Handle	SCNEVI49E2
1/2 Mini 5™ Gracey Curette, DE EverEdge™ 2.0 #9 Handle	SAS1/29E2
13/14 After 5™ Gracey Curette, DE EverEdge™ 2.0 #9 Handle	SRPG13/149E2
15/16 After 5™ Gracey Curette, DE EverEdge™ 2.0 #9 Handle	SRPG15/169E2

IM-D4355



## DEBRIDEMENT FOR COMPREHENSIVE EVALUATION KIT

Exam 5 Instruments Cassette, Orange	IM6053
Plastic Test Stick — Individual	PTS-IND
Nevi™ 4 Sickle Scaler, DE EverEdge™ 2.0 #9 Handle	SCNEVI49E2
17/18 McCall Universal Curette, DE EverEdge™ 2.0 #9 Handle	SM17/189E2
1/2 Rigid Mini 5™ Gracey Curette, DE EverEdge™ 2.0 #9 Handle	SAS1/2R9E2
13/14 Rigid After 5™ Gracey Curette, DE EverEdge™ 2.0 #9 Handle	SRP13/149RE2
15/16 Rigid After 5™ Gracey Curette, DE EverEdge™ 2.0 #9 Handle	SRP15/16R9E2

IM-D6080



## IMPLANT MAINTENANCE KIT

Exam 5 Instrument Cassette, White	IM6052
Plastic Test Stick — Individual	PTS-IND
Colorvue™ Probe, UNC 12, SE Satin Steel Handle	PCVPROMONC
Titanium Implant Scaler, 204SD Sickle Scaler	IMP204SDT
Titanium Implant Scaler, Langer 1/2 Universal Curette	IMPLG1/2T
Titanium Implant Scaler, 1/2 Gracey Mini Five™ Curette	IMPM1/2T
Implacare II H6/7 Scaler	IMPLH6/76

IM-D4910



## PERIODONTAL MAINTENANCE KIT

Double Decker™, 8 Instrument Cassette, Gray	IM9081
Plastic Test Stick — Individual	PTS-IND
3/4 Langer Universal Curette, DE EverEdge™ 2.0 #9 Handle	SL3/49E2
Nevi™ 2 Sickle Scaler, DE EverEdge™ 2.0 #9 Handle	SCNEVI29E2
Nevi™ 4 Sickle Scaler, DE EverEdge™ 2.0 #9 Handle	SCNEVI49E2
1/2 Pattison Gracey LITE Curette, DE EverEdge™ 2.0 #9 Handle	SG1/29APE2
7/8 Pattison Gracey LITE Curette, DE EverEdge™ 2.0 #9 Handle	SG7/89APE2
11/12 Pattison Gracey LITE Curette, DE EverEdge™ 2.0 #9 Handle	SG11/129APE2
13/14 Pattison Gracey LITE Curette, DE EverEdge™ 2.0 #9 Handle	SG13/149APE2
#15/16 Gracey Curette, EverEdge™ 2.0 #9 Handle	SG15/169E2

## DID YOU KNOW?

Anna Pattison Curettes are modified from the original Gracey Curettes. They have shorter, thinner, smaller blades honed at a 60 degree angle with a modified rigid shank specially designed to allow easier access and easier insertion under tight tissue.

Anna  
Pattison  
Curette



Gracey  
Curette