




STREAMLINE PLUS INSERTS

NEW FROM HU-FRIEDY

Enjoy comfort and enhanced ergonomics with new Streamline Plus ultrasonic inserts. These inserts feature the popular design of our Swivel™ line, loved by clinicians for their comfort and efficiency, but in a non-swivel option. Now you can choose what works best for you, your equipment, and your preferences.



STREAMLINE PLUS INSERTS OFFER:

- **Clear line of sight to treatment area** thanks to tapered cone with no bulk near working end
- **Better ergonomics** with large-diameter, textured-silicone grip for reduced pinch force and comfortable, secure grasp
- **Reduced excess spray** with direct flow water delivery
- **Easy at-a-glance tip identification** with color-coded grips
- **Advanced tip longevity** from premium, medical grade steel
- **Customized patient care** with selection of all popular tip designs



STREAMLINE SELECTIONS

Streamline Plus



Streamline Direct Flow



Streamline Classic



Color Coded for Easy Tip Identification

Grip Material

Silicone

Resin

Resin

Direct Flow Water Delivery Near Tip for Reduced Excess Spray

Premium Level Steel Alloy

Ergonomic, 12mm Wide Diameter Handle

Textured Soft Grip for a Comfortable, Secure Grasp

Tapered Cone for Enhanced Visibility of Treatment Site

#10 Universal (Lavender)

UI30SP10

UI30SF10*

UI1030K

#100 Thin (Black)

UI30SP100

UI30SF100*

UI10030K*

#1000 Triple Bend (Orange)

UI30SP1000

UI30SF1000*

UI100030K*

XT (Purple)

UI30SFXT

XT Triple Bend (Purple)

UI30SPXTTB

UI30SFXTTB

After Five Left (Teal)

UI30SPL

UI30SFL*

After Five Straight (Blue)

UI30SPS

UI30SFS*

After Five Right (Red)

UI30SPR

UI30SFR*

Beavertail (Yellow)

UI330K*

For more information, visit [Hu-Friedy.com/Inserts](https://www.hu-friedy.com/Inserts)

Call 1-800-Hu-Friedy or contact your authorized Hu-Friedy representative for more information.

© 2020 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-294/0919

*Visit Hu-Friedy for 25kHz equivalent part codes

How the best perform

