



How the best perform

Shine reNEW Instrument Wipes

Safety Data Sheet

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022
Issue date: 2026-08-06 Revision date: 2026-08-06 Version: 2.0

SECTION 1 Identification

1.1. Product identifier

Product form : Mixture
Product name : Shine reNEW Instrument Wipes
Product code : IMS-1455

1.2. Other means of identification

No additional information available

1.3. Recommended use of the chemical and restrictions on use

Use of the substance/mixture : Cleaner

1.4. Supplier's details

Manufacturer

Hu-Friedy Mfg. Co., LLC
2501 Sanders Rd.
Northbrook, IL 60062 USA
T 1-800-483-7433

Distributor

Acmedent Corp.
91 Citation Dr, Unit 5, Concord L4K 2Y8
905-761-6850
Central Dental Supply
103-251 Saulteaux Crescent, Winnipeg R3J 3C7
204-885-7530
Henry Schein Canada, Inc.
345 Townline Road, SS4, Niagara-On-The-Lake L0S 1J0
905-646-1711
K-Dental Inc.
750 Cochrane Drive, Markham L3R 8E1
416-293-8365
Patterson Dentaire Canada
1205, Boul.henri-Bourassa West, Montreal H3M 3E6
418-688-6546
Sinclair Dental Company Ltd.
900 Harbourside Drive, North Vancouver V7P 3T8
604-986-1544

1.5. Emergency phone number

Emergency number : 1-800-535-5053

SECTION 2 Hazard identification

2.1. Classification of the substance or mixture

GHS classification

Serious eye damage, Category 1
Skin sensitization, Category 1A

2.2. Label elements

GHS labelling

Hazard pictograms (GHS)



Signal word (GHS) : Danger

Shine reNEW Instrument Wipes

Safety Data Sheet

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022

Hazard statements (GHS)	: May cause an allergic skin reaction Causes serious eye damage
Precautionary statements (GHS)	: Avoid breathing dust, fume, gas, mist, vapours, spray. Contaminated work clothing must not be allowed out of the workplace. Wear protective gloves, protective clothing, eye protection, face protection, and hearing protection. If on skin: Wash with plenty of water. Take off contaminated clothing and wash it before reuse. If skin irritation or rash occurs: Get medical advice or attention. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a poison center or doctor. Dispose of contents and/or container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulations.

2.3. Hazards associated with known or reasonably anticipated uses

No additional information available

2.4. Hazards not otherwise classified

No additional information available

2.5. Unknown acute toxicity

Not applicable

SECTION 3 Composition/information on ingredients

3.1. Substances

Not applicable

Shine reNEW Instrument Wipes

Safety Data Sheet

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022

3.2. Mixtures

Name	Chemical name / Synonyms	Product identifier	Conc. (% w/w)
Alcohols, C12-14-secondary, ethoxylated	Alcohols, C12-14-secondary, ethoxylated Alcohols, C12-14, secondary, ethoxylated / C12-14 sec-Pareth-5 / Softanol 70 / C12-14 SEC-PARETH-3 / Softanol-70 / C12-14 sec-Pareth-9 / C12-14 sec-Pareth-8 / C12-14 sec-Pareth-7 / C12-14 sec-Pareth-50 / C12-14 sec-Pareth-40 / C12-14 sec-Pareth-30 / C12-14 sec-Pareth-3 / C12-14 sec-Pareth-20 / Secondary alcohols, C12-14, ethoxylated / C12-14 SEC-PARETH-12 / C12-14 SEC-PARETH-15 / C12-14 SEC-PARETH-20 / C12-14 SEC-PARETH-30 / C12-14 SEC-PARETH-40 / C12-14 SEC-PARETH-5 / C12-14 SEC-PARETH-50 / C12-14 SEC-PARETH-7 / C12-14 SEC-PARETH-8 / C12-14 SEC-PARETH-9 / Ethoxylated secondary alcohols(C12-14) / c12-14 sec-pareth-3 / C12-14 sec-Pareth / C12-14 SEC-ALKETH-12	CAS-No.: 84133-50-6	1 – 5
Acetic acid, chloro-, reaction products with 2-heptadecyl-4,5-dihydro-1H-imidazole-1-ethanol and sodium hydroxide	Acetic acid, chloro-, reaction products with 2-heptadecyl-4,5-dihydro-1H-imidazole-1-ethanol and sodium hydroxide Acetic acid, Chloro-, reaction products with 2-heptadecyl-4,5-dihydro-1H-imidazole-1-ethanol and sodium hydroxide / Acetic acid, 2-chloro-, reaction products with 2-heptadecyl-4,5-dihydro-1H-imidazole-1-ethanol and sodium hydroxide / Stearic acid, aminoethylethanolamine amide-imidazoline, carboxymethylated sodium salt	CAS-No.: 68608-63-9	0.1 – 1

Shine reNEW Instrument Wipes

Safety Data Sheet

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022

Name	Chemical name / Synonyms	Product identifier	Conc. (% w/w)
3(2H)-Isothiazolone, 2-methyl-	3(2H)-Isothiazolone, 2-methyl- 2-Methyl-2,3-dihydroisothiazol-3-one / 2-Methylisothiazol-3(2H)-one / 3(2H)-Isothiazolon-3-one, 2-methyl- / 2-Methylisothiazolin-3(2H)-one / N-Methyl-isothiazolone / methylisothiazolinone / MIT / METHYLISOTHIAZOLINONE / Methyl-4-isothiazolin-3-one, 2- / Methylisothiazolone / Methylisothiazolinone / 2-Methyl-4-isothiazolone-3-one / 2-Methyl-4-isothiazolin-3-one / 2-Methyl-2H-isothiazol-3-one / 3-Isothiazolone, 2-methyl- / 2-Methyl-3-isothiazolone	CAS-No.: 2682-20-4	0.01 – < 0.1
3(2H)-Isothiazolone, 5-chloro-2-methyl-	3(2H)-Isothiazolone, 5-chloro-2-methyl- 5-Chloro-2-methyl-isothiazolone-3(2H)-one / 2-Methyl-5-chloro-2H-isothiazol-3-one / 3(2H)-Isothiazolon-3-one, 5-chloro-2-methyl- / CIT / 5-Chloro-2-methyl-isothiazolin-3(2H)-one / 5-Chloro-2-methyl-4-thiazoline-3-ketone / 5-Chloro-2-methylisothiazol-3(2H)-one / 2H-Isothiazol-3-one, 5-chloro-2-methyl- / 5-Chloro-2-methylisothiazolone / methylchloroisothiazolinone / 2-Methyl-5-chloroisothiazolin-3-one / 5-Chloro-2-methyl-3(2H)-isothiazolone / METHYLCHLOROISOTHIAZOLINONE / Methylchloroisothiazolinone / 4-Isothiazolin-3-one, 5-chloro-2-methyl- / Isothiazol(2H)-3-one, 5-chloro-2-methyl- / 5-Chloro-2-methyl-4-isothiazolin-3-one / 5-Chloro-2-methyl-2H-isothiazol-3-one / 5-Chloro-2-methyl-3-isothiazolone	CAS-No.: 26172-55-4	0.01 – < 0.1

*Chemical name, CAS number and/or exact concentration have been withheld as a trade secret

Specific concentration limits:		
Name	Product identifier	Specific concentration limits (%)
3(2H)-Isothiazolone, 2-methyl-	CAS-No.: 2682-20-4 EC-No.: 220-239-6 EC Index-No.: 613-326-00-9	(0.0015 ≤ C < 100) Skin Sens. 1A; H317
3(2H)-Isothiazolone, 5-chloro-2-methyl-	CAS-No.: 26172-55-4 EC-No.: 247-500-7	(0.0015 ≤ C ≤ 100) Skin Sens. 1; H317 (0.06 ≤ C < 0.6) Eye Irrit. 2; H319 (0.06 ≤ C < 0.6) Skin Irrit. 2; H315 (0.6 ≤ C ≤ 100) Skin Corr. 1B; H314

Shine reNEW Instrument Wipes

Safety Data Sheet

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022

SECTION 4 First-aid measures

4.1. Description of necessary first-aid measures

First-aid measures after inhalation	: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing. Get medical advice/attention if you feel unwell.
First-aid measures after skin contact	: IF ON SKIN: Wash with plenty of Water. Take off contaminated clothing and wash it before reuse. If skin irritation or rash occurs: Get medical advice/attention.
First-aid measures after eye contact	: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor.
First-aid measures after ingestion	: Do not induce vomiting without medical advice. Never give anything by mouth to an unconscious person. Get medical advice/attention if you feel unwell.

4.2. Most important symptoms/effects, acute and delayed

Symptoms/effects after inhalation	: May cause irritation to the respiratory tract.
Symptoms/effects after skin contact	: May cause skin irritation. Repeated exposure may cause skin dryness or cracking. May cause an allergic skin reaction.
Symptoms/effects after eye contact	: Causes serious eye damage. Symptoms may include discomfort or pain, excess blinking and tear production, with marked redness and swelling of the conjunctiva. May cause burns.
Symptoms/effects after ingestion	: May be harmful if swallowed. May cause gastrointestinal irritation, nausea, vomiting and diarrhea.

4.3. Indication of immediate medical attention and special treatment needed, if necessary

Other medical advice or treatment	: Symptoms may be delayed. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
-----------------------------------	--

SECTION 5: Fire-fighting measures

5.1. Suitable (and unsuitable) extinguishing media

Suitable extinguishing media	: Use extinguishing media appropriate for surrounding fire.
Unsuitable extinguishing media	: Do not use water jet.

5.2. Specific hazards arising from the chemical

Fire hazard	: Products of combustion may include, and are not limited to: oxides of carbon.
-------------	---

5.3. Special protective equipment and precautions for fire-fighters

Protection during firefighting	: Keep upwind of fire. Wear full fire fighting turn-out gear (full Bunker gear) and respiratory protection (SCBA).
--------------------------------	--

SECTION 6 Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures	: Use personal protection recommended in Section 8. Isolate the hazard area and deny entry to unnecessary and unprotected personnel.
------------------	--

For non-emergency personnel

No additional information available

For emergency responders

Environmental precautions	: Prevent entry to sewers and public waters.
---------------------------	--

Shine reNEW Instrument Wipes

Safety Data Sheet

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022

6.2. Methods and materials for containment and cleaning up

- For containment : Pick up large pieces, then place in a suitable container. Do not flush to sewer or allow to enter waterways. Use appropriate Personal Protective Equipment (PPE).
- Methods for cleaning up : Sweep or shovel spills into appropriate container for disposal. Provide ventilation.

For further information refer to section 8: "Exposure controls/personal protection"

SECTION 7 Handling and storage

7.1. Precautions for safe handling

- Precautions for safe handling : Avoid contact with skin. Do not get in eyes. Avoid breathing dust/fume/gas/mist/vapours/spray. Do not swallow. Handle and open container with care. When using do not eat, drink or smoke.
- Hygiene measures : Take off contaminated clothing and wash it before reuse. Contaminated work clothing should not be allowed out of the workplace. Wash hands, forearms and face thoroughly after handling.

7.2. Conditions for safe storage, including incompatibilities

- Storage conditions : Keep out of the reach of children. Store tightly closed in a dry, cool and well-ventilated place.
- Specific end uses : Cleaner.

SECTION 8 Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Appropriate engineering controls

- Appropriate engineering controls : Ensure good ventilation of the work station. Provide readily accessible eye wash stations and safety showers.
- Environmental exposure controls : Avoid release to the environment.

8.3. Individual protection measures, such as personal protective equipment

Hand protection:
Wear suitable gloves resistant to chemical penetration. Consult glove manufacturer's product information on material suitability and material thickness.
Eye protection:
Wear eye protection
Skin and body protection:
Wear suitable protective clothing
Respiratory protection:
In case of insufficient ventilation, wear suitable respiratory equipment. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

Other information:

Handle in accordance with good industrial hygiene and safety procedures. Do not eat, drink or smoke when using this product.

Shine reNEW Instrument Wipes

Safety Data Sheet

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022

SECTION 9 Physical and chemical properties

9.1. Basic physical and chemical properties

Physical state	: Liquid
Appearance	: Wipes saturated with liquid.
Colour	: White
Odour	: Floral
Odour threshold	: No data available
pH	: 5.61
Melting point	: No data available
Freezing point	: No data available
Boiling point	: No data available
Flash point	: > 200 °F / 93.3 °C
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20°C/ 68 °F	: No data available
Relative density	: 1.0593
Solubility	: No data available
Partition coefficient n-octanol/water	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Viscosity, kinematic	: No data available
Explosive limits	: No data available
Particle characteristics	: No data available

3(2H)-Isothiazolone, 2-methyl-

Boiling point	> 130 °C Atm. press.: 16 hPa Decomposition: 'yes' Decomp. temp.: 130 °C
Vapour pressure	0.99 Pa (at 20 °C)
Particle characteristics	No data available

3(2H)-Isothiazolone, 5-chloro-2-methyl-

Boiling point	106.5 °C (at 1000 hPa)
Flash point	> 100 °C (closed cup)
Vapour pressure	20.8 hPa (at 20 °C)
Particle characteristics	No data available

9.2. Data relevant with regard to physical hazard classes (supplemental)

No additional information available

SECTION 10 Stability and reactivity

10.1. Reactivity

No dangerous reactions known under normal conditions of use.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

Shine reNEW Instrument Wipes

Safety Data Sheet

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022

10.4. Conditions to avoid

Heat.

10.5. Incompatible materials

None known.

10.6. Hazardous decomposition products

May include, and are not limited to: oxides of carbon.

SECTION 11 Toxicological information

11.1. Information on toxicological effects

Acute toxicity (oral) : Not classified.
Acute toxicity (dermal) : Not classified.
Acute toxicity (inhalation) : Not classified.

Alcohols, C12-14-secondary, ethoxylated (84133-50-6)

LD50 oral rat	2100 mg/kg (Source: NZ_CCID)
---------------	------------------------------

3(2H)-Isothiazolone, 2-methyl- (2682-20-4)

LD50 oral rat	120 mg/kg (Source: EU_CLH)
LD50 dermal rabbit	200 mg/kg (Source: NLM_HSDB)
LC50 inhalation rat	0.11 mg/l/4h

3(2H)-Isothiazolone, 5-chloro-2-methyl- (26172-55-4)

LD50 oral rat	481 mg/kg (Source: IUCLID)
LC50 inhalation rat	1.23 mg/l/4h

Skin corrosion/irritation : Not classified.
pH: 5.61

3(2H)-Isothiazolone, 2-methyl- (2682-20-4)

pH	2.58 Temp.: 25 °C Concentration: 50 g/L
----	---

Serious eye damage/irritation : Causes serious eye damage.
pH: 5.61

3(2H)-Isothiazolone, 2-methyl- (2682-20-4)

pH	2.58 Temp.: 25 °C Concentration: 50 g/L
----	---

Respiratory or skin sensitisation : May cause an allergic skin reaction.
Germ cell mutagenicity : Not classified.
Carcinogenicity : Not classified.
Reproductive toxicity : Not classified.
STOT-single exposure : Not classified.
STOT-repeated exposure : Not classified.

3(2H)-Isothiazolone, 2-methyl- (2682-20-4)

LOAEL (oral, rat, 90 days)	71.2 mg/kg bodyweight Animal: rat, Guideline: OECD Guideline 407 (Repeated Dose 28-Day Oral Toxicity Study in Rodents), Guideline: other:
----------------------------	---

Aspiration hazard : Not classified.

Shine reNEW Instrument Wipes

Safety Data Sheet

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022

Shine reNEW Instrument Wipes	
Viscosity, kinematic	No data available
Alcohols, C12-14-secondary, ethoxylated (84133-50-6)	
Viscosity, kinematic	No data available

Acetic acid, chloro-, reaction products with 2-heptadecyl-4,5-dihydro-1H-imidazole-1-ethanol and sodium hydroxide (68608-63-9)	
Viscosity, kinematic	No data available

3(2H)-Isothiazolone, 2-methyl- (2682-20-4)	
Viscosity, kinematic	No data available

3(2H)-Isothiazolone, 5-chloro-2-methyl- (26172-55-4)	
Viscosity, kinematic	No data available
Symptoms/effects after inhalation	: May cause irritation to the respiratory tract.
Symptoms/effects after skin contact	: May cause skin irritation. Repeated exposure may cause skin dryness or cracking. May cause an allergic skin reaction.
Symptoms/effects after eye contact	: Causes serious eye damage. Symptoms may include discomfort or pain, excess blinking and tear production, with marked redness and swelling of the conjunctiva. May cause burns.
Symptoms/effects after ingestion	: May be harmful if swallowed. May cause gastrointestinal irritation, nausea, vomiting and diarrhea.
Other information	: Likely routes of exposure: ingestion, inhalation, skin and eye.

SECTION 12 Ecological information

12.1. Ecotoxicity

Ecology - general	: May cause long-term adverse effects in the aquatic environment.
Hazardous to the aquatic environment, short-term (acute)	: Not classified.
Hazardous to the aquatic environment, long-term (chronic)	: Not classified.

Alcohols, C12-14-secondary, ethoxylated (84133-50-6)	
LC50 - Fish [1]	3.2 mg/l (Exposure time: 96 h - Species: Pimephales promelas)
EC50 - Crustacea [1]	3.2 mg/l (Exposure time: 48 h - Species: water flea)

3(2H)-Isothiazolone, 2-methyl- (2682-20-4)	
LC50 - Fish [1]	4.77 mg/l Test organisms (species): Oncorhynchus mykiss (previous name: Salmo gairdneri)
EC50 - Crustacea [1]	1.6 mg/l Test organisms (species): Daphnia magna

3(2H)-Isothiazolone, 5-chloro-2-methyl- (26172-55-4)	
LC50 - Fish [1]	1.6 mg/l (Exposure time: 96 h - Species: Oncorhynchus mykiss [semi-static])
EC50 - Crustacea [1]	4.71 mg/l (Exposure time: 48 h - Species: Daphnia magna)
EC50 - Crustacea [2]	0.12 – 0.3 mg/l (Exposure time: 48 h - Species: Daphnia magna [Flow through])
EC50 72h - Algae [1]	0.11 – 0.16 mg/l (Species: Pseudokirchneriella subcapitata [static])
EC50 96h - Algae [1]	0.03 – 0.13 mg/l (Species: Pseudokirchneriella subcapitata [static])

Shine reNEW Instrument Wipes

Safety Data Sheet

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022

12.2. Persistence and degradability

Shine reNEW Instrument Wipes

Persistence and degradability	Not established.
-------------------------------	------------------

Alcohols, C12-14-secondary, ethoxylated (84133-50-6)

Persistence and degradability	Rapidly degradable
-------------------------------	--------------------

Acetic acid, chloro-, reaction products with 2-heptadecyl-4,5-dihydro-1H-imidazole-1-ethanol and sodium hydroxide (68608-63-9)

Persistence and degradability	Rapidly degradable
-------------------------------	--------------------

3(2H)-Isothiazolone, 2-methyl- (2682-20-4)

Persistence and degradability	Rapidly degradable
-------------------------------	--------------------

3(2H)-Isothiazolone, 5-chloro-2-methyl- (26172-55-4)

Persistence and degradability	Rapidly degradable
-------------------------------	--------------------

12.3. Bioaccumulative potential

Shine reNEW Instrument Wipes

Bioaccumulative potential	Not established.
---------------------------	------------------

3(2H)-Isothiazolone, 2-methyl- (2682-20-4)

Partition coefficient n-octanol/water	-0.26 (at 20 °C (at pH 5)
---------------------------------------	---------------------------

3(2H)-Isothiazolone, 5-chloro-2-methyl- (26172-55-4)

Partition coefficient n-octanol/water	-0.71 – 0.75 (at 20 °C)
---------------------------------------	-------------------------

12.4. Mobility in soil

No additional information available

12.5. Other adverse effects

Ozone	: Not classified.
Fluorinated greenhouse gases	: No
Other information	: No other effects known.

SECTION 13 Disposal considerations

Product/Packaging disposal recommendations	: Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.
Additional information	: The generation of waste should be avoided or minimized wherever possible.

SECTION 14 Transport information

In accordance with DOT / TDG

14.1. UN Number

UN-No. (DOT)	: Not regulated
UN-No. (TDG)	: Not regulated

Shine reNEW Instrument Wipes

Safety Data Sheet

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022

14.2. UN Proper Shipping Name

Proper Shipping Name (DOT) : Not regulated
Proper Shipping Name (TDG) : Not regulated

14.3. Transport hazard class(es)

DOT
Transport hazard class(es) (DOT) : Not regulated

TDG
Transport hazard class(es) (TDG) : Not regulated

14.4. Packing group

Packing group (DOT) : Not regulated
Packing group (TDG) : Not regulated

14.5. Environmental hazards

Other information : No supplementary information available.

14.6. Transport in bulk

Not applicable

14.7. Special precautions for user

DOT
Not regulated

TDG
Not regulated

SECTION 15 Regulatory information

15.1. Federal regulations

All components of this product are listed, or excluded from listing, on the United States Environmental Protection Agency Toxic Substances Control Act (TSCA) inventory.

All components of this product are listed, or excluded from listing, on the Canadian DSL (Domestic Substances List) and NDSL (Non-Domestic Substances List) inventories.

15.2. International regulations

No additional information available

15.3. State regulations

California Proposition 65 - This product does not contain any substances known to the state of California to cause cancer, developmental and/or reproductive harm

SECTION 16 Other Information

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022

Revision date : 2026-08-06
Issue date : 2026-08-06
Other information : None.

Shine reNEW Instrument Wipes

Safety Data Sheet

According to the Hazard Communication Standard (CFR29 1910.1200) HazCom 2024 and the Hazardous Products Regulations (HPR) WHMIS 2022

Prepared by : Nexreg Compliance Inc.
www.Nexreg.com



Indication of changes:
SDS update.

SDS HazCom 2024 - WHMIS 2022 (Nexreg) 2025

Disclaimer: We believe the statements, technical information and recommendations contained herein are reliable, but they are given without warranty or guarantee of any kind. The information contained in this document applies to this specific material as supplied. It may not be valid for this material if it is used in combination with any other materials. It is the user's responsibility to satisfy oneself as to the suitability and completeness of this information for the user's own particular use.