



Black Line

Surgical Instruments

Designed for Performance



“Hu-Friedy’s surgical instruments are simply the best. They have been an innovator in instrument design and have maintained a reputation for providing the highest quality instruments for discriminating dentists for decades.”

— DR. EDWARD P. ALLEN



“SwissPerio Black Line instruments are the only ones that completely fulfill my expectations to visual perception and sensation during periodontal surgical interventions due to better vision and spatial orientation as well as improved cognition in periodontal and implant surgery.”

— DR. RINO BURKHARDT



“The Black Line’s dark surface reduces glare from my headlight and the working edge is easy to visualize which helps in maneuvering in a tortuous socket. The edge remains nice and sharp a long time consistent with Hu-Friedy quality.”

— DR. PETER HAUPERS, JR.



“I am so pleased with the new Black Line, it represents a coming together of smart instrument usage and the newest technology for a “home run” in surgical armamentarium.”

— DR. KARL KOERNER



“During delicate instrumentation I was irritated by the reflection of light on instruments and with the Black Line instruments it solves this problem and I am able to perform without compromises. Not only is the Black Line a good-looking instrument, the coating is also harder resulting in sharper longer lasting edges.”

— DR. YATARO KOMIYAMA, DDS, PHD



“Quality of life is achieved balancing your abilities with your target setting and efforts. Choosing the right instruments is always a good idea. Using the Black Line instruments of Hu-Friedy is a step further.”

— DR. STARVOS PELEKANOS



“I love to work with the Black Line instruments! I like the design, the precision and the balance in my hands. Since I introduced them in my everyday practice the work flows much smoother due to reduced light reflection. I like above all the visibility while using this instruments during the tunneling procedures.”

— PROF. GIULIO RASPERINI



“I was amazed to find the cutting edge of the Black Line instruments maintain their acuity throughout the entire periodontal surgery.”

— DR. JONSUZUKI



“I love the Black Line instruments. In addition to their perfect quality, they are supercool.”

— DR. ISTVAN URBAN


Introducing the Hu-Friedy

BLACK LINE SCISSORS COLLECTION


Designed for performance, the Hu-Friedy Black Line is engineered to optimize clinical outcomes by delivering efficiency throughout the entire perio and surgical procedure. From reducing light reflection to enhancing visibility at the surgical site, the Black Line was designed with the clinician in mind. In addition to the surgical hand instrument portfolio, the Black Line now includes 7 Super-Cut scissor designs delivering superior performance through enhanced precision and innovation.

Features and Benefits of Black Line Scissors:

- Super-Cut design featuring a **razor sharp blade and serrated edge for reduced tissue slippage**
- Performance engineered coating for **optimal edge retention**
- Distinct black finish for **enriched contrast and visual acuity at the surgical site and underlying tissue**
- **Reduced light reflection** afforded by a matte black finish
- Gold screw for easy identification of **Super-Cut design**



Black Line Iris Curved Scissors cutting PGA violet suture material after closing a simple loop suture.




SUPER-CUT
DESIGN FEATURING
A RAZOR SHARP BLADE
AND SERRATED EDGE FOR
REDUCED TISSUE SLIPPAGE

GOLD SCREW FOR
EASY SUPER-CUT
IDENTIFICATION

PERFORMANCE
ENGINEERED COATING
FOR REDUCED LIGHT
REFLECTION AND
ENHANCED CONTRAST





BLACK LINE SCISSORS

Black Line LaGrange Scissors cutting PGA Violet suture material after closing a free gingival tissue graft with a crisscross suture technique.



LaGrange Double-Curved Scissors, Super-Cut
S14SCX



Castro Microsurgical,
Super-Cut
SPVX



Goldman-Fox,
Curved,
Super-Cut
S16SCX



Goldman-Fox,
Straight,
Super-Cut
S20SCX



Iris,
Straight/Delicate,
Super-Cut
S17SCX



Iris,
Curved/Delicate,
Super-Cut
S18SCX



Metzenbaum,
Curved,
Super-Cut
S19SCX

BLACK LINE INSTRUMENT COLLECTION


Features and Benefits of the Black Line:

- Performance engineered coating for a harder, smoother surface for **optimal edge retention and enhanced lubricity**
- Distinct black finish for **enriched contrast and visual acuity at the surgical site and underlying tissue**
- **Reduced light reflection** afforded by a matte finished handle and black working ends
- Unique, smooth, large diameter lightweight handle for **maximum comfort, reduced hand fatigue and increased control**
- Handcrafted with Immunity Steel® for **optimal strength and corrosion resistance**

DISTINCT BLACK FINISH FOR
REDUCED LIGHT REFLECTION
AND ENHANCED CONTRAST

PERFORMANCE
ENGINEERED COATING
FOR ENHANCED
LUBRICITY AND
VISUAL ACUITY AT THE
UNDERLYING TISSUE

SMOOTH, LARGE DIAMETER
HANDLE FOR MAXIMUM COMFORT
AND REDUCED HAND FATIGUE



PERIOSTEALS

*A full thickness mucoperiosteal pouch
raised by blunt subperiosteal dissection*




Allen Periosteal Elevator, Anterior

PPAELAX


Short Handle

PPAELASX



3 Prichard Periosteal


PPR3X



9 Molt Periosteal

Elevator

P9X




Allen Periosteal Elevator

PPAELX

Short Handle


PPAELSX



Allen Periosteal Elevator,

Posterior

PPAELPX



Buser Periosteal

PPBUSERX

Short Handle

PPBUSERSX



NEW

Modified
Labanca
Periosteal


PPLABANC3X



NEW


24G
Periosteal, Short
Handle

P24GSPX



14
Goldman-Fox
Periosteal

P14X



152
Kramer-Nevins
Periosteal


PKN152X
Short Handle
PKN152SX



NEW

Buser/Modified
PH2MBH
Large

PPBUH2MX



NEW

Buser/Modified
PH2MBH
Small

PPBUH2M2X



NEW

Swiss Perio
Modified P24GSP
Periosteal


SPP24GSPX



NEW

Swiss Perio
Modified P24GSIC
Periosteal


SPP24GSICX



NEW

Swiss Perio
Modified P24GSP
Periosteal 2


SPP24GSP2X



NEW

Swiss Perio
Modified Schlee
Periosteal

SPPSCHLEEX



NEW

Swiss Perio
Modified PH2
Periosteal

SPPH2X



NEW

Urban Mini Me
Periosteal

PFIWDS1MKX



LUXATING HYBRIDS

A unique cross between a periosteal elevator and a luxating elevator, Luxating Hybrids are used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



NEW

Luxating Hybrid Elevator, Straight Spade

EPTSSPX



NEW

Luxating Hybrids,
Small Straight

EPTSMSX



NEW

Luxating Hybrids,
Small Curved

EPTSMCX



NEW

Luxating Hybrids,
Small Distal

EPTSMDX



NEW

Luxating Hybrids,
Small Mesial

EPTSMMX



NEW

Luxating Hybrids,
Large Distal

EPTLGDX



Luxating Hybrids,
Large Mesial

EPTLGMX

SYNDESMOTOMES



6 Right Bernard
Syndesmotomes

BER6X



NEW

7 Left Bernard
Syndesmotomes

BER7X



NEW

11 Straight Bernard
Syndesmotomes

BER11X

TUNNELING

Used for cutting and recontouring soft tissue during periodontal surgery.




Tunneling Knife 1

TKN1X

Short Handle

TKN1SX




Tunneling Knife 2

TKN2X

Short Handle

TKN2SX



GDQ Giles de Quincey 1
Tunneling Instrument

FQUINCEY1X



GDQ Giles de Quincey 2
Tunneling Instrument

FQUINCEY2X

Short handle instrument designs are made to securely fit within a DIN cassette.

PERIODONTAL KNIVES




Incision made within the sulcus to sever the attachment from the base of the sulcus to the osseous crest.




Allen End-Cutting Intrасulcular Knife

KPAX




1/2
Allen Orban
Knife
KO12KP03AX




1/2
Allen Orban
Knife, Round
KO12KP03RX




Allen
Arrowhead
Knife
KO12KP3R4X



15/16
Kirkland Knife
KK15/16X
Short Handle
KK15/16SX



1/2 Orban Knife
KO1/2X



11
Goldman-Fox
Knife
KGF11X


PERIODONTAL BONE SCRAPER

Bone Scraper recontouring the alveolar crest prior to implant placement.



1/2 Buser Modified Periodontal Scraper 3mm/4mm

PPBUSER12X



2/3 Buser Modified Periodontal Scraper 4mm/5mm

PPBUSER23X



1/3 Buser Modified Periodontal Scraper 3mm/5mm

PPBUSER13X



13K/TG Periodontal Chisel 3.3mm

S13K/TGX



13K/13KL Periodontal Chisel 3.3mm/5mm

S13K/13KLX



Palacci Chisel 4mm/5mm

SPAL13KX




SINUS LIFT INSTRUMENTS


*Kramer-Nevins Sinus Lift instrument
delicately separating the Schneiderian
Membrane from the sinus floor.*




Kramer-Nevins, Obtuse
IMP6578SX




Kramer-Nevins, Acute
IMP6577SCX



Sinus Curette 1
SINC1X
Short Handle
SINC1SX



Sinus Curette 2
SINC2X
Short Handle
SINC2SX




300
IMP300X

SURGICAL CURETTES

Molt curette debriding the osteotomy



2/4 Molt Curette, Double-Ended
CM2/4X



84 Lucas
Surgical Curette,
Double-Ended
2mm
CL84X



85 Lucas
Surgical Curette,
Double-Ended
2.5mm
CL85X



86 Lucas
Surgical Curette,
Double-Ended
2.8mm
CL86X



87 Lucas
Surgical Curette,
Double-Ended
3.5mm
CL87X



88 Lucas
Surgical Curette,
Double-Ended
4.7mm
CL88X

ROOT TIP PICK



13/14 Heidbrink Root Tip Pick
EHB13/14X

NEW



IMPLANTOLOGY

*Rompen Mesial/Distal Plugger
condensing bone between an implant.*



At immediate implant placement after extraction, filling the buccal gap with a biomaterial can significantly improve the long-term esthetic success. The special concave design of these pluggers highly facilitates and improves the precise placing of the chosen material.



– PROF. ERIC ROMPEN



.75mm




.75mm

NEW

Rompen Plugger, Mesial/Distal
PLGROMP2X


BONE GRAFT CONDENSERS / PLUGGERS

3.3mm




Grisdale Bone Condenser
PLGGR1X

3 x 1.5mm



Labanca Plugger, 3 x 1.5mm
PLGLABANCX

3mm



Jovanovic Plugger
PLGJO/4X


.75mm



Rompen Plugger,
Buccal/Lingual
PLGROMP1X

NEW

MEMBRANE REPLACEMENT



Membrane Placement Instrument
PMPIX

SPATULA



Palti Bone Graft Scoop 10x25mm 5x21mm
IMP6522X

NEW


PERIOTOMES

Periodontal ligament is severed from the root surface of the tooth.



Periotome, Anterior, Double-Ended

PT2X



Periotome, Posterior, Double-Ended
PT1X



Periotome, Right, Single-Ended
PT3X




Periotome, Left, Single-Ended
PT4X




Periotome, Anterior Angled, Single-Ended
PT5X



Periotome, Anterior Straight, Single-Ended
PT6X



Wise Periotome
PEREURX



LUXATING ELEVATORS

Luxating elevator inserted into the periodontal ligament space to sever the ligament and luxate the tooth.



Luxating Elevator, 3mm Straight

EL3SX



Luxating Elevator, 2mm
Straight

EL2SX



Luxating Elevator, 2mm
Curved

EL2CX



Luxating Elevator, 3mm
Curved

EL3CX



Luxating Elevator, 4mm
Straight

EL4SX



Luxating Elevator, 4mm
Curved

EL4CX



Luxating Elevator, 5mm
Straight

EL5SX




Luxating Elevator, 5mm
Curved

EL5CX

RETRACTOR

Minnesota Retractor used to reflect to mucoperiosteal during osteotomy debridement.



NEW

Minnesota Retractor
CRMX

BLACK LINE PART CODES

BLACK LINE SCISSORS

Castro Microsurgical, Super-Cut	SPVX
LaGrange Double-Curved, Super-Cut	S14SCX
Goldman-Fox, Curved, Super-Cut	S16SCX
Goldman-Fox, Straight, Super-Cut	S20SCX
Iris, Straight/Delicate, Super-Cut	S17SCX
Iris, Curved/Delicate, Super-Cut	S18SCX
Metzenbaum, Curved, Super-Cut	S19SCX

PERIOSTEALS

3 Prichard Periosteal	PPR3X
9 Molt Periosteal Elevator	P9X
Allen Periosteal Elevator, Anterior	PPAELAX
Allen Periosteal Elevator, Anterior, Short Handle	PPAELASX
Allen Periosteal Elevator	PPAELX
Allen Periosteal Elevator, Short Handle	PPAELSX
Allen Periosteal Elevator, Posterior	PPAELPX
Buser Periosteal	PPBUSERX
Buser Periosteal, Short Handle	PPBUSERSX
152 Kramer-Nevins Periosteal	PKN152X
152 Kramer-Nevins Periosteal, Short Handle	PKN152SX
Modified Labanca Periosteal	PPLABANC3X
24G Periosteal, Short Handle	P24GSPX
Goldman-Fox Periosteal	P14X
Buser/Modified PH2MBH Large	PPBUH2MX
Buser/Modified PH2MBH Small	PPBUH2M2X
Swiss Perio Modified P24GSP Periosteal	SPP24GSPX
Swiss Perio Modified P24GSIC Periosteal	SPP24GSICX
Swiss Perio Modified P24GSP Periosteal 2	SPP24GSP2X
Swiss Perio Modified Schlee Periosteal	SPPSCHLEEX
Swiss Perio Modified PH2 Periosteal	SPPH2X
Urban Mini Me Periosteal	PFIWDS1MKX

LUXATING HYBRIDS / SYNDESMOTOMES

Luxating Hybrids, Straight Spade	EPTSSPX
Luxating Hybrids, Small Straight	EPTSMSX
Luxating Hybrids, Small Curved	EPTSMCX
Luxating Hybrids, Small Distal	EPTSMDX
Luxating Hybrids, Small Mesial	EPTSMMX
Luxating Hybrids, Large Distal	EPTLGDX
Luxating Hybrids, Large Mesial	EPTLGMX
Luxating Hybrid Kit, Black Line	EPTKITX
6 Right Bernard Syndesmotomes	BER6X
7 Left Bernard Syndesmotomes	BER7X
11 Straight Bernard Syndesmotomes	BER11X

TUNNELING INSTRUMENTS

Tunneling Knife 1	TKN1X
Tunneling Knife 1, Short Handle	TKN1SX
Tunneling Knife 2	TKN2X
Tunneling Knife 2, Short Handle	TKN2SX
Giles de Quincey Tunneling Instrument 1	FQUINCEY1X
Giles de Quincey Tunneling Instrument 2	FQUINCEY2X

PERIODONTAL KNIVES

1/2 Allen Orban Knife	KO12KPO3AX
1/2 Allen Orban Knife, Round	KO12KPO3RX
Allen Arrowhead Knife	KO12KP3R4X
Allen End-Cutting Intralucular Knife	KPAX
15/16 Kirkland Knife	KK15/16X
15/16 Kirkland Knife, Short Handle	KK15/16SX
1/2 Orban Knife	KO1/2X
11 Goldman-Fox Knife	KGF11X

PERIODONTAL BONE CHISELS	
1/2 Buser Modified Periodontal Chisel, 3mm/4mm	PPBUSER12X
2/3 Buser Modified Periodontal Chisel, 4mm/5mm	PPBUSER23X
1/3 Buser Modified Periodontal Chisel, 3mm/5mm	PPBUSER13X
13K/TG Periodontal Chisel	S13K/TGX
13K/13KL Periodontal Chisel	S13K/13KLX
Palacci Chisel	SPAL13KX

SINUS LIFT INSTRUMENTS	
Kramer-Nevins Sinus Lift, Obtuse	IMP6578SX
Kramer-Nevins Sinus Lift, Acute	IMP6577SCX
Sinus Curette 1	SINC1X
Sinus Curette 1, Short Handle	SINC1SX
Sinus Curette 2	SINC2X
Sinus Curette 2, Short Handle	SINC2SX
300 Sinus Lift Instrument	IMP300X

SURGICAL CURETTES	
2/4 Molt Curette, Double-Ended	CM2/4X
84 Lucas Surgical Curette, Double-Ended	CL84X
85 Lucas Surgical Curette, Double-Ended	CL85X
86 Lucas Surgical Curette, Double-Ended	CL86X
87 Lucas Surgical Curette, Double-Ended	CL87X
88 Lucas Surgical Curette, Double Ended	CL88X

IMPLANTOLOGY	
Membrane Placement Instrument	PMPIX
Grisdale Bone Condenser	PLGGR1X
Labanca Plugger, 3 x 1.5mm	PLGLABANC3X
Jovanovic Plugger	PLGJO/4X
Rompen Plugger, Buccal/Lingual	PLGROMP1X
Rompen Plugger, Mesial/Distal	PLGROMP2X
Palti Bone Graft Scoop	IMP6522X

PERIOTOMES	
Periotome, Posterior, Double-Ended	PT1X
Periotome, Anterior, Double-Ended	PT2X
Periotome, Right, Single-Ended	PT3X
Periotome, Left, Single-Ended	PT4X
Periotome, Anterior Angled, Single-Ended	PT5X
Periotome, Anterior Straight, Single-Ended	PT6X
Wise Periotome	PEREURX

LUXATING ELEVATORS	
Luxating Elevator, 2mm, Straight	EL2SX
Luxating Elevator, 2mm, Curved	EL2CX
Luxating Elevator, 3mm, Straight	EL3SX
Luxating Elevator, 3mm, Curved	EL3CX
Luxating Elevator, 4mm, Straight	EL4SX
Luxating Elevator, 4mm, Curved	EL4CX
Luxating Elevator, 5mm, Straight	EL5SX
Luxating Elevator, 5mm, Curved	EL5CX

RETRACTOR	
Minnesota Retractor	CRMX

ROOT TIP PICK	
13/14 Heidbrink Root Tip Pick	EHB13/14X

**The
Best In
Practice**

Learn at **Hu-Friedy.com/Black-Line**

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.

©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-330/0421