

CASSETTE CAPACITY CHART

CLEANING EQUIPMENT



IMS™ CASSETTES				RESIN CASSETTES													
MANUFACTURER PHONE	MODEL	TANK DIMENSIONS L x W x D (INCHES)	CAPACITY IN GALLONS	EXAM & HANDPIECE 3 x 8 x 1.25	SMALL 10 INSTRUMENT 5.5 x 8 x 1.25	TWO-TIER 8-INSTRUMENT 3 x 8 x 1.5	SPACE SAVER SEVEN 4.5 x 8 x 1	TWO-TIER 14-INSTRUMENT 4.5 x 8 x 1.5	SIGNA STAT 12 INSTRUMENT 6.5 x 10.5 x 1.25	TWO-TIER 20-INSTRUMENT 6 x 9.7 x 1.5	LARGE 8 x 11 x 1.25	ORAL SURGERY 8 x 14.5 x 1.25	SMALL THIN ORTHO 7 x 8 x 0.75	LARGE THIN ORTHO 8 x 11 x 0.75	TWO-TIER ORTHO 20-INSTRUMENT 6 x 8 x 1.5	SMALL & TWO-TIER 4.5 x 8 x 1.5	LARGE 8 x 11 x 1.5
BrandMax® Tri-Clean™ § 1-844-874-2532	U-10LH / U-10LHREC	12 x 9.5 x 6	2.64	8	3	9	8	6	N/C	3	N/C	N/C	6	N/C	3	3	N/C
	U-13LH / U-13LHREC✓	13 x 11.8 x 6	3.43	9	6	9	8	6	3	3	3	N/C	6	6	3	6	3
	U-19LH/U-19LHREC	13 x 11.8 x 8	4.94	12	10	12	12	8	4	4	5	N/C	8	8	4	8	4
	U-20LH / U-20LHREC✓	20 x 12 x 6	5.06	18	12	15	16	12	6	6	6	3	12	6	9	12	6
Biotrol™ 1-800-822-8550	Purit™ 300*✓	12.25 x 11.25 x 5.5	3	12	6	12	9	7	3	3	3	N/C	5	5	3	7	3
	Purit™ 150*	11 x 5.25 x 5.625	1.5	4	1	4	4	3	N/C	N/C	N/C	N/C	N/C	N/C	N/C	3	N/C
Coltene Whaledent 1-800-221-3046	BioSonic® UC300 (cassette rack accessory P/N UC320 or UC325) ✓	14 x 9 x 6	3	9	6*	9	6	6	3	3	3	N/C	3	3	3	6	3
	BioSonic® UC300R (cassette rack accessory P/N UC320 or UC325) ✓	14 x 9 x 6	3	9	6*	9	6	6	3	3	3	N/C	3	3	3	6	3
	BioSonic® UC150 (Accessory) with basket only	10.4 x 6.2 x 4.4	1.5	4	3	4	4	2	N/C	N/C	N/C	N/C	N/C	N/C	N/C	2	N/C
	BioSonic® UC125/UC100XD (UC164XD)	11.25 x 5.62 x 5	1.25	1	N/C	1	3	1	N/C	N/C	N/C	N/C	N/C	N/C	N/C	1	N/C
	BioSonic® UC125/UC100XD (UC151XD - with original basket) *	11.25 x 5.62 x 5	1.25	4	N/C	3	3	2	N/C	N/C	N/C	N/C	N/C	N/C	1	2	N/C
	BioSonic® UC125/UC100XD (UC151XD - with updated basket) †	11.25 x 5.62 x 5	1.25	4	N/C	3	3	2	N/C	N/C	N/C	N/C	N/C	N/C	1	2	N/C
	BioSonic® UC100 *	11.5 x 5.25 x 4	1	1	N/C	1	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C
Esma, Inc. 1-800-276-2466	E789✓	13 x 11 x 11	6.5	20	12	12	12	9	8	8	6	N/C	10	10	8	9	6
	E789-oms	15 x 10 x 10	6.5							8	6	6	10	10	10	9	6
	E397	12 x 10 x 8	3.5	12	7	9	8	6	5	4	3	N/C	7	7	5	6	3
	E289	12 x 10 x 8	3.5	12	7	9	8	6	5	4	3	N/C	7	7	4	6	3
	E888R	12 x 10 x 8	3.5	12	7	9	8	6	5	4	3	N/C	7	7	4	6	3
	E888	12 x 10 x 8	3.5	12	7	9	8	6	5	4	3	N/C	7	7	4	6	3
Health Sonics 1-800-342-3096	T19.9*✓	20.8 x 13.5 x 12.2	5	17	9	17	16	8	3	2	3	3	10	5	4	8	3
	T13.7*✓	14 x 14 x 10.7	3.5	12	6	12	9	8	3	2	3	N/C	5	5	2	8	3
	T10.4*	13 x 11.5 x 10.3	2.7	10	3	10	8	6	N/C	3	N/C	N/C	5	N/C	3	6	N/C
	T10.4 C/S*	13 x 11.5 x 12.3	3.2	12	7	12	12	8	N/C	4	N/C	N/C	8	N/C	4	8	N/C
	T9*	20.8 x 7.5 x 11.5	2.3	6	N/C	6	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C
	T5.7*	12.8 x 8 x 10.5	1.5	4	N/C	4	4	3	N/C	N/C	N/C	N/C	N/C	N/C	N/C	3	N/C
	T4.4*	10.5 x 7.5 x 10	1.1	1	N/C	1	1	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C
	T3.3*	10.5 x 7.5 x 8.5	0.8	2	N/C	2	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C

All capacity quantities and dimensions are provided by listed equipment manufacturers. They are measured utilizing wrapped or unwrapped cassettes, if your office uses sterilization pouches the capacity may be lower than listed. This data does not take into account weight limits, for which the manufacturer should be consulted. Hu-Friedy provides this for your convenience, but we encourage you to verify with equipment manufacturer before purchase. We are happy to provide sample cassettes to confirm fit.

N/C Not Compatible
■ Must be Hu-Friedy part code
Infinity Series — IMN416 (color) Q - 8 x 10.5
Signature Series — IM416 (color) Q - 8 x 10.5
to accommodate rounded corner of tank

~ According to equipment Manufacturer, when the UC300/UC300R is used with Hu-Friedy 10-Instrument Cassettes, in order to obtain full cassette coverage, it may be necessary to fill the tank beyond the fill line. Ensure the water level does not exceed the lip of the plastic bezel inside the tank. This action will not void the standard warranty.

+ Recessed model
* Discontinued model
✓ Model most compatible with Hu-Friedy Cassettes
• Original style basket with integrated instrument holders on bottom
‡ New style basket with flat bottom

† Must order 2 x 4XL Cassette Racks (Full Size) — part #01-109965S
§ In order to obtain full cassette coverage, it may be necessary to fill tank beyond the fill line. Please ensure water does not extend above top of tank. Capacity show with optimal instrument cassette basket.

CASSETTE CAPACITY CHART

CLEANING EQUIPMENT (CONTINUED)



IMS* CASSETTES																	RESIN CASSETTES	
MANUFACTURER PHONE	MODEL	TANK DIMENSIONS L x W x D (INCHES)	CAPACITY IN GALLONS	EXAM & HANDPIECE 3x8x1.25	SMALL 10 INSTRUMENT 5.5x8x1.25	TWO-TIER 8-INSTRUMENT 3x8x1.5	SPACE SAVER SEVEN 4.5x8x1	TWO-TIER 14-INSTRUMENT 4.5x8x1.5	SIGNA STAT 12 INSTRUMENT 6.5x10.5x1.25	TWO-TIER 20-INSTRUMENT 6x9.7x1.5	LARGE 8 x 11x 1.25	ORAL SURGERY 8x14.5x1.25	SMALL THIN ORTHO 7x8x0.75	LARGE THIN ORTHO 8x11x0.75	TWO-TIER ORTHO 20-INSTRUMENT 6x8x1.5	SMALL & TWO-TIER 4.5x8x1.5	LARGE 8x11x1.5	
L&R 201-991-5330	SweepZone® 200 (No Rack – Use Basket)	11.75 x 6 x 4	1.06	2	N/C	1	1	1	N/C	N/C	N/C	N/C	N/C	N/C	N/C	2	N/C	
	SweepZone® 310R / SweepZone 310 / Quantrex® 310 (Rack #26665)+✓	15.5 x 9 x 6	3.25	12	6	8	6	4	3	2	3	N/C	8	4	4	4	2	
	SweepZone® 360R / Quantrex® 360 Rack #19368+	11.5 x 9.25 x 8	3.6	12	4	12	10	8	4	4	N/C	N/C	7	N/C	3	6	N/C	
	SweepZone® 650R / Quantrex® 650 Rack #26847 (Holds 2 Racks)+✓	19.75 x 11.75 x 6.5	6.5	18	6	18	16	12	6	6	6	N/C	10	10	6	12	6	
	SweepZone® 650R / Quantrex® 650 Rack #19369 (Oral Surgery)	19.75 x 11.75 x 6.5	6.5	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	3	N/C	N/C	N/C	N/C	N/C	
Midmark® QuickClean® 1-800-MIDMARK	QuickClean® QC1	11.8 x 5.9 x 5.9	1.2	3	0	3	1	1	0	0	0	0	0	0	0	1	0	
	QuickClean® QC3/QC3R✓	13 x 11.8 x 5.9	3.3	8	4	6	6	4	3	2	2	0	4	4	2	4	2	
	QuickClean® QC6/QC6R✓	19.9 x 11.8 x 7.8	6.6	20	12	15	15	12	4	8	4	4	12	6	6	12	4	
Midmark® Soniclean® 1-800-643-6275	Soniclean® M550 (Rack #9A407001)*✓	19.5 x 11.5 x 6	5.5	15	9	15	16	12	6	6	6	3	10	10	9	12	6	
	Soniclean® M550 with basket *			15	9	15	16	12	6	6	3	3	10	5	9	12	3	
	Soniclean® M250 (Rack #9A295001)*	11.5 x 9.5 x 6	2.5	9	3	9	8	4	3	3	3■	N/C	6	N/C	3	6	3■	
	Soniclean® M250 with basket *			9	3	9	8	4	3	3	0	0	6	5	3	6	1	
	Soniclean® M150*	11.5 x 6 x 6	1.5	3	0	3	4	3	N/C	N/C	N/C	N/C	N/C	N/C	N/C	3	N/C	
Miele Canada 1-888-325-3957 USA 1-80-991-9380	G7881/G7891 (Rack #4xE523)✓	19.6 x 21 x 19.6	N/A	55	35	55	35	20	11	55	11	11	24	11	24	35	11	
Quala® 1-800-929-4ADC	Sweep 5300✓ (Rack #26665 for large cassettes)	15.5 x 9 x 6	3.25	12	6	8	8	6	3	2	3	N/C	5	5	6	6	3	
	Sweep 5200	11.75 x 6 x 4	1.06	2	N/C	2	2	0	N/C	N/C	N/C	N/C	N/C	N/C	N/C	2	N/C	
	5002	9.4 x 5.4 x 4	0.84	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	
SciCan 1-800-667-7733	Hydrim® L110W G4✓ (*w/ modified racks)	23 x 15.5 x 18	N/A	30	20	30	15	15	5/**10	5/**10	**10	**8	10	**10	10	15	**10	
	New Hydrim® C61W G4✓ (Std rack / 6 rack)	12 x 19.5 x 16	N/A	12/24	6/12	12/24	9/18	9/18	3/6	3/6	3/6	3/6	6/12	3/6	6/12	9/18	3/6	
Sonix IV 1-843-554-0239 1-800-878-2769 Ortho Rack# (SS6A-TUCH128) Standard Rack# (SS6A-TUCH266)	SR118/ST118✓	19.5 x 11.5 x 8	7.8	24	12	20	20	16	4	10	4	4	12	12	10	16	4	
	SR116/ST116	19.5 x 11.5 x 6	5.5	18	9	15	14	12	3	6	3	3	8	8	6	12	3	
	SR236/ST236	11.5 x 12.8 x 6	3.5	9	6	9	10	7	3	3	3	N/C	4	4	3	7	3	
	SR128/ST128✓	11.5 x 9.5 x 8	3.5	12	8	12	10	8	4	4	4	N/C	6	N/C	4	8	N/C	
	SR126/ST126	11.5 x 9.5 x 6	2.7	6	3	9	6	6	3	3	N/C	N/C	4	N/C	3	6	N/C	
	SR238/ST238	11.5 x 12.8 x 8	4.5	9	6	9	10	8	4	4	4	N/C	6	6	4	8	4	
Sultan 1-800-637-8582	ProSonic 2000 (2 Racks #60230)*✓	19.5 x 11.5 x 6	5.5	18	12	18	16	12	6	6	6	4	10	10	6	12	6	
	ProSonic 1000 (1 Rack #60230)*✓	11.5 x 9.5 x 6	2.5	12	4	9	8	6	3	3	3■	N/C	5	N/C	3	6	3	
	ProSonic 600*	11.5 x 6 x 6	1.5	4	4	3	4	3	N/C	3	N/C	N/C	N/C	N/C	N/C	3	N/C	
	ProSonic 300*	9.5 x 5.5 x 4	0.75	3	N/C	2	N/C	2	N/C	N/C	N/C	N/C	N/C	N/C	N/C	2	N/C	
	Resurge™ Ultrasonic Cleaner	9.8 x 6.9 x 3	0.8	2	10	1	2	1	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	N/C	
	Resurge™ Ultrasonic Cleaner	13 x 9.4 x 4.1	2	6	4	6	6	4	2	2	2	N/C	4	4	2	4	2	
Tuttnauer 1-800-624-5836	CSU3✓	15.5 x 9 x 6	3	9	6	6	6	6	4	3	3	N/C	5	5	6	6	3	
	CSU1	9.38 x 5.38 x 4	0.8	2	N/C	2	2	1	N/C	N/C	N/C	N/C	N/C	N/C	N/C	1	N/C	

All capacity quantities and dimensions are provided by listed equipment manufacturers. They are measured utilizing wrapped or unwrapped cassettes, if your office uses sterilization pouches the capacity may be lower than listed. This data does not take into account weight limits, for which the manufacturer should be consulted. Hu-Friedy provides this for your convenience, but we encourage you to verify with equipment manufacturer before purchase. We are happy to provide sample cassettes to confirm fit.

* Discontinued models
N/C Not Compatible
■ Must be Hu-Friedy part code
Infinity Series — IMN416 (color) Q - 8 x 10.5
Signature Series — IM416 (color) Q - 8 x 10.5
to accommodate rounded corner of tank

^ According to equipment Manufacturer, when the UC300/UC300R is used with Hu-Friedy 10-Instrument Cassettes, in order to obtain full cassette coverage, it may be necessary to fill the tank beyond the fill line. Ensure the water level does not exceed the lip of the plastic bezel inside the tank. This action will not void the standard warranty.

+ Recessed model
* Discontinued model
✓ Model most compatible with Hu-Friedy Cassettes
• Original style basket with integrated instrument holders on bottom
▸ New style basket with flat bottom

† Must order 2 x 4XL Cassette Racks (Full Size) — part #01-109965S
§ In order to obtain full cassette coverage, it may be necessary to fill tank beyond the fill line. Please ensure water does not extend above top of tank. Capacity show with optimal instrument cassette basket.

CASSETTE CAPACITY CHART

STERILIZER EQUIPMENT



IMS™ CASSETTES															RESIN CASSETTES	
MANUFACTURER PHONE	MODEL	CHAMBER DIMENSIONS L x W x D (INCHES)	EXAM & HANDPIECE 3 x 8 x 1.25	SMALL 10 INSTRUMENT 5.5 x 8 x 1.25	TWO-TIER 8-INSTRUMENT 3 x 8 x 1.5	SPACE SAVER SEVEN 4.5 x 8 x 1	TWO-TIER 14-INSTRUMENT 4.5 x 8 x 1.5	SIGNA STAT 12 INSTRUMENT 6.5 x 10.5 x 1.25	TWO-TIER 20-INSTRUMENT 6 x 9.7 x 1.5	LARGE 8 x 11 x 1.25	ORAL SURGERY 8 x 14.5 x 1.25	SMALL THIN ORTHO 7 x 8 x 0.75	LARGE THIN ORTHO 8 x 11 x 0.75	TWO-TIER ORTHO 20-INSTRUMENT 6 x 8 x 1.5	SMALL & TWO-TIER 4.5 x 8 x 1.5	LARGE 8 x 11 x 1.5
A-dec® / W&H 1-800-547-1883	LEXA MN-111	11x18	16	8	8	12	6^	4	3	4	4	8	4	4	6^	2^
	LISA MB 17✓	9.8 x 13.8	9	6	N/C	6	N/C	N/C	N/C	3	N/C	3	3	N/C	N/C	N/C
	LISA MB 22 (Available in Canada)✓	9.8 x 17.7	15	6	N/C	9	N/C	3	N/C	3	3	6	3	N/C	N/C	N/C
Midmark® 1-800-MIDMARK 1-800-643-6275	M11w/ Horizontal Rack #9A215001✓	11x18	16	8	16	12	12	N/C	N/C	4	4	8	4	8	12	4
	M11w/ Vertical Rack #9A215002✓	11x18	16	8	16	12	12	4	4	4	4	8	4	4	12	4
	M11 with trays	11x18	12	6	12	8	8	3	4	2	2	6	3	6	8	2
	M9 standard rack with trays	9 x 15	6	3	6	4	3	3	3	N/C	N/C	4	N/C	3	3	N/C
	M3 with tray	7.6 x 1.6 x 12.1	2	1	N/C	1	N/C	1	N/C	N/C	N/C	1	N/C	N/C	N/C	N/C
Miele Canada 1-888-325-3957 USA 1-80-991-9380	PS1201B, Class B✓	10x16	9	6	6	6	3	3	3	3	2	6	3	3	3	N/C
Mocom +39 051703168	B Futura/Classic 28✓	11 x 17.5	12	8	12	12	8	4	N/C	4	4	8	4	8	12	4
Pelton & Crane 1-800-659-6560	Magnaclave✓	15.25x15x26	28	28	28	28	28	14	14	14	3	14	14	14	28	14
	Delta XL10 w/ Rack #024338✓	9.8 x 17.5	15	9	12	9	9	3	3	3	3	4	3	9	9	3
	Delta XL8	8.25 x 14	8	1	8	1	N/A	0	0	N/C	N/C	2	N/C	1	4	N/C
	Delta Q 10 Rack #024338*✓	9.8 x 17.5	15	9	12	9	9	3	3	3	3	4	3	9	9	3
	Delta Q 8*	8.25 x 14	8	1	8	1	N/A	0	0	N/C	N/C	2	N/C	1	4	N/C
SciCan Canada 1-800-667-7733 USA 1-800-572-1211	Bravo™17✓	10 x 13.5	9	6	9	6	6	3▣	3	3	N/A	5	3	3	6	3
	Bravo™17V✓	10 x 13.5	9	6	9	6	6	3▣	3	3	N/A	5	3	3	6	3
	Bravo™21V✓	10 x 17.5	12	6	12	9	9	3▣	3	3	3	6	3	6	9	3
	Quantim 16✓	9.8 x 13	9	6	9	6	6	3	3	3	N/C	3	3	3	6	3
	Statim® 7000	8.6x13.5x2.5	4	4	4	4	2	2	1	2	N/C	3	3	3	2	2
	Statim® 5000 G4	7.1x15 x 3.2	2	2	2	2	1	2	1	N/C	N/C	3	N/C	N/A	1	N/C
	Statim® 2000 G4	7.1x11 x 1.6	2	1	1	1	1	1	N/C	N/C	N/C	1	N/C	N/A	N/C	N/C
	StatClave® G4	11x15	16	8	16	12	12	4	4	4	4	8	4	8	12	4
Tuttnauer 1-800-624-5836	3870*✓	15x30	31	20	23	20	20	11	11	10	7	17	14	14	35	10
	EZ™10✓	10 x 19	15	10	11	10	10	5	5	5	3	8	7	7	17	5
	2540M✓	10 x 19	15	10	11	10	10	5	5	5	3	8	7	7	17	5
	EZ™9	9 x 18	12	8	9	8	8	4	4	4	3	7	5	5	14	4
	2340M	9 x 18	12	8	9	8	8	4	4	4	3	7	5	5	14	4
	EZ™9Plus	9x19.8	16	10	12	10	20	6	4	5	3	9	7	7	17	5
	Valueklave™ 1730	7 x 12	5	3	4	3	3	0	0	0	0	0	0	0	6	0
	EZ™11Plus✓	11x19.8	26	17	20	17	17	10	9	8	6	15	10	12	30	8
	Elara 11 Class B✓	11x19.8	31	20	23	20	20	11	11	10	7	17	14	14	35	10

All capacity quantities and dimensions are provided by listed equipment manufacturers. They are measured utilizing wrapped or unwrapped cassettes, if your office uses sterilization pouches the capacity may be lower than listed. This data does not take into account weight limits, for which the manufacturer should be consulted. Hu-Friedy provides this for your convenience, but we encourage you to verify with equipment manufacturer before purchase. We are happy to provide sample cassettes to confirm fit.

N/A Not Available
N/C Not Compatible
■ Must be Hu-Friedy part code
Infinity Series — IMN416 (color) Q - 8 x 10.5
Signature Series — IM416 (color) Q - 8 x 10.5
to accommodate rounded corner of tank

* Discontinued model
✓ Model most compatible with Hu-Friedy Cassettes
** Optional Rack required for these cassettes
▣ Loose on the rails

♦ Does not fit within rail
^ Insert tray on bottom shelf and then place ortho boxes
° Data not available from equipment manufacturer
Please work with the manufacturer for capacity information.