



Introducing Hu-Friedy's New **ATRAUMAIR FORCEPS**

Hu-Friedy's new line of Atraumair Forceps are designed for optimal access and enhanced crown and root engagement.

Each forceps is meticulously crafted to the highest quality standard, ensuring superior performance and longevity. Atraumair Forceps are available in 13 unique patterns that provide a complete line for a variety of extraction cases.



ATRAUMAIR POINTS OF PERFORMANCE:

- Enhanced beak geometries and longitudinal beak serrations provide a **firm, secure grasp of the crown and root reducing the risk of crushing or fracturing the tooth**
- Laser marked hinge for **greater visibility when bagged for sterilization**
- Matte finish for **enhanced contrast and reduced light reflection**
- Unique hole design reduces weight of the forceps, providing **increased comfort and control**
- Meticulously crafted from the finest surgical grade stainless steel

How the best perform

Hu-Friedy

ATRAUMAIR FORCEPS

67 Forceps
Upper Molars
| **F67XS**



1 Standard Forceps
Upper Incisors
| **F1XS**



1 Apical Forceps
Upper Incisors
| **FAF1XS**



35 Forceps, Upper
Canines & Premolars
| **F35XS**



35 Apical Forceps, Upper
Canines & Premolars
| **FAF35XS**



51 Root Forceps
Upper Roots, Serrated
| **FX51XS**



36 Forceps
Lower Incisors
& Canines
| **F36XS**



36 Apical Forceps
Lower Premolars,
Canines & Incisors
| **FAF36XS**



13 Forceps
Lower Premolars
| **F13XS**



79 Forceps
Lower Molars
| **FX79XS**



73 Forceps
Lower Molars
| **FX73XS**



23 Cowhorn Forceps
Lower Molars
| **F23XS**



45 Root Forceps
Lower Roots, Serrated
| **FX45XS**