

POWER PLUS INSERTS


REACH FURTHER FOR POCKET DEBRIDEMENT

New **POWER PLUS** inserts from Hu-Friedy feature **POWER** for effective removal of tenacious deposits and heavy extrinsic stain **PLUS** extended length to reach deep into periodontal pockets. With a 25% longer tip design*, a variety of tip sizes, and a focus on clinician comfort, Power PLUS inserts offer an advanced solution for debridement of tenacious deposits.



FEATURES:

- **2+ mm longer tip*** for extended reach into periodontal pockets
- **Tip options** to accommodate diverse pocket anatomy
 - Thinner tip for tenacious deposits in narrow pockets
 - Standard size tip to reach heavy debris in deeper pockets
- **Tip shape similar to periodontal probe** for reach, familiarity and ease of use
- **Special blend of stainless steel** for durability and tip longevity
- **Improved clinician comfort**
 - Popular Swivel™ handle for excellent control and maneuverability
 - Ergonomic wide diameter handle for reduced finger pinching
 - Textured silicone grip for secure grasp
 - Uniform diameter handle for ease of movement when shifting grip location



Power PLUS Standard
| UI30PPSC

Power PLUS Thin
| UI30PPTC

	Tip Description	Tip Cross-Section	Recommendation	Part Code
Power PLUS Standard	Standard diameter shank extending to a tapered tip	Round	Moderate to heavy deposits Tenacious deposits Orthodontic cement Wide, deep distended pockets	UI30PPSC
Power PLUS Thin	Thinner diameter shank extending to a tapered tip	Round	Light to moderate deposits Narrow, deep pockets Crowded supragingival areas	UI30PPTC