



EMPOWERMENT KIT

The perfect set of scalers to enhance the effectiveness of your Ultrasonic and AIRFLOW treatment protocol

COLOR CODED

ACCESSORY AREA


ON AVERAGE
10% THINNER
BLADES

How the best perform


Hu-Friedy

Easy to Identify

The Empowerment Kit **instruments are color coded** to easily identify scalers and curettes that correspond to the anatomical areas of the mouth and tooth surfaces. This will help maximize efficiency by reducing the time spent identifying different working ends, allowing the clinician to continue moving at the same fast pace they are accustomed to when using power therapies.



UNIVERSAL SCALERS



Also included in the Empowerment Kit are Universal Scalers. The versatility of the Universal Scaler enables it to be used on any area of the mouth.

Access tight, hard to reach places with an average of 10% thinner blades.

The most common place where calculus is left behind is the cervical area of the tooth^[1]. The 10% thinner blades^[2] of the Empowerment Kit are specifically designed for easier access to tight interproximal spaces and allow for easier insertion subgingivally. Combined with the clinically proven sharp cutting edge of EverEdge® 2.0 scalers^[2], clinicians can effectively remove the calculus subgingivally, reducing the risk of leaving behind calculus.




[1] Cobb, C. M. (1996). Non-surgical pocket therapy: mechanical. Annals of Periodontology, 1, 443-490.

[2] Based on average difference when compared to the standard Hu-Friedy instrument. Data on file, available upon request.

Customize Your Kit


The Empowerment Kit includes seven scalers. We also offer four alternative scalers to accomodate all preferences.

The Nevi 1 Empowerment Scaler has an elongated disk end for anterior lingual stain and calculus removal. This is paired with a super thin, anterior curved sickle scaler.




Nevi 1
Empowerment Scaler
| SCNEVI1LPC9E2

Terminal Shank: standard
Blade Length: standard
Blade Width: -11%




#13/14
Empowerment Gracey
| SG13/14PC9E2

Terminal Shank: +20%
Blade Length: -12%
Blade Width: -18%




ALTERNATIVE
AVAILABLE



#11/12
Empowerment Gracey
| SG11/12PC9E2

Terminal Shank: +1%
Blade Length: -24%
Blade Width: -18%




ALTERNATIVE
AVAILABLE




#17/18 Mini 5
Empowerment Gracey
| SAS17/18PC9E2

Terminal Shank: +2%
Blade Length: -39%
Blade Width: -9%



#15/16 After 5
Empowerment Gracey
| SRPG15/16PC9E2


Terminal Shank: +32%
Blade Length: standard
Blade Width: -9%






**#7/8
Empowerment Gracey**
| SG7/8PC9E2

Terminal Shank: +11%
Blade Length: -7%
Blade Width: -18%




**Nevi 4
Empowerment Scaler**
| SCNEVI4LPC9E2

Terminal Shank: standard
Blade Length: standard
Blade Width: -9%



**#1/2
Empowerment Gracey**
| SG1/2PC9E2

Terminal Shank: +56%
Blade Length: -24%
Blade Width: -18%




ALTERNATIVE
AVAILABLE



**#7/8 Empowerment
Younger Good Scaler**
| SYG7/8PC9E2

Terminal Shank: standard
Blade Length: standard
Blade Width: standard




ALTERNATIVE
AVAILABLE




**#1/2 After 5
Empowerment Gracey**
| SRPG1/2PC9E2

Terminal Shank: +67%
Blade Length: standard
Blade Width: -9%



**#3/4
Empowerment Langer**
| SL3/4PC9E2

Terminal Shank: standard
Blade Length: standard
Blade Width: standard



EMPOWERMENT KIT**| POWCOMPKIT**

Infinity 8 Instrument Cassette, Purple

| IMN5084

#11/12 Old Dominion University Explorer

| EXD11/126

■ Nevi 4 Empowerment Scaler

| SCNEVI4LPC9E2

■ #7/8 Empowerment Gracey

| SG7/8PC9E2

■ #13/14 Empowerment Gracey

| SG13/14PC9E2

■ #11/12 Empowerment Gracey

| SG11/12PC9E2

■ #1/2 Empowerment Gracey

| SG1/2PC9E2

■ #7/8 Empowerment Younger Good Scaler

| SYG7/8PC9E2

■ Nevi 1 Empowerment Scaler

| SCNEVI1LPC9E2**Additional Scaler Options**

■ #17/18 Mini 5 Empowerment Gracey

| SAS17/18PC9E2

■ #15/16 After 5 Empowerment Gracey

| SRPG15/16PC9E2


■ #1/2 After 5 Empowerment Gracey

| SRPG1/2PC9E2

■ #3/4 Empowerment Langer

| SL3/4PC9E2**Accessory Area**

The Empowerment Kit is the perfect complement to any power therapy. The Infinity Cassette comes with an accessory area to conveniently hold your Ultrasonic or AIRFLOW® handpiece.




PERIODONTAL



RESTORATIVE




ORTHODONTIC

INSTRUMENT
MANAGEMENT

ENDODONTIC



SURGICAL

POWER
SCALING

For more, visit

HU-FRIEDY.COM/EmpowermentKit

Call 1-800-Hu-Friedy or contact your authorized Hu-Friedy representative for more information.

©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-302/0120

