



HU-FRIEDY SURGICAL KIT

LUXATING HYBRID KIT, BLACK LINE



PART CODE	DESCRIPTION
IMN5081	8 Instrument Infinity Series™ Cassette, Gray
EPTSSPX	Luxating Hybrid Elevator, Straight Spade, Black Line
EPTSMSX	Luxating Hybrid Elevator, Small Straight, Black Line
EPTSMCX	Luxating Hybrid - Small Curved, Black Line
EPTSMDX	Luxating Hybrid - Small Distal, Black Line
EPTSMMX	Luxating Hybrid - Small Mesial, Black Line
EPTLGDX	Luxating Hybrid - Large Distal, Black Line
EPTLGMX	Luxating Hybrid - Large Mesial, Black Line
EPTKITX	Complete Luxating Hybrid Kit, Black Line

LUXATING HYBRID KIT, BLACK LINE



LUXATING HYBRID KIT, BLACK LINE



- A IMN5081**
8 Instrument Infinity Series™ Cassette, Gray
IMS® cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



- B EPTSSPX**
Luxating Hybrid Elevator, Straight Spade, Black Line
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- C EPTSMSX**
Luxating Hybrid Elevator, Small Straight, Black Line
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- D EPTSMCX**
Luxating Hybrid - Small Curved, Black Line
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- E EPTSMDX**
Luxating Hybrid - Small Distal, Black Line
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- F EPTSMMX**
Luxating Hybrid - Small Mesial, Black Line
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- G EPTLGDX**
Luxating Hybrid - Large Distal, Black Line
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- H EPTLGMX**
Luxating Hybrid - Large Mesial, Black Line
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.