



HU-FRIEDY SURGICAL KIT

BONE GRAFTING KIT



PART CODE	DESCRIPTION
IM4128	12 Instrument Signature Series® Cassette, Blue
10-130-05	5 Scalpel Handle, Round Straight
CSEUWES6	Euwe Bone Scoop/Plugger, #6 Handle
SP20	20 Corn Suture Pliers
PPBUSER6	Buser Periosteal, #6 Handle
P96	9 Molt Periosteal Elevator
PPBUSE1/36	PJO/1 & PJO/3 Combined 13KL, #6 Handle
PPBUSE2/36	PJO/2 & PJO/3 Combined 13KL, #6 Handle
CL866	86 Lucas DE Surgical Curette, #6 Handle
SPR1/26	PR1/2 Prichard Curette, #6 Handle
MCUPE	Measuring cup 25ml with Marking
S16	16 Goldman-Fox Scissors
NHHB	Hegar-Baumgartner Needle Holder
PMPI	Membrane Placement Instrument
IM1002	6-Pack Signature Series® Instrument Clip
IMBGKIT	Complete Bone Grafting Kit

BONE GRAFTING KIT



BONE GRAFTING KIT



A IM4128
12 Instrument Signature Series® Cassette, Blue

IMS® cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



B 10-130-05
5 Scalpel Handle, Round Straight

Used for cutting gingival tissue and making surgical incisions. Provides greater rotational movement with fingertip control and enhances dexterity. Handle measures 6" long. Compatible with standard scalpel blade sizes.



C CSEUWES6
Euwe Bone Scoop/Plugger, #6 Handle

Designed for easy transport and non-stick compaction of bone graft material. 3.0mm/7.0mm.



D SP20
20 Corn Suture Pliers

Ideal for use with guided tissue regeneration procedures. 15cm



E PPBUSER
Buser Periosteal, #6 Handle

Designed for reflecting and retracting the muco-periosteum after incisions of the gingival tissue. 3.8mm/4.3mm



F P96
9 Molt Periosteal Elevator

Designed for reflecting and retracting the muco-periosteum after incisions of the gingival tissue. 4.3mm/8.4mm



G PPBUSE1/36
PJO/1 & PJO/3 Combined 13KL, #6 Handle

Used to remove autologous bone for transplantation material with 4mm and 6mm wide working ends.



H PPBUSE2/36
PJO/2 & PJO/3 Combined 13KL, #6 Handle

Used to to remove autologous bone for transplantation material with 5mm and 6mm wide working ends.



I CL866
86 Lucas DE Surgical Curette, #6 Handle

Designed for curettage, cyst removal and tooth socket debridement. 2.8mm



J SPR1/26
PR1/2 Prichard Curette, #6 Handle

Larger and heavier curette for the removal of granulation tissue and tenacious subgingival deposits measuring 7.5" long with 1.5mm wide working ends.



K MCUPE
Measuring Cup 25ml with Marking



L S16
16 Goldman-Fox Scissors

Used for trimming tissue or cutting sutures. 12.5cm with 1 serrated blade.




M NHHB
Hegar-Baumgartner Needle Holder

Used to grasp and guide the needle when suturing. 14 cm. Recommended suture sizes: 1-0, 2-0, 3-0, 4-0.



N PMPI
Membrane Placement Instrument

Micropoint tip for precise placement and manipulation of membranes.



O IM1002
6-Pack Signature Series® Instrument Clip

Secure and maintain open position of scissors, hemostats, and other hinged instruments in lid of cassette.