



HU-FRIEDY SURGICAL KIT

CRESTAL SINUS LIFT KIT



PART CODE	DESCRIPTION
IM4162	16 Instrument Signature Series® Cassette, White
10-130-03	3 Scalpel Handle
P96	9 Molt Periosteal Elevator
CM2/4S6	2/4 Molt DE Curette, #6 Handle
TP43	43 Adson-Brown Tissue Pliers
CSEUWES6	Euwe Bone Scoop/Plugger, #6 Handle
PCP65006	20mm Probe 8-10-12-14-16-18-20, #6 Handle
OSTPAL1	1 Palti Osteotome 2.7mm
OSTPAL2	2 Palti Osteotome 3.2mm
OSTPAL3	3 Palti Osteotome 3.7mm
OSTPAL4	4 Palti Osteotome 4.2mm
PFIWDS1	1 Woodson DE Plastic Filling Instrument
PPR3S6	Prichard Periosteal, #6 Handle
IMCRESTKIT	Complete Crestal Sinus Lift Kit

CRESTAL SINUS LIFT KIT



CRESTAL SINUS LIFT KIT



A IMN4162
16 Instrument Signature Series® Cassette, White

IMS® cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



B 10-130-03
3 Scalpel Handle

Used for cutting gingival tissue and making surgical incisions. Handle is 4.9" long and 13mm wide at base. Compatible with standard scalpel blade sizes.




C P96
9 Molt Periosteal Elevator

Designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. 4.3mm/8.4mm



D CM2/4S6
2/4 Molt DE Curette, #6 Handle

Designed for curettage, cyst removal and tooth socket debridement with 3.6mm and 7mm working ends.



E TP43
43 Adson-Brown Tissue Pliers

Used to grasp and stabilize soft tissue during suturing. 12mm



F CSEUWES6
Euwe Bone Scoop/Plugger, #6 Handle

Designed for easy transport and non-stick compaction of bone graft material. 3.0mm/7.0mm.



G PCP65006
20mm Probe 8-10-12-14-16-18-20, #6 Handle

A double-ended instrument with an explorer on one end and a probe on the other, excellent for dental examinations and charting.



H OSTPAL1
1 Palti Osteotome 2.7mm

Used for preparing bone. 7-10-13-15-18mm



I OSTPAL2
2 Palti Osteotome 3.2mm

Used for preparing bone. 7-10-13-15-18mm



J OSTPAL3
3 Palti Osteotome 3.7mm

Used for preparing bone. 7-10-13-15-18mm



K OSTPAL4
4 Palti Osteotome 4.2mm

Used for preparing bone. 7-10-13-15-18mm



L PFIWDS1
1 Woodson Plastic Filling Instrument

Designed to reflect and retract the mucoperiosteum after incisions of the gingival tissue with 2.2mm and 3.8mm width.



M PPR3S6
Prichard Periosteal, #6 Handle

Secure and maintain open position of scissors, hemostats, and other hinged instruments in lid of cassette.