



HU-FRIEDY SURGICAL KIT

MICROSURGICAL KIT



PART CODE	DESCRIPTION
10-130-70	7 Mini Scalpel Blade Holder
S13K/13KL6	13K/13KL DE Surgical Chisel
C36/376	36/37 Back-Action DE Chisel
8-905DD	Precision Dressing Forceps, Diamond Dusted
NH5024CORT	Cortellini Micro Castro Needle Holder
PPBUSER6	Buser Periosteal, #6 Handle
SPV	Microsurgical Scissors 18cm
IM8120	12 Instrument Signature Series® Cassette, Red
PCPUNC156	CC SE UNC Probe, 1-15, #6 Handle
MIR5	5 FS CS Mouth Mirror
MH6	Cone Socket Mirror Handle, #6 Handle
IMMICSUR	Complete Microsurgical Kit

MICROSURGICAL KIT



MICROSURGICAL KIT



- A 10-130-70**
7 Mini Scalpel Blade Holder
Designed to hold micro scalpel blades for minimally invasive procedures and microsurgery



- B S13K/13KL6**
13K/13KL DE Surgical Chisel
Used to recontour periodontal bone or to remove remaining tooth fibers during periodontal surgery with curved working ends 3.2mm and 4.9mm wide.



- C C36/376**
36/37 Back-Action DE Chisel
Designed for removing bone adjacent to the tooth without causing trauma and is especially useful on the distal of last molars with working ends 3.8mm wide.



- D 8-905DD**
Precision Dressing Forceps, Diamond Dusted
Used to grasp and stabilize soft tissue, 7" length.



- E NH5024CORT**
Cortellini Micro Castro Needle Holder
Used to grasp and guide the needle when suturing.



- F PPBUSER6**
Buser Periosteal, Satin #6 Handle
Designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. 3.8mm/4.3mm



- G SPV**
Microsurgical Scissors
Microsurgical scissors unique design features allow extremely smooth cutting of fine and coarse tissue with reduced trauma. 18cm



- H IM8120**
12 Instrument Signature Series® Cassette, Red
IMS® cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



- I PCPUNC156**
CC SE UNC Probe, 1-15, #6 Handle
Single-ended, color-coded probe. ergonomic, Satin Steel™ handle. University of North Carolina.



- J MIR5**
5 FS CS Mouth Mirror
Single-sided Cone Socket mirror. 15/16 in. (24mm) diameter.



- K MH6**
Cone Socket Mirror Handle, #6 Handle
Used with Single and Double-Sided Cone Socket Mirrors. The blunt end of the mirror handle is ideal for testing sensitivity to percussion and assessing tooth mobility.