



CHECK IT →

SWISS INSTRUMENTS^{PM}

PIEZON® INSTRUMENT TESTER →

EVEN THE BEST ULTRASONIC INSTRUMENTS ARE SUBJECT TO NORMAL WEAR

- > 2 mm shorter = 50% performance loss

REGULAR REPLACEMENT DELIVERS MEASURABLE TREATMENT EFFICACY

- > Easy removal of calculus = time savings
- > Effortless treatment = more clinician comfort
- > Gentle on the tissue = increased patient comfort

PLACE TIP ON CORRESPONDING OPENING TO CHECK FOR WEAR →

