

Water Usage Chart

DentaPure™ DP365B

Independent Water Bottle Cartridge



INSTRUCTIONS

Each box below represents one liter of water. Check off a box every time a liter of water is consumed. When the chart is full, it is time to replace your DentaPure™ Cartridge.

The DentaPure™ DP365B Cartridge is EPA registered to deliver compliant* treatment water for 365 days, at which time you should change your cartridge.

However, if your practice prefers to keep water usage records, the DentaPure DP365B Cartridge will treat 240 liters of water. The chart below is a convenient tool to track your water consumption.

One chart should be used per operatory.

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40
41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60
61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80
81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
101	102	103	104	105	106	107	108	109	110	111	112	113	114	115	116	117	118	119	120
121	122	123	124	125	126	127	128	129	130	131	132	133	134	135	136	137	138	139	140
141	142	143	144	145	146	147	148	149	150	151	152	153	154	155	156	157	158	159	160
161	162	163	164	165	166	167	168	169	170	171	172	173	174	175	176	177	178	179	180
181	182	183	184	185	186	187	188	189	190	191	192	193	194	195	196	197	198	199	200
201	202	203	204	205	206	207	208	209	210	211	212	213	214	215	216	217	218	219	220
221	222	223	224	225	226	227	228	229	230	231	232	233	234	235	236	237	238	239	240

Date Installed: _____ Installed by: _____

*Compliant with the current EPA standard for potable water of <500 CFU/mL.
All claims made when used with potable water.
Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted.
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. DDOC00598 RevB 0522