

TECHNICAL DATA SHEET

<u>PRODUCT</u>	Crosstex Self-Sealing Autoclave Bags	
<u>MATERIALS</u>	Paper	Medical grade paper
	Adhesive Tape	Polyester film tape, water based adhesive
	Indicators	Color process indicators – steam directly on the bag. Water based flexographic ink.
<u>PRODUCT/ PACKAGING REQUIREMENTS</u>	Unique hospital double-fold seal with pull-off adhesive tab. Gusset feature provides 50% more inner space. Item # OC – 2.5" x 1.5" Gusset x 10.5" – White/Blue – 1000/Carton Item # OCS – 2.5" x 1.5" Gusset x 8.75" – White/Pink – 1000/Carton	
<u>STORAGE REQUIREMENTS</u>	Store in a Cool Dry place	
<u>SHELF LIFE</u>	N/A	