

# Lab Transport Bags

**REF BAGS**

## INSTRUCTIONS FOR USE

1. Insert lab work into longer rear pouch.
2. Insert Rx into shorter front pouch.
3. Remove liner and carefully press along adhesive to close.

*This product conforms to OSHA guidelines for containing and transporting laboratory cases.*



**CROSSTEX**



Manufactured for: Crosstex International, Inc.  
6789 W. Henrietta Road, Rush, NY 14543 USA  
585.359.0130 | 631.582.6777 | crosstex.com

Made in USA

BAGS-IFU / REV B - 09/2021

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.

©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved.



Not made with natural rubber latex.

**HuFriedyGroup**  
The Best In Practice