


Safe-Flo™ Adapter

REF M1020ADAP

INSTRUCTIONS FOR USE

1. Insert the end of the saliva ejector into the blue end of the adapter.



2. Insert the valve/saliva ejector combination onto the open end of the high volume evacuation suction line.
3. Bend the saliva ejector tube to the desired position and then place into the mouth.
4. After single patient use, dispose of the valve/ejector combination in the standard trash.

CROSSTEX



Crosstex International, Inc.
6789 W. Henrietta Road, Rush, NY 14543 USA
585.359.0130 | 631.582.6777 | crosstex.com



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



GH REP STERIS GmbH c/o BDO AG,
Langfeldweg 116A,
CH-2504 Biel-Bienne,
SWITZERLAND



For professional use only.
Not made with natural rubber latex.

Made in USA

M1020ADAP-IFU-EN / REV D - 03/2022

In case of any serious incident that has occurred and involves our devices within the European Union, please report it to Hu-Friedy Mfg. Co., LLC and your national competent authority of the member state.

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.

US Pat. No. 8,256,464, 8,714,200, 8,752,579 and foreign counterparts.

©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved.