


SPARKLE™ AIR/WATER SYRINGE TIPS



Locking indicator and groove securely installs the tip




Beveled edge increases patient comfort



Quad Core™ Air Channel ensures air/water separation



Rounded non-scratch end prevents o-ring damage



Available in 6 colors

Reduce patient risk and increase workflow productivity.

Studies have shown that the irregular inner surfaces of traditional metal air/water syringe tips allow for the accumulation of deposits and possibly bioburden.¹ This poses a challenge for reliable sterilization as the inside channel is hard to reach and cannot be visually inspected; potentially putting patients at risk for cross-contamination. Single-use Sparkle™ Air/Water Syringe Tips help eliminate the risks associated with reusable syringe tips as well as improve daily workflow and patient comfort.

- **Save time:** No need to reprocess metal syringe tips between patients.
- **Single-use:** Disposable tips eliminate the risk of cross-contamination.
- **No adapter required:** Locking groove ensures tip is securely installed.
- **Efficient:** Ensures distinct air/water separation for every procedure.
- **Durable:** Rigid construction prevents bending during cheek retraction.
- **Translucent colors:** Allow for visualization of water flow and/or blockages.

SKU	Description	Quantity
BCSAWSCL	Clear	250/Bag
BCSAWSWH	White (opaque, not translucent.)	250/Bag
BCSAWS	Assorted (Blue, Purple, Yellow and Amber)	250/Bag
BCSAWSB	Assorted Bulk (Blue, Purple, Yellow and Amber)	1500/Bag

Learn more at [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

1) Puttaiah, P., Cottone, J., Guildersleeve, J., Azmoudeh, A., Tenney, J. "Rational for Using Single-Use Disposable Air/Water Syringe Tips." Compendium, November 1999, Vol. 20 No.11:1-7.

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted.

©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. T01290 Rev A 0322

Crosstex is now a proud member of

HuFriedyGroup
The Best In Practice