

Comfort Plus™ Premium Saliva Ejectors




Make your patient's dental experience more comfortable.

Comfort Plus™ Premium Saliva Ejectors are engineered with long, rounded edges, smooth surfaces, and strategically placed suction ports with curved openings to protect delicate mucosal tissue and optimize fluid and light debris removal.

- **Increases Patient Comfort:** Suction ports positioned along the side rather than wrapping around the end of the tip to minimize trauma to mucosal tissue.
- **Prevents Blocked Suction:** Engineered to optimize suction for uninterrupted fluid and light debris removal.
- **Sturdy and Robust:** Holds shape when formed and placed in mouth.



Perfectly placed suction ports take in fluid, not tissue.



Comfort Plus™ Premium Saliva Ejectors

ZWWCP	White with White Tip	100/Bag
ZCWCP	Clear with White Tip	100/Bag
ZWGCPM	White with Green Tip Mint Scented	100/Bag
ZWPCPBG	White with Pink Tip Bubble Gum Scented	100/Bag

Call your preferred distributor to order, or learn more at [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

Not made with natural rubber latex.

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, their affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.

©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-842 / 0923

Crosstex is now a proud member of

 **HuFriedyGroup**
The Best In Practice